



Antibody & TDP-43 RNA Aptamer Dual IHC Stain_TS_v1.0

Day1: Immunohistochemistry – Antibody staining (Novolink Polymer Detection kit)

Things to do while you're waiting for your timer

- ✓ Check DPX and put a pipette in ready for later
- ✓ Get the kit from the cold room

- ✓ Do you need to make antigen retrieval solution?
Check you have enough.
- ✓ If using, dilute citric acid retrieval 1:10 (stock: dH₂O)

- ✓ Make up 1:10 TBS (Tris:saline)
Tris is in the fridge, saline is in the cupboard

- ✓ Make up DAB
50µl chromogen + 1ml Dab substrate buffer

The steps

Deparaffinise slides

Xylene	3 min	<input type="checkbox"/>
Xylene	3 min	<input type="checkbox"/>
Alcohol	2 min	<input type="checkbox"/>
Alcohol	2 min	<input type="checkbox"/>

Place slides in tap water (but see below.

NB: Do not place slides in tap water if employing optional picric acid step below

(Optional step) Remove Formalin Pigment

Saturated Alcoholic Picric acid	15 min	<input type="checkbox"/>
Wash in warm, running tap water	15 min	<input type="checkbox"/>

Antigen retrieval

Put 500ml dH ₂ O in the pressure cooker before starting	15 min	<input type="checkbox"/>
--------------------------------------------------------------------	--------	--------------------------

Coverplate

Drops

dH ₂ O	5 min	<input type="checkbox"/>
Peroxidase block (3 drops)	30 min	<input type="checkbox"/>
Wash with TBS	5 min	<input type="checkbox"/>
Protein block (3 drops)	15 min	<input type="checkbox"/>
Make up antibody (antibody:TBS)		<input type="checkbox"/>
Wash with TBS	5 min	<input type="checkbox"/>
Primary antibody	30 min	<input type="checkbox"/>
Wash with TBS	5 min	<input type="checkbox"/>
Post primary block (2 drops)	30 min	<input type="checkbox"/>
Wash with TBS	5 min	<input type="checkbox"/>
Novolink Polymer (2 drops)	30 min	<input type="checkbox"/>
Wash with TBS	5 min	<input type="checkbox"/>
DAB	5 min	<input type="checkbox"/>
dH ₂ O at room temp	Overnight	<input type="checkbox"/>



Antibody & TDP-43 RNA Aptamer Dual IHC Stain_TS_v1.0

Day 2: Immunohistochemistry – TDP-43 Aptamer

Things to do while you're waiting for your timer

The steps

Antigen retrieval (citric acid pH6)	15 min	<input type="checkbox"/>
Coverplate		<input type="checkbox"/>
dH ₂ O	5 min	<input type="checkbox"/>
Peroxidase block (3 drops)	30 min	<input type="checkbox"/>
Wash with TBS	5 min	<input type="checkbox"/>
Avidin block (4 drops)	15 min	<input type="checkbox"/>
Wash with TBS	5 min	<input type="checkbox"/>
Biotin Block (4 drops)	15 min	<input type="checkbox"/>
Wash with TBS	5 min	<input type="checkbox"/>
Make up aptamer in MilliQ H ₂ O		<input type="checkbox"/>
Wash with dH ₂ O	5 min	<input type="checkbox"/>
Aptamer incubation (1in500)	3 hours	<input type="checkbox"/>
Fix step 4% PFA in dH ₂ O	Overnight	<input type="checkbox"/>

Day3: Immunohistochemistry – TDP-43 Aptamer

	Re-coverplate		<input type="checkbox"/>
	Wash with dH ₂ O	5 min	<input type="checkbox"/>
✓ Make up Fast Red Chromogen (1:60)	Anti-Biotin/Alk Phos 1:200 in MilliQ	30 min	<input type="checkbox"/>
1 µl Fast Red + 59 µl Fast Red substrate buffer So ...	Wash with dH ₂ O	5 min	<input type="checkbox"/>
17 µl chromogen + 983 µl Fast Red substrate buffer	Fast Red (1:60)	5 min	<input type="checkbox"/>
	Wash with dH ₂ O	5 min	<input type="checkbox"/>
	Wash well in running tap water		<input type="checkbox"/>
	Counterstain with haematoxylin		<input type="checkbox"/>
	Blue in lithium carbonate		<input type="checkbox"/>
	Dehydrate, clear and mount		<input type="checkbox"/>

References for citation of this method

To cite this dual staining method, please cite both Rifai et al., 2024 and Waldron et al., 2024 (this SOP) below:

The citation, Rifai et al., 2024, for the first publication employing dual staining for NEK1 and TDP-43 in *Brain Pathology*:

Clinicopathological analysis of NEK1 variants in amyotrophic lateral sclerosis. Olivia M. Rifai, Fergal M. Waldron, Danah Sleibi, Judi O'Shaughnessy, Danielle J. Leighton*, Jenna M. Gregory** (2024). *Brain Pathology* doi: 10.1111/bpa.13287 (in press at the time of publication of this SOP) *equal contributions, †corresponding author

The citation, Waldron et al., 2024, for this SOP published on *protocols.io* is:

Antibody and TDP-43 RNA aptamer dual staining to detect patterns of co-pathology in FFPE-preserved human tissue, as described in Rifai et al., 2024 (Brain Pathology): An SOP and tick-sheet. Fergal M. Waldron, Olivia M. Rifai, Jenna M. Gregory†, (2024). *Protocols.io* 2024; DOI: dx.doi.org/10.17504/protocols.io.14egn6wnm15d/v1, †corresponding author.