Protocol Summary

Procedure: RUO DISCOVERY Universal (v0.00.0394)

DISCOVERY ULTRA

Duke University Medical Center, 100 Trent Drive Durham, NC 27710

Validated: No Active: Yes

Protocol No Protocol Name Version Creation Date
181 EZH2HDAC2 Std CC1 Rb HQ manual HRP DAB 1 04/19/2023 3:33:23 PM

- 1 First Sequence [Selected]
- 2 Deparaffinization [Selected]
- 3 Warmup Slide to [70 Deg C], and Incubate for [8 Minutes] (Cycle 1)
- 4 Incubate for [8 Minutes] (Cycle 2)
- 5 Incubate for [8 Minutes] (Cycle 3)
- 6 Pretreatment [Selected]
- 7 Cell Conditioning [Selected]
- 8 CC1 Reservoir [Selected]
- 9 Warmup Slide to [99 Deg C], and Incubate for 4 Minutes (Cell Conditioner #1)
- 10 CC1 8 Min [Selected]
- 11 CC1 16 Min [Selected]
- 12 CC1 24 Min [Selected]
- 13 CC1 32 Min [Selected]
- 14 CC1 40 Min [Selected]
- 15 CC1 48 Min [Selected]
- 16 CC1 56 Min [Selected]
- 17 Inhibitor [Selected]
- 18 DISCOVERY Inhibitor [Selected]
- 19 Apply One Drop of DISC Inhibitor, and Incubate for [8 Minutes]
- 20 Inhibitor CM [Selected]
- 21 Apply One Drop of Inhibitor CM, and Incubate for [8 Minutes]
- 22 Antibody [Selected]
- 23 1st Antibody Manual Application [Selected]
- 24 Warmup Slide to [36 Deg C] from Very Low Temperatures (Primary Antibody)
- 25 Hand Apply (Primary Antibody), and Incubate for [60 Minutes]
- 26 Linking Antibody [Selected]
- 27 2nd Antibody [Selected]
- 28 Warmup Slide to [36 Deg C] from Very Low Temperatures (2nd Antibody)
- 29 Apply One Drop of [Anti-Rabbit HQ] (Detection #1), and Incubate for [0 Hr 12 Min]
- 30 Enzyme conjugate [Selected]
- 31 Apply One Drop of [Anti-HQ HRP] (Conjugate #1), and Incubate for [12 Minutes]
- 32 DAB [Selected]
- 33 Counterstain [Selected]
- 34 Use RB for Counterstain [Selected]
- 35 Apply One Drop of [HEMATOXYLIN] (Counterstain), and Incubate for [4 Minutes]
- 36 Post Counterstain [Selected]
- 37 Use RB for Post Counterstain [Selected]
- 38 Apply One Drop of [BLUING REAGENT] (Post Counterstain), and Incubate for [4 Minutes]