

## Activation Induced Marker (AIM) Staining Protocol

### Materials

- Brilliant Stain Buffer Plus
- CD40 Antibody, anti-human, pure functional grade

### Flow Antibodies Needed:

Membrane Antibody	Fluorochrome	Clone/Vendor/Catalog #	Vol. Per Test (µl)
CXCR5	BUV395	RF8B2/BD/740266	1
CD8	BUV496	RPA-T8/BD/612942	2
CD3	BUV805	UCHT1/BD/612895	2
CD45RA	BV421	HI100/Biolegend/304130	2
LIVE/DEAD	eFluor 506	eBioscience/65-0866-18	0.5
CD14	V500	M5E2/BD/561391	2
CD19	V500	HIB19/BD/561121	2
CD4	BV605	RPA-T4/BD/562658	4
CD38	BV786	HIT2/BD/563964	4
CCR7	FITC	G043H7/Biolegend/353216	2
CD40L	PerCP-ef710	24-31/eBioscience/46-1548-42	4
CD69	PE	FN50/BD/555531	10
PD-1	PE-Dazzle594	EH12.2H7/Biolegend/329940	2
OX40	PE-Cy7	Ber-ACT35/Biolegend/350012	2
CD137	APC	4B4-1/Biolegend/309810	4
HLA-DR	AF700	LN3/ebioscience/56-9956-42	4

CD40	X (unconjugated)	(Milty Biotech, 130-094-133)	1.5
Brilliant Stain Buffer Plus		BD Horizon/566385	10

Prepared Individual Peptides or Peptide Pools

Stimuli	Stock Concentration	Final Concentration
Peptide Pool	1mg/mL	Assay Dependent
DMSO (Negative Control)		Same concentration as peptide
PHA (Positive Control)	1mg/mL	1ug-20ug/mL

## Protocol

### Peptide Stimulation Solution

1. Label U-bottom plate with donor, stimulation solution, name and date.
2. Prepare PHA and DMSO mix separately
3. Prepare and arrange the remaining stimulation solution. Mix thoroughly by pipetting up and down before adding to the experimental plate.
4. Add appropriate stimulus solution to each well in 96-well U-bottom plates.
5. After adding stimulation solutions, prepare an additional solution of anti-CXCR5 antibody as described in below table.

Antibody	Fluorochrome	Clone/vendor/catalog	Amount per well(50ul) (uL)
CXCR5	BUV395	RF8B2/BD/740266	1
HR5			49

6. Add anti-CXCR5 antibody solution to all wells already containing stimulus.
7. Keep plate in incubator at 37°C until cells are ready to be added.

### PBMC Counting and Stimulus Preparation

8. Obtain indicated number of vial(s) of PBMCs.
9. For each donor, prepare sterile 50ml tubes with 10ml HR5 and 20μL Benzonase per vial to be thawed.
10. Thaw PBMC vials.
11. Centrifuge @ 1200 rpm for 7 min.
12. Resuspend cells in HR5 and determine cell number.
13. Centrifuge @ 1200rpm for 7 min.
14. While sample(s) spinning, prepare the CD40 antibody solution the stock for all donors.
15. Resuspend each donor at 1.5ug/mL per 10 million cells per ml in prepared 1.5μg/ml CD40 antibody solution.

16. Incubate the tube for 15 minutes at 37°C/5% CO<sub>2</sub>.
17. Add 100µL of CD40 antibody-treated PMBCs to each well already containing stimulus.
18. Incubate plate for a total of 20-24 hours at 37°C/5% CO<sub>2</sub>.
19. After incubation, spin plate at 1400rpm/4°C/2min.
20. Wash plate by adding 200µL PBS and spinning at 1400rpm/4°C for 2 min.
21. Resuspend cells in 100µL of antibody mix and incubate at 4°C for 30 minutes, protected from light. Wrap plate in aluminum foil and place in fridge.
22. After incubation, add 100µL MACS buffer and spin plate at 1400rpm/4°C/2 min.
23. Wash 1X plate using 200 µL MACs buffer at 1400 rpm/4C/2 min.
24. Wash 1X plate using PBS µL MACs buffer at 1400 rpm/4C/2 min.
25. Resuspend in 120uL PBS.
26. Wrap in foil and store at 4°C until analysis.