

User Guide | CG000582 | Rev E

Xenium In Situ Gene Expression

Probe Hybridization, Ligation & Amplification

For use with:

Xenium Slides & Sample Prep Reagents (2 slides, 2 rxns) PN-1000460

Xenium Decoding Consumables (1 run, 2 slides) PN-1000487

Xenium Pre-Designed Gene Expression Panel - (2 rxns)*

Xenium Add-on Custom Gene Panel - (4 & 16 rxns)*

Xenium Standalone Custom Gene Panel - (4 & 16 rxns)*

Xenium Instrument Accessory Kit Module A PN-1000530

*Refer to the 10x Genomics Support Website for the most current list of available panels and part numbers.

Notices

Document Number

CG000582 | Rev E

Legal Notices

© 2023 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: www.10xgenomics.com/trademarks. 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those thirdparty marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: www.10xgenomics.com/patents. The use of products described herein is subject to 10x Genomics Terms and Conditions of Sale, available at www.10xgenomics.com/legal-notices, or such other terms that have been agreed to in writing between 10x Genomics and user. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

Instrument & Licensed Software Updates Warranties

Updates to existing Instruments and Licensed Software may be required to enable customers to use new or existing products.

Support

Email: support@10xgenomics.com 10x Genomics 6230 Stoneridge Mall Road Pleasanton, CA

Document Revision Summary

Document Number

CG000582

Title

Xenium In Situ Gene Expression User Guide

Revision

Rev E

Revision Date

October 13, 2023

Specific Changes

- Added new Xenium Add-On Custom Gene Panel and Xenium Standalone Custom Gene Panel reaction sizes (4 &16 rxns) (pages 1 and 9).
- Updated Probe Panel Handling guidance (page 21).
- Updated Slide Incubation Guidance (page 28).
- Updated Probe Hybridization Get Started table with information for Xenium Add-On Custom Gene Panels and Xenium Standalone Custom Gene Panels (page 31).
- Added new section for Custom Probe Preparation (page 34).
- Updated resuspension volumes for custom probes (page 34).
- Updated Probe Hybridization Mix volumes for add-on custom and standalone custom probes (page 36).
- Updated thermal cycler lid temperature instructions (pages 44, 48, and 54).
- Updated handling guidance for Xenium Amplification Mix and Xenium Amplification Enzyme (pages 47-48).
- Updated slide storage guidance during long-term incubations (pages 50, 55, and 56).
- Updated Xenium Cassette Lid Cleaning guidance (page 64).

General Changes

Updated for general minor consistency of language and terms throughout.

Table of Contents

Introduction	
Reagent Kits	7
Recommended Thermal Cyclers	11
Additional Kits, Reagents & Equipment	12
Protocol Steps & Timing	14
Stepwise Objectives	15
Tips & Best Practices	
lcons	2
General Reagent Handling	2
Pipette Calibration	2
Probe Panel Handling	2
Xenium Slide Handling	22
Processing a Single Xenium Slide	23
Reagent Addition to Wells	24
Reagent Removal from Wells	24
Xenium Cassette Lid Application & Removal	25
Xenium Cassette Storage	26
Slide Incubation Guidance	27
Tissue Detachment on Xenium Slides	29
Step 1: Probe Hybridization	
1.0 Get Started	3
1.1 Buffer Preparation	33
1.2 Custom Probe Preparation (optional)	34
1.3 Probe Hybridization	35
Step 2: Post Hybridization Wash	
2.0 Get Started	39
2.1 Post Hybridization Wash	40
Step 3: Ligation	
3.0 Get Started	43
3.1 Ligation	44

Table of Contents

Step 4: Amplification	
4.0 Get Started	47
4.1 Amplification	48
4.2 Post Amplification Wash	50
Step 5: Autofluorescence Quenching	
5.0 Get Started	52
5.1 Autofluorescence Quenching	53
5.2 Nuclei Staining	56
Troubleshooting	
Appendix	
Probe Panel Selection	68
Sample Shipping	68

Introduction

Reagent Kits	7
Recommended Thermal Cyclers	11
Additional Kits, Reagents & Equipment	12
Protocol Steps & Timing	14
Stepwise Objectives	15

Reagent Kits

Xenium In Situ Gene Expression Reagent Kits

Refer to SDS for handling and disposal information.

Xenium Slides & Sample Prep Reagents - (2 Slides, 2 Rxns) PN-1000460

Xenium Slides & Sample Prep Reagents (2 slides, 2 rxns), PN-1000460 Store at -20°C				
		#	PN	
•	Xenium Probe Hybridization Buffer	1	2000390	
\circ	Xenium Post Hybridization Wash Buffer	1	2000395	
	Xenium Ligation Buffer	1	2000391	
	Xenium Ligation Enzyme A	1	2000397	
	Xenium Ligation Enzyme B	1	2000398	
•	Xenium Amplification Mix	1	2000392	
•	Xenium Amplification Enzyme	1	2000399	
0	Reducing Agent B	1	2000087	
	Xenium Autofluorescence Mix	1	2000753	
	Xenium FFPE Tissue Enhancer	1	2000798	
	Xenium Nuclei Staining Buffer	1	2000762	
	Perm Enzyme B	1	3000553	
	Xenium Slides (2 pack)	1	3000941	
			10x	


All items, except Xenium FFPE Tissue Enhancer (PN-2000798) and Perm Enzyme B (PN-3000553), are needed for this workflow.

Xenium Decoding Consumables - (1 Run, 2 Slides) PN-1000487

Xenium Decoding Consumables (1 run, 2 slides), PN-1000487	
Store at ambient temperature	
	# PN
Xenium Cassette Kit (2 cassettes + 16 lids)	1 1000566
Extraction Tip	1 2000757
Pipette Tips	1 3000866
Xenium Buffer Cap	4 3000949
Xenium Objective Wetting Consumable	1 2000749
Deionized Water (bottle)	1 3001198
Xenium Sample Wash Buffer A (bottle)	1 3001199
Xenium Sample Wash Buffer B (bottle)	1 3001200
Xenium Probe Removal Buffer (bottle)	1 3001201
	10x GENOMICS

Only the Xenium Cassette Kit (2 cassettes + 16 lids) is needed for this workflow.

Xenium Pre-Designed Gene Expression Panel - (2 Rxns)




Example Xenium Pre-Designed Gene Expression Panel kit. Pre-designed panel kits contain one or two tubes - Xenium Probes with or without Xenium Probe Dilution Buffer.

*Refer to the 10x Genomics Support Website for the most updated list of available panels and part numbers.



Use TE Buffer instead of Xenium Probe Dilution Buffer (to be discontinued) without any impact on assay performance. Relevant guidance is provided in the workflow steps.

Xenium Add-on Custom Gene Panel - (4 & 16 rxns)



Example Xenium Add-on Custom Gene Panel kit. Add-on custom panel kits contain one or two tubes - Xenium Probes with or without Xenium Probe Dilution Buffer.

*Refer to the 10x Genomics Support Website for the most updated list of available panels and part numbers.



Use TE Buffer instead of Xenium Probe Dilution Buffer (to be discontinued) without any impact on assay performance. Relevant guidance is provided in the workflow steps.

Xenium Standalone Custom Gene Panel - (4 & 16 rxns)



Example Xenium Standalone Custom Gene Panel kit. Standalone custom panel kits contain one or two tubes - Xenium Probes with or without Xenium Probe Dilution Buffer.

*Refer to the 10x Genomics Support Website for the most updated list of available panels and part numbers.



Use TE Buffer instead of Xenium Probe Dilution Buffer (to be discontinued) without any impact on assay performance. Relevant guidance is provided in the workflow steps.

Xenium Instrument Accessory Kit Module A PN-1000530

Xenium Instrument Accessory Kit Module A PN-1000530 Store at ambient temperature		
,	#	PN
Waste Bottle	1	3000955
Xenium Waste Tip Tray	1	3000957
Xenium Thermocycler Adaptor	1	3000954
		10x GENOMICS

Only the Xenium Thermocycler Adaptor (PN-3000954) is needed for this workflow.

Recommended Thermal Cyclers

Supplier	Description	Part Number
Bio-Rad	C1000 Touch Thermal Cycler with 96-Deep Well Reaction Module (discontinued)	1851197
Analytik Jena	Biometra TAdvanced 96 SG with 96-well block (silver, 0.2 mL) and gradient function	846-x-070-241 (where x=2 for 230 V; 4 for 115 V; 5 for 100 V, 50-60 Hz)
VWR	Gradient thermal cycler, XT ⁹⁶ Gradient, with 96- well gradient block and standard lid	76452-153
Marshall Scientific	MJ Research PTC-200 Thermal Cycler (discontinued)	05434-05

Additional Kits, Reagents & Equipment

10x Genomics has tested all items listed below. These items perform optimally with the assay. **Substituting materials may adversely affect system performance.** For items with multiple options listed, choose option based on availability and preference. Refer to the manufacturer's website for regional part numbers.

Item	Description	Supplier	Part Number (US)
Plastics			
1.5 ml tubes	DNA LoBind Tubes, 1.5 ml	Eppendorf	022431021
	Low DNA Binding Tubes, 1.5 ml	Sarstedt	72.706.700
15 ml tubes	15 ml PP Centrifuge Tubes	Corning	430791
50 ml tubes	Self-Standing Polypropylene Centrifuge Tubes (50 ml), sterile	Corning	430921
Pipette tips	Tips LTS 20UL Filter RT-L20FLR (or equivalent)	Rainin	30389226
	Tips LTS 200UL Filter RT-L200FLR (or equivalent)	Rainin	30389240
	Tips LTS 1ML Filter RT-L1000FLR (or equivalent)	Rainin	30389213
Kits & Reagents			
Nuclease-free Water	Nuclease-free Water (not DEPC treated)	Thermo Fisher Scientific	AM9937
TE Buffer	TE Buffer, TRIS-EDTA, 1X Solution, pH 8.0 (nuclease-free)	Thermo Fisher Scientific	BP24731
PBS	PBS - Phosphate Buffered Saline (10X) pH 7.4, RNase-free	Thermo Fisher Scientific	AM9624
10% Tween-20	Tween 20 Surfact-Amps Detergent Solution (10% solution) (not 100% Tween diluted to 10%)	Thermo Fisher Scientific	28320
	10% Tween 20	Bio-Rad	1662404
Ethanol	Ethyl alcohol, Pure (200 Proof, anhydrous)	Millipore Sigma	E7023-500ML
Equipment			
Pipettes	Pipet-Lite LTS Pipette L-20XLS+	Rainin	17014392
	Pipet-Lite LTS Pipette L-200XLS+	Rainin	17014391
	Pipet-Lite LTS Pipette L-1000XLS+	Rainin	17014382
Mini centrifuge	VWR Mini Centrifuge (or any equivalent mini centrifuge)	VWR	76269-064
Thermomixer			5382000023
Thermoblock	Eppendorf SmartBlock 2.0 mL (or any equivalent Thermoblock)	Eppendorf	5362000035
Additional Materials			
Waterbath (bead bath untested)			

Item	Description	Supplier	Part Number (US)		
Thermal Cycler (see Reco.	mmended Thermal Cyclers)				
Ice bucket					
Vortex					
PCR Sealing Film (optiona	l, if storing slides long-term)				
	for Water Bath (recommended), apure Water System or equivalent				

This list may not include some standard laboratory equipment.

Protocol Steps & Timing

Steps	Timing	Stop & Store
Day 1		
Step 1: Probe Hybridization		
1.1 Buffer Preparation1.2 Custom Probe Preparation (optional)1.3 Probe Hybridization	20 min 10 min 16-24 h (overnight)	
Day 2		
Step 2: Post Hybridization Wash		
2.1 Post Hybridization Wash	35 min	
Step 3: Ligation		
3.1 Ligation	~2 h	
Step 4: Amplification		
4.1 Amplification	~2 h	
4.2 Post Amplification Wash	15 min	4°C overnight or ≤4 days
Step 5: Autofluorescence Quenching		
5.1 Autofluorescence Quenching	45 min	4°C overnight or ≤4 days (in the dark)
5.2 Nuclei Staining	10 min	4°C overnight or ≤4 days (in the dark)*

^{*}Alternatively, slides can be stored for up to 24 days after Nuclei Staining. For long-term storage, the Xenium Cassette Lid should be replaced with a slide seal and slides should be stored at 4°C in the dark. PBS-T storage buffer should be exchanged with fresh PBS-T every 3-4 days. Xenium Cassettes must be covered with a Xenium Cassette Lid or slide seal during storage.



Storing slides for more than 24 days may carry the risk for a lower number of genes or transcripts detected per cell, changes in tissue morphology over time, and microbe growth. These risks are dependent on many factors including input tissue quality and how nuclease-free/microbe-free the workflow is.

Stepwise Objectives

Xenium In Situ Gene Expression assays RNA at the subcellular level by using targeted probes in formalin fixed & paraffin embedded (FFPE) or fresh frozen (FF) tissue sections. FFPE tissue sections placed on Xenium Slides are deparaffinized and decrosslinked as described in Xenium In Situ for FFPE - Deparaffinization & Decrosslinking (Demonstrated Protocol – CG000580). FF tissue sections placed on Xenium slides are fixed and permeabilized as described in Xenium In Situ for Fresh Frozen - Fixation & Permeabilization (Demonstrated Protocol – CG000581).


Pre-designed, add-on custom, or standalone custom probe panels are then added to the tissue. Each circularizable DNA probe contains two regions that hybridize to the target RNA and a third region that encodes a gene-specific barcode. The two ends of the probes bind the target RNA and are ligated to generate a circular DNA probe. Following ligation, the circularized probe is enzymatically amplified, generating multiple copies of the gene-specific barcode for each RNA target.

Xenium slides containing FFPE or FF tissue sections are then loaded for imaging and analysis on the Xenium Analyzer instrument for high-throughput, automated in situ analysis. Fluorescently labeled oligos bind to the amplified DNA probes. Cyclical rounds of fluorescent probe hybridization, imaging, and removal generate optical signatures specific for each barcode, which are converted into a gene identity. Identified transcripts can be visualized using Xenium Explorer software.

This document outlines the protocol for generating Xenium In Situ Gene Expression data from FFPE and FF tissue sections placed on Sample Areas of a Xenium slide.


Xenium Slide

The Xenium slide has one Sample Area measuring $10.45 \times 22.45 \text{ mm}$ and is defined by a fiducial frame. The imageable area, measuring $12 \text{ mm} \times 24 \text{ mm}$, includes the area within the Sample Area + fiducial frame. FFPE or FF tissue sections are placed within the Sample Area for downstream processing and analysis.




Step 1: Probe Hybridization

Pre-designed, add-on custom, or standalone custom probe panels are added to the FFPE or FF tissue sections. The DNA probes are flanked by two regions that independently hybridize to the target RNA and also contain a gene-specific barcode sequence. The probes hybridize to their complementary target RNA in an overnight incubation.




Step 2: Probe Hybridization Wash

Excess, unbound probes are washed away in the post hybridization wash step.




Step 3: Ligation

After removal of unbound probes, a ligase is added to seal the junction between the probe regions that have hybridized to RNA. Ligation of the probe ends on the targeted RNA region generates a circular DNA probe. This ligation ensures a unique level of probe specificity to the target region.




Step 4: Amplification

The ligation products are enzymatically amplified. Hundreds of copies of the gene-specific barcode are generated during the amplification process.



Step 5: Autofluorescence Quenching

Autofluorescence Quenching diminishes unwanted autofluorescence and enhances signal-to-noise ratio in the treated FFPE and FF tissue sections. Next, nuclei are stained with DAPI (derived from Xenium Nuclei Staining Buffer) to assist in identification of tissue or regions of interest during an instrument overview scan. Finally, tissue sections on Xenium slides assembled into Xenium Cassettes are loaded into the Xenium Analyzer for imaging and decoding.



Tips & Best Practices



Icons



Tips & Best Practices section includes additional guidance



Signifies critical step requiring accurate execution



Troubleshooting section includes additional guidance

General Reagent Handling

- Fully thaw reagents at indicated temperatures. Thoroughly mix reagents before use.
- When pipette mixing reagents, unless otherwise specified, set pipette to 75% of total volume.
- Keep all enzymes and Master Mixes on ice during setup and use, unless otherwise stated.
- Promptly move reagents back to the recommended storage.

Pipette Calibration


• Follow manufacturer's calibration and maintenance schedules.

Probe Panel Handling

- 10x Genomics provides the following types of probe panels: pre-designed, add-on custom, and standalone custom. Add-on custom panels are used to supplement pre-designed panels. Standalone custom probe panels are used alone and do not require pre-designed panels.
- Pre-designed probes are good for two Xenium slides. Add-on and Standalone custom probes are good for either four or sixteen Xenium slides.
- Add-on and Standalone custom probes are delivered lyophilized at room temperature and should be stored at -20°C upon resuspension.
- Custom probes must be resupsended prior to use. See Custom Probe Preparation for more details.
- Record the Custom Panel Design ID and Slide Number before starting the workflow. This information is critical for identifying the correct electronic decode file when setting up the Xenium Analyzer in downstream steps.

Xenium Slide Handling

- Always wear gloves when handling slides.
- The bottom of the slide is indicated by the etched label, which should be readable when in the proper position.
- The tissue sections should always be placed within the Sample Area on etched label side of the slide.
- Hold the slide on the label. DO NOT touch the tissue sections or near fiducials.
- Minimize exposure of the slides to sources of particles and fibers.
- Keep the slide cassette flat on the bench when adding reagents to the Sample Area.
- Ensure that no absorbent surface is in contact with the reagents on the slide during incubation.
- When pipetting reagent onto a slide, avoid generating bubbles. Avoid pipetting directly onto the tissue.
- After aspirating reagent from a slide, pipette new reagent onto same slide before moving onto aspiration of second slide.



Processing a Single Xenium Slide

• Xenium reagent kits are sufficient for two reactions, and for optimal Xenium Analyzer throughput, two slides should be run at the same time.



- It is possible to perform the Xenium In Situ Gene Expression workflow with a single slide. To do this, ensure the following best practices are followed for optimal assay performance:
 - Assemble a mock Xenium Cassette using a blank slide and a cassette from the Xenium Cassette Kit (2 cassettes), PN-1000566.
 - Insert the blank slide into the Xenium Cassette. Cassettes should be assembled following the instructions in Troubleshooting for Xenium Cassette Assembly.
 - Attach a Xenium Cassette Lid from the Xenium Cassette Kit (2 cassettes), PN-1000566 to the cassette containing the blank slide following Tips & Best Practices for Xenium Lid Application. It is not necessary to add liquid to the slide well before adding the lid.
 - For all incubation steps with the thermal cycler lid closed, ensure the mock slide cassette is placed alongside the Xenium slide cassette containing tissue on the Thermocycler Adaptor.

Reagent Addition to Wells

- Place assembled cassette flat on a clean work surface.
- Dispense and remove reagents along the side of the well without touching the tissue sections and without introducing bubbles.



• Always cover the Sample Area completely when adding reagents to the well. A gentle tap may help spread the reagent more evenly.



Reagent Removal from Wells

- Place assembled cassette flat on a clean work surface.
- Slightly tilt the cassette while removing the reagent.
- Place the pipette tip on the bottom edge of the well.
- Remove reagents along the side of the well without touching the tissue sections.
- Remove all liquid from the well in each step.



Xenium Cassette Lid Application & Removal

Application

- Place the Xenium Cassette flat on a clean work surface.
- Hold the Xenium Cassette Lid with index and middle finger on two upper tabs and thumb on the lower clip.
- Align the Xenium Cassette Lid with the surface of the Xenium Cassette. Hook the two upper clips into the two holes on the top of the cassette.
- Push the lid down until the lower clip clicks into place.
- Inspect the lid to confirm placement.

Removal

- Place the Xenium Cassette flat on a clean work surface.
- Push on the top of the two upper tabs with index and middle fingers.
- Use thumb to push in on the lower clip.
- While maintaining inward pressure, pull upward with thumb until the lower clip disengages.
- Ensure that no liquid splashes out of the well.



Xenium Cassette Lids are a single use item and should be discarded after each use. PBS-T washes DO NOT require sealing of the cassette with a lid.

Xenium Cassette Storage

- Store cassettes sealed with a Xenium Cassette Lid at the indicated stopping points listed throughout the protocol and as outlined in the Protocol Steps & Timing on page 14.
- Cassettes should always be stored hydrated with recommended reagent and at 4°C.



Slide Incubation Guidance

Incubation at a specified temperature

• Position a Xenium Thermocycler Adaptor on a thermal cycler that is set at the incubation temperature. Ensure thermal cycler has reached appropriate temperature prior to starting incubation.



- Ensure that the Thermocycler Adaptor is fully inserted into the thermal cycler and is in contact with the thermal cycler block surface uniformly.
- When incubating a slide, position the slide on the Thermocycler Adaptor with the tissue side facing up.





- Ensure the Sample Area is aligned with the corresponding area on the Thermocycler Adaptor. DO NOT close the lid.
- When incubating a slide encased in a cassette, place the assembled unit on the Thermocycler Adaptor with the well facing up. Ensure the cassette is in complete contact with the Thermocycler Adaptor. The cassette should always be sealed with a Xenium Cassette Lid when on the Thermocycler Adaptor unless indicated otherwise.




• After each thermal cycler incubation, at least 400 μ l should be aspirated from cassette well. If less than 400 μ l is recovered, contact support@10xgenomics.com.

Tightening the thermal cycler lid

- Thermal cycler lid contact with the Xenium Cassette Lid is critical for assay performance.
- For thermal cyclers with adjustable lids, tighten the lid until an audible click is heard. Tightening past the click risks breaking the slide.

Incubation at room temperature

- Place the assembled cassette on a flat, clean, non-absorbent work surface.
- Ensure that no absorbent surface is in contact with the reagents on the slide/cassette during incubation.



Tips & Best Practices 28

Tissue Detachment on Xenium Slides



- Monitor section adhesion on Xenium slides throughout the workflow.
- Tissue detachment during the workflow can negatively impact performance. If observed, contact support@10xgenomics.com.
- For more information, refer to Troubleshooting.

Tips & Best Practices 29

Step 1:

Probe Hybridization

1.0 Get Started	31
1.1 Buffer Preparation	33
1.2 Custom Probe Preparation (optional)	34
1.3 Probe Hybridization	35

1.0 Get Started

Items			10x PN	Preparation & Handling	Storage
Equilibrate	to room	temperature			
	•	Xenium Probe Hybridization Buffer	2000390	Thaw at room temperature for 15 min or until completely thawed. Check for precipitate and invert until clear. Maintain at room temperature after thawing. After use, return to -20°C.	−20°C
	•	Xenium Probe Dilution Buffer (to be discontinued)	2000393	Thaw at room temperature. Vortex and centrifuge briefly. Maintain at room temperature after thawing. After use, return to -20°C. Use TE Buffer instead of Xenium Probe Dilution Buffer without any impact on assay performance. See details below.	-20°C
	•	Xenium Pre- Designed Gene Expression Probes*	-	Thaw at room temperature. See Probe Hybridization for additional handling instructions.	-20°C
	0	Xenium Add-on Custom Probes*	-	Resuspend add-on custom probes according to Custom Probe Preparation. For additional handling instructions, see Probe Hybridization.	-20°C (after resuspension)
	0	Xenium Standalone Custom Probes*	-	Resuspend standalone custom probes according to Custom Probe Preparation. For additional handling instructions, see Probe Hybridization.	-20°C (after resuspension)
Obtain					
		Assembled cassettes containing FFPE or FF tissue samples	-	Consult Xenium in Situ for FFPE - Deparaffinization & Decrosslinking (Demonstrated Protocol CG000580) or Xenium in Situ for Fresh Frozen - Fixation & Permeabilization (Demonstrated Protocol CG000581), respectively.	-
		Nuclease-free Water	-	-	Ambient

Items		10x PN	Preparation & Handling	Storage
	10X PBS, pH 7.4	-	-	Ambient
	10% Tween-20	-	-	Ambient
	Heatblock or waterbath	-	Preheat to 95°C.	Ambient
	Xenium Cassette Lids (16 ct)	3001046	See Tips & Best Practices.	Ambient
	Xenium Thermocycler Adaptor	3000954	See Tips & Best Practices.	Ambient
	TE Buffer, TRIS- EDTA, 1X Solution, pH 8.0 (nuclease-free)	-	The pH of the stock solution should be 8.0. Readjusting the pH is NOT recommended.	Ambient



*Thaw appropriate probe panels based on experimental needs. Refer to the 10x Genomics Support Website for the most updated list of available panels and part numbers.

Step 1: Probe Hybridization 32

1.1 Buffer Preparation

Prepare the following buffers fresh before starting the Xenium In Situ Gene Expression workflow. The volumes of each buffer are sufficient for washes in all subsequent steps.

a. Prepare 1X PBS according to the table below before use and maintain at **room temperature.** Add reagents in the order listed. Invert gently to mix.

1X PBS	Stock	Final	1X+10% (ml)	2X+10% (ml)
Nuclease-free Water	-	-	13.5	27.0
10X PBS, pH 7.4	10X	1X	1.5	3.0
Total	-	-	15.0	30.0

b. Using 1X PBS from step 1.1a, prepare PBS-Tween Buffer (PBS-T) according to the table below before use and maintain at **room temperature.** Add reagents in the order listed. Invert gently to mix.

PBS-T	Stock	Final	1X+10% (ml)	2X+10% (ml)
1X PBS (prepared at Step 1.1a)	-	-	9.95	19.9
10% Tween-20	10%	0.05%	0.05	0.1
Total	-	-	10.0	20.0

Step 1: Probe Hybridization 33

1.2 Custom Probe Preparation (optional)

Proceed to Probe Hybridization, step 1.3, directly if using pre-designed probes only. Add-on and standalone custom probes are delivered lyophilized and must be resuspended before use. Resuspend add-on or standalone custom probes according to the instructions below before proceeding with Probe Hybridization.



Confirm the number of reactions provided for the add-on or standalone custom probes prior to resuspension.

- a. Centrifuge custom probe panel tube briefly.
- **b.** Resuspend lyophilized custom probes in TE Buffer according to the following table.

Custom Probe Resuspension	10x PN	TE Buffer (μl)
4 reactions/kit		
Xenium Add-on	varies	140
Xenium Standalone	varies	140
16 reactions/kit		
Xenium Add-on	varies	700
Xenium Standalone	varies	700

Xenium Probe Dilution Buffer may be used interchangeably with TE Buffer if sufficient volume is available.

- **c.** Replace the cap firmly and agitate on vortex mixer at **room temperature** for **5 min**.
- **d.** Centrifuge custom probe panel tube briefly and maintain at **room temperature**.



If custom probes are already resuspended, thaw at **room temperature** prior to starting Probe Hybridization.

Custom probes resuspended in Xenium Probe Dilution Buffer or TE Buffer can be stored at -20°C until the expiration date specified on the kit.

1.3 Probe Hybridization



Before starting this protocol, ensure that tissue sections have been appropriately deparaffinized and decrosslinked if working with FFPE tissues. Ensure that tissue sections have been appropriately fixed and permeabilized if working with fresh frozen tissues. Consult Xenium In Situ for FFPE - Deparaffinization & Decrosslinking Demonstrated Protocol (CG000580) or Xenium In Situ for Fresh Frozen - Fixation & Permeabilization Demonstrated Protocol (CG000581), respectively, for more information.



During reagent removal steps, ensure that **ALL** the liquid is removed from the wells. See Tips & Best Practices for guidance on Reagent Removal.

a. Prepare a thermal cycler with the following incubation protocol and start the program.

Lid Temperature	Reaction Volume	Run Time
50°C	100 μΙ	-
Step	Temperature	Time hh:mm:ss
Pre-equilibrate	50°C	Hold
Probe Hybridization	50°C	Overnight (16 - 24 h)
Hold	50°C	Hold

- **b.** Obtain probes that have been thawed or equilibrated to **room temperature**. Pulse spin the probe tube and remove an aliquot appropriate for the number of desired Xenium slides (**34 µl** per slide).
- **c.** Preheat probes by incubating for **2 min** at **95°C** in a heatblock or waterbath, followed by **1 min** on **ice**.
- **d.** Equilibrate to room temperature, centrifuge briefly, and maintain at **room temperature**.
- **e.** Prepare Probe Hybridization Mix according to the three options below. Preparation instructions for each option can be found on the following page.
 - Pre-designed probe panels only
 - Add-on custom probe panels used with pre-designed probe panels
 - Standalone custom probe panels only

Prepare Probe Hybridization Mix shortly before use and maintain at **room temperature**. Add reagents in the order listed. Pipette mix and centrifuge briefly.

Option 1: Probe Hybridization Mix (pre-designed probe panels only)

Probe Hybridization Mix (pre-designed probe panels only)	10x PN	1X+5% (μl)	2X+5% (μl)
Xenium Probe Hybridization Buffer	2000390	315.0	630.0
TE Buffer	-	177.0	354.0
Xenium Pre-Designed Gene Expression Probes*	-	33.0	66.0
Total	-	525.0	1,050.0

Option 2: Probe Hybridization Mix (add-on custom probe panels used with pre-designed probe panels)

	Probe Hybridization Mix (add-on custom probe panels used with pre- designed probe panels)	10x PN	1X+5% (µl)	2X+5% (µl)
	Xenium Probe Hybridization Buffer	2000390	315.0	630.0
	TE Buffer	-	144.0	288.0
	Xenium Pre-Designed Gene Expression Probes*	varies	33.0	66.0
\circ	Xenium Add-on Custom Probes*†	varies	33.0	66.0
	Total	-	525.0	1,050.0

Option 3: Probe Hybridization Mix (standalone custom probe panels only)

	Probe Hybridization Mix (standalone custom probe panels only)	10x PN	1X+5% (μl)	2X+5% (μl)
	Xenium Probe Hybridization Buffer	2000390	315.0	630.0
	TE Buffer	-	177.0	354.0
0	Xenium Standalone Custom Probes*†	varies	33.0	66.0
	Total		525.0	1,050.0

^{*}Refer to the 10x Genomics Support Website for the most updated list of available panels and part numbers.

[†]Custom probes resuspended in Xenium Probe Dilution Buffer or TE Buffer can be stored at -20°C until the expiration date specified on the kit.



Record the Custom Panel Design ID and Slide Number before starting workflow. This information is critical for identifying the correct electronic decode file when setting up the Xenium Analyzer in downstream steps.

Continued from step e of Probe Hybridization 1.3

- **f.** Retrieve the assembled Xenium Cassette containing FFPE or fresh frozen tissue sections.
- **g.** Remove all PBS-T from FFPE or fresh frozen tissues as prepared according to Xenium in Situ for FFPE Deparaffinization & Decrosslinking (CG000580) or Xenium in Situ for Fresh Frozen Fixation & Permeabilization (CG000581) Demonstrated Protocols, respectively.
- **h.** Add **500 μl** room-temperature Probe Hybridization Mix along the side of the well to uniformly cover the tissue sections, without introducing bubbles.
- i. Apply a new Xenium Cassette Lid on the Xenium Cassette and place on the Xenium Thermocycler Adaptor on the pre-heated thermal cycler. Tightly close the thermal cycler lid until an audible click is heard.



- Audible click will be heard in thermal cyclers with an adjustable lid (i.e. Bio-Rad C1000 Touch Thermal Cycler).
- **j.** Skip Pre-equilibrate step to initiate Probe Hybridization.
- k. After Probe Hybridization is complete, immediately proceed to next step.

Step 1: Probe Hybridization 10xgenomics.com 37

Step 2:

Post Hybridization Wash

2.0 Get Started	39
2.1 Post Hybridization Wash	40



2.0 Get Started

Items			10x PN	Preparation & Handling	Storage
Equilibrate	to room	temperature			
		Xenium Post Hybridization Wash Buffer	2000395	Thaw at room temperature for 30 min or until thawed completely. Vortex and centrifuge briefly. Keep the buffer at room temperature after thawing.	−20°C
Obtain					
		PBS-T	-	Prepared at Step 1.1.	Ambient

2.1 Post Hybridization Wash

a. Remove the Xenium Cassette from the Thermocycler Adaptor and place on a flat, clean work surface. DO NOT let the cassette cool down before proceeding to PBS-T washes.



Fluid on the Thermocycler Adaptor may indicate a reagent leak from the cassette. See Troubleshooting for more details.

b. Remove the Xenium Cassette Lid. Using a pipette, remove all Probe Hybridization Mix from well corners. Discard old Cassette Lids.

If processing two slides at a time, remove and add reagent from first slide before proceeding to second slide. Ensure tissue sections are always covered with reagent in between removal and addition steps to avoid drying out of tissue samples. At least 400 µl of liquid should be aspirated from each slide.

c. Immediately add 500 µl PBS-T prepared at step 1.1 along the side of the well to uniformly cover the tissue sections, without introducing bubbles. Removal and addition of buffers should be done quickly to prevent drying of tissue sections.



Small bubbles on the surface of the slide are normal and unlikely to compromise assay performance. DO NOT aspirate or pop bubbles, as this can lead to detachment or scratching of the tissue.

- **d.** Incubate for **1 min** at **room temperature**.
- e. Prepare a thermal cycler with the following incubation protocol and start the protocol.

Lid Temperature	Reaction Volume	Run Time
37°C (lid may be set to the lowest temperature if the instrument does not enable 37°C)	100 μΙ	-
Step	Temperature	Time hh:mm:ss
Pre-equilibrate	37°C	Hold
Post Hybridization Wash	37°C	00:30:00
Hold	37°C	Hold

- **f.** Using a pipette, remove all PBS-T from well corners.
- **g.** Add **500 μl** PBS-T.
- **h.** Incubate for **1 min** at **room temperature**.
- i. Remove all PBS-T.

- j. Add $500 \mu l$ Xenium Post Hybridization Wash Buffer to the well.
- k. Apply a new Xenium Cassette Lid on the Xenium Cassette and place on the Thermocycler Adaptor on the pre-heated thermal cycler. Close the thermal cycler lid.
- 1. Skip Pre-equilibrate step to initiate Post Hybridization Wash.



- Start thawing Ligation reagents during Post Hybridization Wash incubation as outlined in the Get Started table in step 3.0.
- m. After the Post Hybridization Wash is complete, immediately proceed to the next step.

Step 3:

Ligation

3.0 Get Started	43
3.1 Ligation	44



3.0 Get Started

Items		10x PN	Preparation & Handling	Storage
Equilibrate to room	temperature	•		
	Xenium Ligation Buffer	2000391	Thaw at room temperature for 15 min or until completely thawed. Vortex and centrifuge briefly. Maintain at room temperature after thawing.	-20°C
Place on ice				
	Xenium Ligation Enzyme A	2000397	Pipette mix and centrifuge briefly. Maintain on ice until ready to use.	-20°C
	Xenium Ligation Enzyme B	2000398	Pipette mix and centrifuge briefly. Maintain on ice until ready to use.	-20°C
Obtain				
	PBS-T	-	Prepared at Step 1.1.	Ambient

Step 3: Ligation 10xgenomics.com 43

3.1 Ligation

a. Prepare Ligation Mix shortly before using. Add reagents in the order listed. Pipette mix 10X and centrifuge briefly. Maintain on ice.

Ligation Mix	10x PN	1X+10% (μl)	2X+10% (μl)
Xenium Ligation Buffer	2000391	481.2	962.5
Xenium Ligation Enzyme A	2000397	13.8	27.5
Xenium Ligation Enzyme B	2000398	55.0	110.0
Total	-	550.0	1,100.0

- **b.** Remove the Xenium Cassette from the Thermocycler Adaptor and place on a flat, clean work surface.
- **c.** Remove the Xenium Cassette Lid. Using a pipette, remove all Xenium Post Hybridization Wash Buffer from the well. Discard old Cassette Lids.

If processing two slides at a time, remove and add reagent from first slide before proceeding to second slide. Ensure tissue sections are always covered with reagent in between removal and addition steps to avoid drying out of tissue samples. At least 400 μ l of liquid should be aspirated from each slide.

- **d.** Immediately add 500 μ l PBS-T prepared at step 1.1 to the well. Removal and addition of buffers should be done quickly.
- e. Incubate at room temperature for 1 min.
- **f.** Prepare a thermal cycler with the following incubation protocol. Place a Thermocycler Adaptor on the thermal cycler and start the program.

Lid Temperature	Reaction Volume	Run Time
37°C (lid may be set to the lowest temperature if the instrument does not enable 37°C)	100 μΙ	-
	_	Time
Step	Temperature	hh:mm:ss
Step Pre-equilibrate	37°C	hh:mm:ss Hold

- g. Using a pipette, remove all PBS-T from well corners.
- **h.** Add **500 μl** PBS-T.
- i. Incubate at room temperature for 1 min.

Step 3: Ligation 10xgenomics.com 44

- j. Remove all PBS-T.
- **k. Repeat** steps h-j one more time.
- **1.** Add **500** μ **1** Ligation Mix to the well.
- **m.** Apply a new Xenium Cassette Lid on the Xenium Cassette and place on the Thermocycler Adaptor on the pre-heated thermal cycler. Close the thermal cycler lid.
- **n.** Skip Pre-equilibrate step to initiate Ligation.



Start thawing Amplification reagents (except enzymes) during Ligation incubation as outlined in the Get Started table in step 4.0.

o. After Ligation is complete, **immediately** proceed to next step.

Step 3: Ligation 10xgenomics.com 45

Step 4:

Amplification

4.0 Get Started	47
4.1 Amplification	48
4.2 Post Amplification Wash	50



4.0 Get Started

Item		10x PN	Preparation & Handling	Storage
Place on ice				
	Xenium Amplification Mix	2000392	Immediately after initiating Ligation step, thaw on ice. Vortex and centrifuge briefly. Ensure completely thawed before use, with no precipitate remaining.	-20°C
	Xenium Amplification Enzyme	2000399	Transfer to ice before use. Pipette mix and centrifuge briefly. Maintain on ice until ready to use.	-20°C
Obtain				
	PBS-T	-	Prepared at Step 1.1.	Ambient
	TE Buffer, TRIS- EDTA, 1X Solution, pH 8.0 (nuclease- free)	-	-	Ambient

Step 4: Amplification 47

4.1 Amplification

a. Prepare Amplification Master Mix on ice shortly before use. Add reagents in the order listed. Pipette mix 10X and centrifuge briefly. Maintain on ice.

Amplification Master Mix	10x PN	1X +10% (μl)	2X +10% (μl)
Xenium Amplification Mix (Thaw completely before use)	2000392	495.0	990.0
Xenium Amplification Enzyme	2000399	55.0	110.0
Total	-	550.0	1,100.0

- **b.** Remove the Xenium Cassette from the Thermocycler Adaptor and place on a flat, clean work surface.
- **c.** Remove the Xenium Cassette Lid. Using a pipette, remove all Ligation Mix from the well. Discard old Cassette Lids.

If processing two slides at a time, remove and add reagent from first slide before proceeding to second slide. Ensure tissue sections are always covered with reagent in between removal and addition steps to avoid drying out of tissue samples. At least 400 μ l of liquid should be aspirated from each slide.

- **d.** Immediately add 500 μ l PBS-T prepared at step 1.1 to the well.
- e. Incubate for 1 min at room temperature.
- **f.** Prepare a thermal cycler with the following incubation protocol. Place a Thermocycler Adaptor on the thermal cycler and start the program.

Lid Temperature	Reaction Volume	Run Time
30°C (lid may be set to the lowest temperature if the instrument does not enable 30°C)	100 μΙ	-
Step	Temperature	Time hh:mm:ss
Pre-equilibrate	30°C	Hold
Amplification	30°C	02:00:00

- g. Using a pipette, remove all PBS-T from well corners.
- **h.** Add **500 µl** PBS-T.
- i. Incubate for 1 min at room temperature.
- j. Remove all PBS-T.

Step 4: Amplification 10xgenomics.com 48

- **k. Repeat** steps h-j one more time.
- 1. Immediately add 500 μl Amplification Master Mix to the well.
- **m.** Apply a new Xenium Cassette Lid on the Xenium Cassette and place on the Thermocycler Adaptor on the thermal cycler. Close the thermal cycler lid.
- **n.** Skip pre-equilibrate step to initiate Amplification.



Start thawing Autofluorescence Quenching reagents during Amplification incubation as outlined in the Get Started table in step 5.0.

o. After Amplification is complete, immediately proceed to next step.

Step 4: Amplification 10xgenomics.com 49

4.2 Post Amplification Wash

- **a.** Remove the Xenium Cassette from the Thermocycler Adaptor and place on a flat, clean work surface.
- **b.** Remove the Xenium Cassette Lid and using a pipette, remove all Amplification Mix from the well. Discard old Cassette Lids.
- c. Add 500 μ l TE Buffer to the well.
- d. Incubate 1 min at room temperature.
- e. Remove all TE buffer.
- **f. Repeat** steps c-e one more time.
- **g.** Add **500** μ **l** TE Buffer to the well.



h. Store slides in TE Buffer overnight or for ≤4 days at 4°C with a new Xenium Cassette Lid applied on the Xenium Cassette. Alternatively, proceed to next step. If storing slides, DO NOT discard the lid after use and instead save for step 5, Autofluorescence Quenching.



Xenium Cassettes must be covered with a Xenium Cassette Lid during storage.

Step 4: Amplification 10xgenomics.com 50

Step 5:

Autofluorescence Quenching

5.0 Get Started	52
5.1 Autofluorescence Quenching	53
5.2 Nuclei Staining	56



5.0 Get Started

Items			10x PN	Preparation & Handling	Storage			
Equilibrate to room temperature								
		Xenium Autofluorescence Mix	2000753	Thaw in a thermomixer (with 2.0-ml thermoblock) for 15 min at 37°C, 300 rpm with shaking. Cool to room temperature for 5 min. Vortex for 30 sec and centrifuge briefly. Alternatively, thaw in a waterbath for 15 min at 37°C. Cool to room temperature for 5 min. Vortex for 30 sec and centrifuge briefly.*	-20°C			
	\bigcirc	Reducing Agent B	2000087	Thaw at room temperature. Vortex and centrifuge briefly.	-20°C			
	•	Xenium Nuclei Staining Buffer	2000762	Thaw at room temperature. Vortex and centrifuge briefly. Keep in the dark until ready to use.	-20°C			
Obtain								
		Nuclease-free Water	-	-	Ambient			
		1X PBS	-	Prepared at Step 1.1.	Ambient			
		PBS-T	-	Prepared at Step 1.1.	Ambient			
		100% Ethanol	-	-	Ambient			



*Pre-heat thermomixer or waterbath to 37°C in advance of intended use.

5.1 Autofluorescence Quenching

- a. Prepare the following for Autofluorescence Quenching:
 - i. Prepare diluted Reducing Agent B. Add reagents in the order listed. Maintain at room temperature.

	Diluted Reducing Agent B	10x PN	Stock	Final	1X+10% (µl)	2X+10% (μl)
	1X PBS (prepared at Step 1.1)	-	-	-	544.5	1,089.0
0	Reducing Agent B	2000087	-	-	5.5	11.0
	Total	-	-	-	550.0	1,100.0

ii. Prepare 70% Ethanol. Add reagents in the order listed. Maintain at room temperature.

70% Ethanol	10x PN	Stock	Final	1X+10% (µl)	2X+10% (μl)
Nuclease-free Water	-	-	-	330.0	660.0
100% Ethanol	-	100%	70%	770.0	1,540.0
Total	-	-	-	1,100.0	2,200.0

iii. Prepare Autofluorescence Solution using thawed Xenium Autofluorescence Mix prepared according to step 5.0. Add reagents in the order listed and vortex to mix. Maintain at room temperature in the dark until ready to use.

Autofluorescence Solution	10x PN	Stock	Final	1X+10% (μl)	2X+10% (μl)
100% Ethanol	-	100%	-	544.5	1,089.0
Xenium Autofluorescence Mix	2000753	-	-	5.5	11.0
Total	-	-	-	550.0	1,100.0

b. Retrieve the Xenium Cassette from step 4.2h and place on a flat, clean work surface.

- **c.** If stored, remove the Xenium Cassette Lid and using a pipette, remove all TE Buffer from the well. **Save lid** for use in following indicated steps.
 - If processing two slides at a time, remove and add reagent from first slide before proceeding to second slide. Ensure tissue sections are always covered with reagent in between removal and addition steps to avoid drying out of tissue samples. At least 400 μ l of liquid should be aspirated from each slide.
- **d.** Add **1,000 μl** 1X PBS prepared at step 1.1 to the well and incubate for **1 min** at **room temperature**.
- e. Remove all 1X PBS.
- **f. Repeat** steps d-e two more times.
- g. Add 500 μl Diluted Reducing Agent B prepared at step 5.1ai to the well.
- h. Apply a Xenium Cassette Lid on the Xenium Cassette, and incubate for 10 min at room temperature. Use lid from step 5.1c if previously stored. Apply new lid if not stored.
- i. Remove the Xenium Cassette Lid. Using a pipette, remove all Diluted Reducing Agent B from the well. Discard old Cassette Lids.
- j. Add 1,000 µl 70% Ethanol prepared at step 5.1aii. Wait 1 min.
- k. Remove all 70% Ethanol.
- **l.** Add **1,000 μl** 100% Ethanol. Wait **1 min**.
- m. Remove all 100% Ethanol.
- **n.** Repeat steps 1-m for a total of two washes.
- o. Immediately add $500 \mu l$ Autofluorescence Solution prepared at step 5.1aiii. Pipette mix thoroughly before dispensing onto sample to prevent settling of reagent.
- **p.** Apply a new Xenium Cassette Lid on the Xenium Cassette, and incubate for **10 min** at **room temperature in the dark**.
- **q.** Prepare a thermal cycler with the following incubation protocol. Place a Thermocycler Adaptor on the thermal cycler and start the program.

Lid Temperature	Reaction Volume	Run Time
37°C (lid may be set to the lowest temperature if the instrument does not enable 37°C)	100 μΙ	-
Step	Temperature	Time hh:mm:ss
Pre-equilibrate	37°C	Hold
Drying	37°C	00:05:00

- r. Remove the Xenium Cassette Lid. Using a pipette, remove all Autofluorescence Solution. Discard old Cassette Lids.
- s. Add 1,000 μl 100% Ethanol. Wait 2 min.
- t. Remove all 100% Ethanol.
- u. Repeat steps s-t two more times.
- v. Place Xenium Cassette without lid on the Thermocycler Adaptor on the thermal cycler to dry. DO NOT close the thermal cycler lid.
- w. Skip pre-equilibrate step to initiate Drying.
- **x. Immediately** remove the Xenium Cassette from the Thermocycler Adaptor and place on a flat, clean work surface.
- y. Add 1,000 µl 1X PBS prepared at step 1.1 to rehydrate the tissue and incubate for 1 min at room temperature in the dark.
- z. Remove all 1X PBS.
- aa. Add 1,000 µl PBS-T and incubate for 2 min at room temperature in the dark.



Optional: photograph the slide against a white background. This image can be used for comparison purposes to identify tissue detachment downstream in the workflow. See Troubleshooting for more details.



ab. Store slides for ≤4 days at 4°C in the dark with a new Xenium Cassette Lid applied on the Xenium Cassette. Alternatively, proceed to next step. If storing slides, DO NOT discard the lid after use and instead save for Step 6, Nuclei Staining.



Xenium Cassettes must be covered with a Xenium Cassette Lid during storage.

5.2 Nuclei Staining

- a. Retrieve thawed Xenium Nuclei Staining Buffer prepared according to the Get Started table in step 5.0.
- **b.** Retrieve the Xenium Cassette from step 5.1ab and place on a flat, clean work surface.
- c. If stored, remove the Xenium Cassette Lid and using a pipette, remove all PBS-T from the well. **Save lid** for use in following indicated steps.
- d. Add 500 μl Xenium Nuclei Staining Buffer and incubate 1 min at room temperature in the dark.
- e. Remove all Nuclei Staining Buffer.
- **f.** Add **1,000 μl** PBS-T prepared at step 1.1 to the well.
- g. Incubate for 1 min at room temperature in the dark.
- **h.** Remove all PBS-T.
- **i.** Repeat steps f-h two more times.
- **i.** Add **1,000 ul** PBS-T.



k. Store slides for ≤4 days* at 4°C in the dark with a Xenium Cassette Lid applied on the Xenium Cassette. Alternatively, proceed to the Xenium Analyzer User Guide (CG000584). Use lid from step 5.2c if previously stored. Apply new lid if not stored.

*Alternatively, slides can be stored for **24 days** at **4°C in the dark** with a slide seal post-Nuclei Staining. PBS-T storage buffer should be exchanged with fresh PBS-T every 3-4 days.



Xenium Cassettes must be covered with a Xenium Cassette Lid or slide seal during storage.



Storing slides for more than 24 days may carry the risk for a lower number of genes or transcripts detected per cell, changes in tissue morphology over time, and microbe growth. These risks are dependent on many factors including input tissue quality and how nuclease-free/microbe-free the workflow is.


Troubleshooting



Tissue Detachment and Folding

Tissue detachment may result in a lack of decodable data in the region where detachment occurred. If the tissue has folded on itself, this may also cause elevated signal in the overlapping areas. Inspect images carefully to identify these areas. If tissue detachment is observed during this workflow, contact support@10xgenomics.com

Tissue Detachment in Human Breast as viewed on Xenium Analyzer Overview Scan



Percentages represent tissue detachment/"area that cannot be analyzed" at the end of the Xenium Analyzer workflow. White circles indicate areas of tissue detachment.

Tissue Detachment on Xenium Slides

Tissue sections may detach from Xenium slides during on-slide workflows. Tissue adhesion to the slide is impacted both by tissue processing and tissue quality.

Careful tissue preparation is critical for adhesion in on-slide workflows. Consult Xenium In Situ for Fresh Frozen - Tissue Preparation Guide (Document CG000578) and Xenium In Situ for FFPE - Tissue Preparation Guide (Document CG000579) for tissue QC and sectioning best practices. Below are some additional best practices for minimizing detachment during on-slide workflows:

- Do not pipette directly onto the tissue.
- Gently add and remove reagents from the well. Forceful addition or removal of reagents can agitate tissue and lead to detachment.
- Avoid touching in and around the Sample Area of the Xenium slide.
- Work quickly and carefully during reagent addition and removal.

In addition to following best practices, it is possible to monitor section adhesion on Xenium slides throughout the workflow. Taking a photograph of the slide at the beginning of the on-slide workflow and comparing with post-assay workflow images can help identify whether tissue shape has changed significantly, an indication of detachment. Steps when slide photos can be taken are noted in the protocol. These QC images can be compared with the DAPI overview scan as part of the Web Summary file to see whether tissue morphology has changed in the workflow.

If tissue detachment occurs, send pictures to support@10xgenomics.com for further assistance.

Bubbles during Workflow

Bubbles may occur throughout the Xenium In Situ Gene Expression workflow, including after Probe Hybridization and Ligation, and during PBS-T washes. Bubbles floating on the surface of the slide are unlikely to compromise assay performance. However, bubbles that are in contact with the tissue during a Xenium Analyzer run may result in a lack of decodable data in the tissue area where the bubbles occurred.

Avoid generating bubbles during reagent dispensing by pipetting slowly and avoiding expelling air from the pipette tip. Gently tap or rock the cassette after reagent dispension and inspect the cassette carefully to ensure liquid is fully covering the tissue. DO NOT aspirate or pop the bubbles as this could lead to tissue detachment or scratching of the tissue. Ensure there are no bubbles on the assembled cassette before loading it into the Xenium Analyzer.

Number of Washes

Post Hybridization and post-Ligation washes are critical for assay performance. Failure to perform the correct number of washes can reduce the fraction of usable decodable data. A similar effect is observed when washing for less than the recommended time, or when reagent is carried over during the washes. Remove all liquid from the well when washing, and refer to User Guide for correct number of washes and incubation times.

Samples Dry Out

Drying of tissue samples may lead to decreased decoding efficiency and unusable data. Work quickly and ensure reagents are dispensed evenly across tissues during incubation and wash steps throughout the workflow to prevent drying out of tissues. If processing two slides at a time, remove and add reagent from first slide before proceeding to second slide. Ensure tissue sections are always covered with reagent in between removal and addition steps. Note that there are no safe stopping points except for those described in the protocol and outlined specifically in the Protocol Steps & Timing.

Cassette Assembly Failure

Incorrect assembly of the Xenium Cassette with a Xenium slide can negatively impact assay performance. Always dry the front and the back of the slide completely using a lint-free laboratory wipe while avoiding touching or damaging of the tissue sections. Inspect the slide carefully to ensure it is seated fully within the cassette before assembly.

Example scenarios that may indicate incorrect Xenium Cassette assembly are described below:

- If a gap appears between the two halves of the cassette after assembly.
- If the cassette does not click shut or appears domed after assembly.
- If the Xenium Thermocycler Adaptor is wet following removal from the thermal cycler, indicating reagent leakage from the cassette.

If the cassette is incorrectly assembled, disassemble and reassemble the cassette as instructed in the following pages.

Incorrect cassette assembly as indicated by a gap between the two halves of the cassette

Correct cassette assembly

Incorrect cassette assembly







Inspect gasket during cassette assembly. Damaged gaskets can lead to leaks in the cassette.

Xenium Cassette Assembly



Exercise caution when handling slide edges to prevent injury.

Place top and bottom halves of cassette on bench



Press slide down into grooves of the bottom half of the cassette until it sits firmly in place



Apply even pressure on top of cassette until all clips click shut. Verify that clips are completely secured over tabs



Place Xenium slide with tissue side facing upwards into bottom half of cassette; ensure



Secure clips of top half with tabs of bottom half (on both sides)



Slides in images are representative.




Once cassette is assembled, DO NOT remove slide until after Xenium Analyzer imaging and decoding for optional H&E staining step.

Xenium Cassette Removal

Pull inner clips from inner tabs to detach top and bottom halves of cassette



Open cassette by continuing to lift inner clips



Hold slide by the label and lift slide out from bottom half



Slides in images are representative.

Xenium Cassette and Lid Cleaning

Xenium Cassettes and Lids are a single use item and are to be discarded after use. Cassettes (prior to cassette assembly) or lids that are accidentally dropped may be reused after thorough cleaning. Note that PBS-T washes DO NOT require sealing of the cassette.

Cleaning Procedure:

- Rinse the lid under running Milli-Q Water
- Spray with 70% isopropanol
- Rinse under running Milli-Q Water
- Spray with 70% isopropanol a second time
- Rinse under running Milli-Q Water
- Air dry

Incorrect Autofluorescence Quenching


Variation in stain color is normal and tissue-type dependent in tissue sections correctly stained with Autofluorescence Solution. Incorrect staining scenarios are listed below:

- Uneven staining with Autofluorescence Solution may be visible as a non-uniform stain across a tissue section.
- Overquenching can cause tissue to overheat on the Xenium Analyzer, and data generated in the overheated spots may be compromised or missing.

Improper Autofluorescence Quenching risks lower transcript quality scores and reduced median transcripts per cell but depends on the sample type.

Ensure Autofluorescence Solution is well mixed and dispensed uniformly across the tissue sections to avoid uneven staining. Autofluorescence Quenching has been optimized for a large number of tissue types and modifying the dilution listed in the User Guide is not recommended. Cassette should be sealed properly and firmly during incubation to prevent reagent evaporation.

Incorrect or insufficient Autofluorescence Quenching may cause variation in staining of tissues (Mouse Brain pictured below)




Incorrect Nuclei Staining


Incorrect staining of nuclei may lead to poor image quality and an inability to easily identify tissue or regions of interest when selecting areas to image during a Xenium Analyzer overview scan. Follow the Nuclei Staining protocol as instructed using the Xenium Nuclei Staining Buffer provided in the Xenium Slides & Sample Prep Reagents Kit - (2 slides, 2 rxns), PN-1000460. Confirm Xenium Nuclei Staining Buffer is well mixed and applied uniformly across tissue sections. All incubations with Xenium Nuclei Staining Buffer should be performed in the dark. If an alternate staining protocol or buffer is used, lower quality images may be obtained.

Incorrect or insufficient nuclei staining may impact image quality and region of interest selection (as viewed on a Xenium Analyzer overview scan)

Correct Nuclei Staining Protocol



No Nuclei Staining Performed



Appendix

Probe Panel Selection	68
Sample Shipping	68

Probe Panel Selection

Ensure that a compatible gene panel has been selected prior to executing the Xenium In Situ Gene Expression workflow. 10x Genomics provides the option of using pre-designed gene panels, pre-designed panels that are customized by adding genes of interest, and standalone custom gene panels.

Pre-designed Gene Panels

Refer to the 10x Genomics Support Website for the most current list of available panels and part numbers.

Custom Gene Panels

Contact your 10x Genomics Sales Executive for information about designing add-on custom gene panels that are compatible with pre-designed panels or standalone custom gene panels. If you do not know your Sales Executive, please contact customerservice@10xgenomics.com.

If utilizing a custom panel, the Design ID on the label of the tube containing the custom panel should match with the first portion of the custom gene panel electronic file name.

Sample Shipping

Processed Xenium slides may be shipped following the Xenium In Situ Gene Expression workflow. After Nuclei Staining, remove all PBS-T from last step, disassemble the Xenium Cassette, and place slides in a mailer filled to capacity with PBS-T. Ship the slide mailer containing processed Xenium slides in a container with ice packs. Place no more than two slides per mailer. Note that assay performance may be compromised post-shipping and handling.

Appendix 10xgenomics.com 68