

## z-Movi® Cell Avidity Analyzer

### Chip Cleaning Protocol (Version 6)

After finishing an experiment, ensure that chips are thoroughly cleaned according to protocol.

**Improper chip cleaning can impact monolayer adherence and the ability to run acoustic force ramps.** To achieve best cleaning results, start to clean a chip immediately after its final run.

1. Pull through 400  $\mu$ L dH<sub>2</sub>O. Repeat.
2. Pull through 400  $\mu$ L 10% bleach. Scrub back and forth several times. Repeat scrubbing with 400  $\mu$ L more bleach.
3. Pull in 400  $\mu$ L bleach and let incubate for 30 minutes at room temperature. Pull completely through once incubation time is finished.
4. Repeat step 3.

*Note: If there is cell debris after completely pulling through bleach, repeat step 4 until all debris is clear. You can also actively pull through new bleach and scrub in between incubations for more thorough cleaning.*

5. Pull through 400  $\mu$ L H<sub>2</sub>O five times.
6. Pull through 400  $\mu$ L 1M NaOH. Scrub back and forth several times. Repeat scrubbing with 400  $\mu$ L more NaOH.
7. Incubate with **1M NaOH for one hour at RT**. Screw the cap on the chip to prevent evaporation.

*\*\*\*Note: longer NaOH incubations can be detrimental to the chip, do not extend the NaOH incubation time\*\**

8. Pull NaOH completely through once incubation time is finished. The fluid should run evenly and smoothly (low glass surface hydrophobicity). This is a good indication of glass surface cleanliness and uniformity.
9. Pull through 400  $\mu$ L dH<sub>2</sub>O three times.
10. Pull through air multiple times until dry.

*Note: Make sure that the chip is completely dry before long-term storage. Use a larger (20 mL) syringe for faster drying. The chip is dry once you leave it in the dry incubator at 37C for 5 min, then pull air and do not see any solution flowing through.*

11. Chips are ready for the next avidity experiment. For long-term storage, place the clean, dry chip in a dry incubator at 37C.

Note: Improper chip cleaning (debris/ clogging) is not covered under warranty. Please take care to clean the chips after each experiment.

For help, please email your LUMICKS field application scientist for support.