

Protocol: Primary neuron culture for live-imaging of axonal cargoes

C. Alexander Boecker¹, Erika L.F. Holzbaur^{2,3}

¹Department of Neurology, University Medical Center Goettingen, 37077 Goettingen, Germany;

²Department of Physiology, Perelman School of Medicine, University of Pennsylvania, Philadelphia, PA 19104, USA; ³Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, USA

Key words

Dissection, mouse, primary neuron, live-imaging

Abstract

This protocol describes the preparation and culture of mouse primary cortical neurons for live-imaging experiments. Cortices were dissected from mouse embryos at day 15.5. Cortical neurons were isolated by digestion with 0.25% Trypsin and trituration with a serological pipette. Neurons were plated on glass-bottom imaging dishes in Attachment Media. After 5 hours in culture, Attachment Media was replaced with Maintenance Media, and AraC was added on the next day to prevent glia cell proliferation. Neurons were transfected 16-24 hours before imaging using Lipofectamine 2000.

Materials

Dissecting microscope

Dumont #5 Mirror Finish Forceps (Fine Science Tools, Cat# 1125123)

Micro spring scissors (Fine Science Tools)

Hemocytometer or automated cell counter (Countess 3, Thermo Fisher)

35 mm glass-bottom imaging dishes (MatTek, Cat# P35G-1.5-20-C)

15 mL conical tubes

10 cm cell culture dish

Reagents

Poly-L-lysine, mol wt 70,000 – 150,000 (Sigma, Cat# P1274-100MG)

10x HBSS (Thermo Fisher, Cat# 14185052)

1M HEPES (Thermo Fisher, Cat# 15630080)

2.5% Trypsin (Thermo Fisher, Cat# 15090046)

Minimum Essential Medium, MEM (Thermo Fisher, Cat# 11095072)

Horse Serum, heat inactivated (Thermo Fisher, Cat# 16050122)

D-(+)-Glucose Glucose solution, 45% (Sigma, G8769-100ml)

Sodium Pyruvate, 100 mM (Thermo Fisher, Cat# 11360070)

Trypan Blue Stain 0.4% (Thermo Fisher, Cat# T10282)

Neurobasal Medium (Thermo Fisher, Cat# 21103049)

B-27 Supplement (Thermo Fisher, Cat# 17504044)

GlutaMAX (Thermo Fisher, Cat# 35050061)

Penicillin-Streptomycin (Thermo Fisher, Cat# 15140122)

AraC (Sigma, Cat# C6645)

Lipofectamine 2000 Transfection Reagent (Thermo Fisher, Cat# 11668019)

Hibernate E Low Fluorescence Imaging Medium (BrainBits, Cat# HELF)

Safety warnings

Take necessary precautions with sharp objects during dissection. Follow institutional recommendations for disposal of animal tissue and biohazardous materials.

Day before dissection:

- Coat glass-bottom imaging dishes with PLL: hydrate 100 mg PLL (Sigma) in 50 mL 0.1 M borate buffer, pH 8.5. Store PLL stock solution (2 mg/mL) in 1 mL aliquots at -80°C. On the day before neuron dissection, dilute PLL in ddH₂O 1:20 to a final concentration of 100 µg/mL. Add 1 mL PLL to each glass-bottom imaging dish (MatTek) and incubate overnight at 37°C. Only coat the glass center with PLL.

Note: for easy handling, we find it helpful to place imaging dishes in 10 cm or 15 cm cell culture dishes.

- Prepare HBSS, attachment media and maintenance media: for 500 mL 1x HBSS, combine 50 mL 10x HBSS with 5 mL 1M HEPES, fill up to 500 mL with ddH₂O and filter-sterilize. Store 1x HBSS at 4°C and use within one month. For 50 mL attachment media, combine 5 mL heat-inactivated horse serum, 500 µL 100 mM sodium pyruvate, and 660 µL 45% Glucose; add MEM up to 50 mL. For 50 mL maintenance media, combine 500 µL GlutaMAX, 500 µL Penicillin/Streptomycin,

660 μ L 45% Glucose, 1 mL B-27, and add Neurobasal up to 50 mL. Store attachment media and maintenance media at 4°C.

Note: Maintenance Media should be used within 7 days. Attachment media can be kept at 4°C for 3-4 weeks.

Dissection of cortical neurons

- In the morning of the day of dissection: wash PLL-coated imaging dishes twice with sterile ddH₂O. Add 2 mL attachment media per imaging dish and leave dishes at 37°C in cell culture incubator. Warm required amount of attachment media and 1x HBSS (4.5 mL for one dissection) in 37°C water bath. Aliquot maintenance media into 10 cm cell culture dish to equilibrate in 37°C / 5% CO₂ cell culture incubator. Let 2.5% trypsin aliquots thaw at room temperature.
- Sacrifice pregnant mouse, dissect embryos, and place embryonic brains in HBSS on ice. Using a dissecting microscope, remove meninges from brain hemispheres with fine forceps. Isolate cortices using fine forceps and small spring scissors. Transfer dissected cortices into a 15 mL conical tube filled with 5 mL HBSS and keep on ice until all cortices are collected.

Note: Use clean and sterile equipment for all dissection steps to prevent bacterial contamination of neuron cultures.

Note: We find that using ice-cold HBSS helps preventing the tissue from getting sticky during the dissection. If HBSS gets too warm during the dissection, replace with fresh cold HBSS.

- Perform all following steps under a sterile tissue culture hood. Once all cortices are collected, remove HBSS from 15 mL conical tube and add 4.5 mL warm (37°C) HBSS and 0.5 mL 2.5% trypsin. After adding trypsin, invert the tube to mix. Then incubate for 10 minutes in a 37°C water bath.
- Remove HBSS-trypsin solution with a 5 mL serological pipette. Wash thrice with 7 mL attachment media: add attachment media, then wait until cortex tissue has settled at the bottom of the conical and remove attachment media with a serological pipette to repeat the washing step.

Note: we do not recommend using a vacuum aspirator for removing HBSS and attachment media, instead use a 10 mL serological pipette.

- Add 5 mL attachment media after the last washing step, then triturate cortices by pipetting up and down forcefully with a 5 mL serological pipette 10 – 15 times. Trituration is complete when no tissue clumps are visible and attachment media turns turbid. Let media with triturated tissue settle for 1-2 minutes, then transfer top 4.5 mL to a new tube to remove any remaining cell clumps.
- Mix 10 μ L cell suspension with 10 μ L 0.5% trypan blue in an Eppendorf tube, then count cells using a hemocytometer or an automated cell counter.
- Dilute cortical neurons to 1,000,000 cells/mL. For transfection and live-imaging, plate 200,000 cells per live-imaging dish. Place imaging dishes in 37°C cell culture incubator.

Note: Take up cells in a 200 μ L pipette, then plate drops of cells in different areas to distribute neurons evenly across the live-imaging dish.

- After 3-4 hours, use an aspirator to remove all attachment media and replace with 2 mL pre-equilibrated maintenance media per imaging dish. Cells should be attached to the glass-bottom dish at this point.

Note: maintenance media must always be pre-equilibrated to 5% CO₂ in 37°C incubator before adding to cells.

Neuronal cell culture

- On the day following the dissection, dilute AraC to 10 μ M in maintenance media and bring to 37°C. Add 200 μ L maintenance media + AraC to each imaging dish for a final AraC concentration of 1 μ M.
- Every 3-4 days, remove 600 μ L maintenance media from each dish and replace with 750 μ L fresh, pre-equilibrated maintenance media.

Note: cultured neurons are sensitive. Try to keep time outside the cell culture incubator to a minimum. If possible, use a separate incubator for primary neurons and keep openings to a minimum.

Transfection

- Transfect primary neurons on DIV6-7, ~ 16 hours before live-imaging.
- Replace conditioned media with fresh, pre-equilibrated maintenance media (2 mL per imaging dish). Save old media = conditioned media in a 10 cm cell culture dish at 37°C in the cell culture incubator.

- For each imaging dish, prepare two tubes with transfection reagents. In tube 1, add plasmid DNA to 200 μ L Neurobasal medium. In tube 2, add Lipofectamine 2000 to 100 μ L Neurobasal medium.

Note: it is important to use non-supplemented Neurobasal and not maintenance medium to set up the transfection reaction.

Note: the amount of Lipofectamine 2000 and plasmid DNA depends on the specific construct(s) used and may require optimization. We find that for transfection with one plasmid, 0.4 μ g DNA and 1 μ L Lipofectamine 2000 works well in most cases.

- Combine contents of tube 1 + 2 and mix by gently pipetting up and down 4-5 times. Incubate mix at room temperature for 10 minutes
- Add Lipofectamine-DNA mix to imaging dishes dropwise. Then incubate for 45 minutes at 37°C in cell culture incubator.
- Remove all transfection media and replace with conditioned media collected earlier. Return cells to incubator and image on the next day.

