



Plasmid Expansion

1. Thaw the competent cells on ice
2. Chill approximately $5\mu\text{g}$ ($2\mu\text{l}$) of Plasmid DNA in 1.5mL microcentrifuge tube
3. Add 50 μl of competent cells to the DNA. Mix gently by pipetting up and down or flicking the tube 4–5 times to mix the cells and DNA. Do not vortex.
4. Place the mixture on ice for 30 minutes. Do not mix.
5. Heat shock at 42°C for 30 seconds*. Do not mix.
6. Add 950 μl of room temperature media* to the tube.
7. Place tube at 37°C for 60 minutes. Shake vigorously (250 rpm) or rotate.
8. Warm selection plates to 37°C .
9. Spread 50–100 μl of the cells and ligation mixture onto the plates.
10. Incubate overnight at 37°C .
 - a. Please note: For the duration and temperature of the heat shock step as well as for the media to be used during the recovery period, please follow the recommendations provided by the competent cells’ manufacturer.

Plasmids and antibiotics

PLVX mCherry (HLA-A2)	Amp 50 $\mu\text{g}/\text{mL}$
Ploc MLANA (MART1)	Blasticidin (50 $\mu\text{g}/\text{mL}$)
psPAX2	Amp 50 $\mu\text{g}/\text{mL}$
VSV.G	Amp 50 $\mu\text{g}/\text{mL}$

Preparation of Liquid Media

LB-Ampicillin Broth and LB-Blasticidin Broth:

1. Add 20g of LB broth Lennox into 1 liter of distilled water in a Erlenmeyer flask
2. Autoclave the LB broth in the Flask on a wet cycle for 30min
3. Storage of LB in 4°C is recommended and is good for up to 1 year
4. Ensure that the LB broth has come to room temperature before you add ampicillin
5. The final concentration of ampicillin in the broth should be 50 $\mu\text{g}/\text{mL}$

6. LB broth containing ampicillin should be stored at 4°C for up to 1 month
-

Preparation of Solid Media

LB-Amp Agar and LB-Blasticidin Agar: Each petri Dish takes about 10mL of LB-Agar so scale the volume accordingly

1 LB agar tablet makes 50mL

1. Take 1 tablet and dissolve it in 50mL of distilled water in an autoclaved flask
 2. Autoclave the LB agar for 30min on a wet cycle
 3. Once the cycle is complete check to make sure that all the agar is dissolved
 4. Allow the LB agar to cool until it is comfortable to the touch (50°C), a water bath set to 50°C is useful for this step
 5. While the agar is cooling label the plates using sharpie for each antibiotics
 - a. Ampicillin = Red
 - b. Blasticidin = Blue
 6. The final concentration of ampicillin and Blasticidin should be 50µg/mL
 7. Make sure that the LB agar has cooled to 50°C as excessive heat will degrade the antibiotics
-