

Detailed operation manual of six requirements

Test environment: A horizontal worktable in a laboratory in accordance with JIS 8703-1983,

ISO 554-1976 (temperature 23° C, relative humidity 50%)

1) Average diameter of filter opening

1 Examine the average diameter of filter opening by referring to the catalog of the sheet

2 Observe under a polarized light microscope for sheets with a diameter of filter opening of

less than 10  $\mu$  m

3 Confirm opening size by using scanning electron microscopy (SEM), if possible

2) Thickness

1 Examine the thickness by referring to the catalog of the sheet

2 Prepare 10 test sheets (cut with scissors to 20 cm x 20 cm square )

3 Put 10 test sheets on a horizontal worktable

4 Lay a glass plate (200  $\pm$  2 mm square, a mass of 82  $\pm$  2 g, a thickness of 0.7 mm) on the

test sheets

5 After about 10 sec, measure the test pieces to 0.1 mm using a vertical ruler

## 20 3) Chemical resistance

21 1 Cut in half the sheet in order to compare each other

22 2 Immerse the whole part of one cut sheet in the test solution (xylene)

23 3 Cover the sheet with a lid and kept at room temperature for about 8 h

24 4 Wash with water, dried, and observe both samples

25 5 If no obvious change is observed, repeat the procedure 1-4 by changing the test solution

26 (formalin, methanol)

27

## 28 4) Water permeability test

## 29 Equipment / Consumables

- 30 • 1 x 50 mL plastic measuring pipette
- 31 • P1000 Pipette Tips (without filter) (Gilson) 1 pc
- 32 • P200 Pipette Tips (without filter) (Gilson) 2 pcs
- 33 • Container 1 (at least 8-10 cm in diameter and 15 cm deep, made of a material that
- 34 can be easily pierced) for example, 1 L-PET bottle.
- 35 • Container 2 (at least 15 cm in diameter and at least 5 cm in depth)
- 36 • Burette stand
- 37 • Palafilm
- 38 • stopwatch (measures up to 2 decimal places)
- 39 • funnel
- 40 • 500mL beaker
- 41 • Scissors and/or cutter
- 42 • plastic cutter
- 43 • drill

44

## 45 Solution

- 46 • 500mL distilled water

Test sheets (cut with scissors and/or cutter to 3cm x 3cm square)

# 1 Make-up for water permeability test device

## 1 Create a water storage tank


Cut off the suction opening of the 50mL disposable measuring pipette with a plastic cutter

## 2 Create drainage tank 1


Drill a hole to make an overflow port at a height of about 12 cm from the bottom of the container 1 (inner diameter of hole, about 10 mm)

## 3 Create a sample reservoir

- ① Cut the attachment part (a) of the tip of the P1000 pipette with scissors from the tip body (b) indicated by red arrow




- ② For one sheet, cut off at 2 cm from the tip of (b). In the case of two or more sheets, cut off 1.5 cm from the tip. (Red line)




- ③ Measure the inner diameter (cm) of the cut side of (b) and calculate the cross-sectional area (inner diameter x inner diameter x 3.14) (Cross-sectional area of test sheet).

- ④ Cut a 2 cm area from the opening of the tip of the P200 pipette with scissors (c)
- ⑤ Cut a 3 cm area from the opening of the other P200 pipette tip with scissors (d)
- ⑥ Connect (b) Connect (c) and (d) (yellow arrow)
- ⑦ Cover (b) with the test sheet (3cm×3cm) on the red line
- ⑧ Fix the sheet by fitting (a) (blue arrow)
- ⑨ Cut the sheet that protrudes from (a)
- ⑩ Seal and fix the connected parts with Parafilm



Complete view of the sample reservoir



## Assembly diagram of permeability test device

- 4 Vertical installation of a water storage tank on the workbench using a Burette stand
- 5 Assemble water storage tank, drainage tank 1, drainage tank 2, according to the reference above diagram
- 6 Adjust the height of the overflow opening to "5mL" on the scale of a 50mL disposable measuring pipette
- 7 Attach the sample reservoir to the end of the reservoir
- 8 Wrap the parafilm to prevent liquid leakage

## 2 Operation Procedure

- 1 Fill water storage tank with water by using a funnel until it flows out of the overflow port
- 2 Measure transit time ( $t_2 - t_1$  value) it takes for the water in a water storage tank to pass from the "40mL" scale to the "20mL" scale using a stopwatch
- 3 Repeat the operation in step 2 and confirm that the  $t_2 - t_1$  value was almost constant
- 4 Measure  $t_2 - t_1$  value 5 times

## 3 Calculation

The permeability coefficient was calculated by the following formula, applying JIS A

1218:2020.

$$k = 2.303 \frac{aL}{A(t_2 - t_1)} \log_{10} \frac{h_1}{h_2}$$

$\kappa$ :	permeability coefficient	cm/s		
a:	cross-section area of storage tank	cm <sup>2</sup>		
L:	thickness of test sheet	cm		
A:	cross-section area of test sheet	cm <sup>2</sup>		
$t_2 - t_1$ :	measurement time	s		
$h_1$ :	water level difference at time $t_1$	cm		
$h_2$ :	water level difference at time $t_2$	cm		

## 109 5) Water retention test

110 1 Prepare 3 sheets cut in 100 mm x 100 mm sized square.

111 2 Measure the mass of a cut sheet in the standard condition to 1 mg-level

112 3 Immerse 3 sheets in the water for at least 15 min

113 4 Remove 3 sheets from the water with tweezers

114 5 Allow the water to drip off for at least one min

115 6 Measure the mass of a cut sheet after wetting and dripping off the water (m2)

116 7 Calculate the water retention rate by the following formula

117 
$$m = (m_2 - m_1) / m_1 \times 100$$

118 m: water retention rate (%)

119 m1: mass of the specimen in the standard condition (mg)

120 m2: mass of the specimen after wetting and dripping off the water (mg)

121 8 Calculate the average value in accordance with JIS A 5209:2014, Rule B (rounding

122 method).

123

## 124 6) Cell transit test

## 125 Equipment / Consumables

126 ● P1000 Pipette Tips (without filter) (Gilson) 1 pc

127 ● P100 Pipette (Gilson) 2 pcs

128 ● P100 Pipette

- 129 ● P20 Pipette
- 130 ● P20 Pipette Tips
- 131 ● 2mL micro tube
- 132 ● 1.5mL/2mL microtube rack
- 133 ● microscope slide
- 134 ● alcohol-resistant marking pen
- 135 ● optical microscope (400x)
- 136 ● aspirator
- 137 ● 50mL plastic centrifuge tube
- 138 ● cell counterplate (Watson, Neubauer Improved)
- 139 ● Palafilm
- 140 ● Scissors

141

## 142 Solution/Reagent

- 143 ● 10% neutral buffered formalin solution
- 144 ● Residual specimen (ascites or pleural effusion) or cultured cells ( $1 \times 10^6$ - $10^7$  cells)
- 145 ● PBS 70mL
- 146 ● 0.4% Trypan blue solution
- 147 ● 95% Ethanol
- 148 ● Papanicolaou staining solution

149

150 1 Preparation of sample reservoir (described in line 58-80)

151 2 Preparation of cell suspension (collection of cells from ascites and/or pleural fluid)

152 1 Centrifugation residual specimens (ascites and/or pleural effusion) at 2000 rpm  
 153 for 5 min

154 2 Gently remove the supernatant with an aspirator

155 3 Collect more than 1mL of cell pellet

156 4 Add 30ml to 35ml of PBS and mix by inverting

157 5 Stand still for 15 min

158 6 Gently remove the supernatant with an aspirator

159 7 Repeat PBS wash (step 4-6)

160 8 Counts all cells in a cell suspension

161 ① Add 90  $\mu\ell$  of 0.4% trypan blue solution to 10  $\mu\ell$  of cell suspension (10-fold  
 162 dilution).

163 ② Gentle vortex for 2-3 sec

- ③ Inject 6  $\mu\ell$  samples with a P20 pipette through the sample inlet of the cell counter plate
- ④ Count the number of cells in 4 compartments using the cell counter plate.
- ⑤ Calculate the number of cells per 1  $\mu\ell$  of stock solution {(average number of cells of 4 compartments)  $\times$  10  $\times$  (dilution factor)}

Note : The count of one plot should be about 100.

If the number of cells is large, dilute the cell suspension 10 to 100 times with PBS.

If the number of cells is small, reduce the amount of 0.4% trypan blue solution added to the cell suspension 10  $\mu\ell$ .

- 7 Adjusted to  $1 \times 10^6$ -  $10^7$  cells/ml in 10% neutral buffered formalin solution

### 3 Operation Procedure

- 1 Apply 600  $\mu\ell$  of cell suspension ( $1 \times 10^6$ - $10^7$  cells/ml) to the sample reservoir (filter assembly) with a P1000 pipette
- 2 Allow to natural filtration
- 3 Measure the number of cells in the filtrate
  - ① Write an enclosure on a glass slide with an alcohol-resistant marking pen
  - ② 15  $\mu\ell$  drops of 95% ethanol in the enclosure, followed promptly by 15  $\mu\ell$  drops of filtrate
  - ③ Allow to natural drying
  - ④ Perform Papanicolaou staining
  - ⑤ Count cells in 10 fields of view at 400x under the microscopy