

## **PBMC Thawing Protocol:**

### **Reagent:**

PBMC washing medium; A: RPMI-1640 with 5 to 10% FBS, 100 U/ml penicillin, 100 U/ml streptomycin, 1 mM sodium pyruvate, 2 mM L-glutamine;

B: 50% X-vivo 15 medium (Lonza) + 25U/ benzonase,

C: 1X CTL-Anti-Aggregate-Wash™ (CTL))

PBMC complete culture medium: RPMI-1640 with 10% FBS, 100 U/ml penicillin, 100 U/ml streptomycin, 1 mM sodium pyruvate, 2 mM L-glutamine.

### **Protocol:**

1. Warm Washing buffer and medium to 37°C in a water bath.
2. Remove vials from liquid nitrogen and transport them to the lab on dry ice.
3. Thaw frozen vials, only 1 vial at a time, in a 37°C water bath. When cells are nearly completely thawed, carry the vials to the hood and swab them with 70% EtOH.
4. Gently remove PBMCs (avoid pipetting up and down, as the cells are very fragile at this stage) and transfer the cells into a 50 mL falcon tube (Fisher scientific #14-432-22) containing 25 mL warmed washing buffer.
5. Use 1mL washing medium to rinse out the cryovial and gently mix the cells by inverting the 15 mL Falcon tube ~5x.
6. Wash 1: Spin the cells: 400xg for 8 min at room temperature. Pour off the supernatant.
7. Wash 2:
  - 7.1 Suspend the cell pellet in 1 mL prewarmed medium (dropped slowly along the side of the tube) and resuspend the cell pellets, add 9 mL complete medium then Spin the cells 400xg for 8 min at room temperature.
  - 7.2 If cells were thawed in the presence of benzonase, perform an additional wash with culture medium in the absence of benzonase.
8. Count cells and determine viability by Trypan blue staining