Immunoblotting

Authors: Will Hancock-Cerutti^{1,2,3}, Pietro De Camilli^{1,3}

¹Departments of Neuroscience and of Cell Biology, Howard Hughes Medical Institute, Program in Cellular Neuroscience, Neurodegeneration and Repair, Yale University School of Medicine, New Haven, Connecticut 06510, USA

²Interdisciplinary Neuroscience Program and MD-PhD Program, Yale University School of Medicine, New Haven, Connecticut 06510, USA

³Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815

Abstract

This protocol describes collection of protein from cultured cells and immunoblotting, including immunoblotting of large proteins using a proprietary tris-acetate buffer system.

Keywords

Immunoblot, Western blot, tris-glycine, tris-acetate, siRNA, cGAMP, HT-DNA

Solutions to prepare

DMEM containing 10% FBS, 100 U/ml penicillin, 100 mg/ml streptomycin, and 2 mM L-glutamine (all from Gibco).

RIPA buffer: 150 mM NaCl, 10 mM Tris, 0.5 mM EDTA, 0.5% NP40 supplemented immediately before use with Protease Inhibitor Cocktail (Roche) and PhosStop phosphatase inhibitor (Roche)

TBS: 50 mM Tris-Cl, 150 mM NaCl, adjust pH to 7.5

TBST: TBS with 0.1% TWEEN-20 (Sigma-Aldrich)

3x Laemmli buffer: 188 mM Tris-HCl, 3% SDS, 30% glycerol, 0.01% bromophenol blue, 15% β-mercaptoethanol

Tris-glycine running buffer: 25 mM tris base, 192 mM glycine, 0.1% SDS in milliQ water

Tris-glycine transfer buffer: 25 mM Tris, 192 mM glycine, 20% methanol in milliQ water. Chill to 4 °C.

A. Cell culture and treatments

- 1. HeLa-M cells were cultured at 37°C in 5% CO₂ and DMEM containing 10% FBS, 100 U/ml penicillin, 100 mg/ml streptomycin, and 2 mM L-glutamine (all from Gibco).
- 2. For any given experiment, cells were plated at such density so as to be approximately 90% confluent at the time of lysis.
- 3. For experiments using siRNA, 60 pmols of the indicated siRNA was transfected using 6μL Lipofectamine RNAiMax (ThermoFisher) in Opti-MEM (Gibco) per well according to manufacturer protocol. Cells were lysed 72 hours after siRNA transfection.
- For experiments using cGAMP, 8 μg/L of cGAMP was transfected using using 6μL Lipofectamine RNAiMax (ThermoFisher) in Opti-MEM (Gibco) per well according to manufacturer protocol.
- 5. For experiments using herring testes (HT)-DNA, 1 ug HT-DNA was transfected using using 3 µL Lipofectamine 2000 (ThermoFisher) in Opti-MEM (Gibco) per well according to manufacturer protocol.

B. Cell lysis and sample preparation

- 1. Supplement RIPA buffer with Protease Inhibitor Cocktail (Roche) and PhosStop phosphatase inhibitor (Roche) and chill on ice.
- 2. Aspirate media from cells and rinse cells with PBS on ice. Aspirate PBS thoroughly.
- 3. Pipette RIPA lysis buffer onto cells and scrape cells using a cell lifter (Corning).
- 4. Pipette lysis buffer containing cell mass into Eppendorf tube.
- 5. Incubate Eppendorf tube on ice for 30 minutes.
- 6. Every 10 minutes, pipette lysis mixture up and down 10 times with a P-200 pipette tip (a total of 3 cycles).

NOTE: Take care not to introduce bubbles.

7. Centrifuge at 15,000 x g for 10 min at 4° C and collect the post-nuclear supernatant in a new Eppendorf tube.

NOTE: Samples can be snap frozen in liquid nitrogen at this step and stored at -70°C.

- 8. Determine protein concentration in sample using Pierce BCA assay (ThermoFisher)
- 9. Prepare samples at desired concentration and add 3x Laemmli buffer.

C. Gel electrophoresis and immunoblotting (Tris-glycine buffer system)

- 1. Incubate samples at 95°C for 5 minutes.
- 2. During this incubation, prepare gel apparatus with Mini PROTEAN TGX 4-20% trisglycine gels (Bio-Rad) and running buffer.
- 3. Load samples into gel and run until dye front reaches bottom (90-120 V).
- 4. Remove gel and set up transfer cassette with nitrocellulose membrane.
- 5. Transfer at 120 V for 1 hour at 4°C in tris-glycine transfer buffer.
- 6. Remove nitrocellulose membrane and stain for total protein with ponceau stain.
- 7. Wash with milliQ water.
- 8. Block membrane with 5% milk in TBST for 1 hour at 22°C.
- 9. Incubate membrane with primary antibodies in 2.5% milk in TBST overnight at $4^{\circ}C$.

NOTE: Optimal primary antibody incubation time and temperature can be determined empirically for a given primary antibody

- 10. Wash membrane for 5 minutes with TBST. Repeat a total of 3 times.
- 11. Incubate membrane with secondary antibodies conjugated to IRdye 800CW or IRdye 680CW (1:10,000, Licor) in 2.5% milk in TBST for 1 hour at 22°C
- 12. Wash membrane for 5 minutes with TBST. Repeat a total of 3 times.
- 13. Wash membrane for 5 minutes with TBS. Repeat a total of 3 times.
- 14. Image membranes using a Licor Odyssey Infrared Imager.

D. Gel electrophoresis and immunoblotting (Tris-acetate buffer system)

- 1. For VPS13C immunoblotting, samples were lysed and post-nuclear supernatant was collected as above.
- 2. Mix post-nuclear supernatant with NuPAGE LDS Sample Buffer and Reducing Agent (Thermofisher) and incubated for 10 minutes at 70 °C.
- 3. During this incubation, prepare gel apparatus with NuPage Tris-Acetate 3-8% gels and NuPage Running Buffer (Thermofisher)
- 4. Remove gel and set up transfer cassette with nitrocellulose membrane.
- 5. Transfer at 0.05 mA for 16 hours at 4°C in NuPage transfer buffer (Thermofisher).
- 6. Remove nitrocellulose membrane and stain for total protein with ponceau stain.
- 7. Wash with milliQ water.
- 8. Block membrane with 5% milk in TBST for 1 hour at 22°C.
- 9. Incubate membrane with primary antibodies in 2.5% milk in TBST for 2 hours at 22°C.

NOTE: Optimal primary antibody incubation time and temperature can be determined empirically for a given primary antibody

- 10. Wash membrane for 5 minutes with TBST. Repeat a total of 3 times.
- 11. Incubate membrane with secondary antibodies conjugated to IRdye 800CW or IRdye 680CW (1:10,000, Licor) in 2.5% milk in TBST for 1 hour at 22°C.
- 12. Wash membrane for 5 minutes with TBST. Repeat a total of 3 times.
- 13. Wash membrane for 5 minutes with TBS. Repeat a total of 3 times.
- 14. Image membranes using a Licor Odyssey Infrared Imager.