

Protocol for hippocampal neuronal cultures

Authors: Andrés Guillén-Samander^{1,2} and Pietro De Camilli^{1,2}

¹Departments of Neuroscience and of Cell Biology, Howard Hughes Medical Institute, Program in Cellular Neuroscience, Neurodegeneration and Repair, Yale University School of Medicine, New Haven, Connecticut 06510, USA

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815

Solutions to prepare

Poly-D-Lysine (PDL)

- Dilute PDL to 0.1mg/ml in Borate Buffer 0.1M, pH 8.5 and filter sterilize

Plating Medium (filter sterilize and keep stable for 3 weeks at 4°C)

- Neurobasal (Gibco): 183ml
- FBS (Gibco): 10ml (to final concentration of 5%)
- Glutamax (Gibco): 2ml
- B27 (Gibco): 4ml
- Penicilin (Gibco): 0.5ml (to final concentration of 50U/ml)
- Streptomycin (Gibco): 0.5ml (to final concentration of 50mg/ml)

Neuronal HBSS (nHBSS; filter sterilize and keep stable for 3 weeks at 4°C)

- HBSS (Gibco): 490ml
- 1M HEPES (Gibco): 6ml
- 100mM Pyruvic Acid (Gibco): 6ml
- H2O: 100ml
- Penicilin (Gibco): 1.5ml (to final concentration of 50U/ml)
- Streptomycin (Gibco): 1.5ml (to final concentration of 50mg/ml)

Neuronal Medium (filter sterilize and keep stable for 3 weeks at 4°C)

- Neurobasal (Gibco): 193ml
- Glutamax (Gibco): 2ml
- B27 (Gibco): 4ml
- Penicilin (Gibco): 1.5ml (to final concentration of 50U/ml)
- Streptomycin (Gibco): 1.5ml (to final concentration of 50mg/ml)

Papain solution (Prepare fresh and adjust pH if needed with 0.1M NaOH, incubate at 37°C for 30 minutes to dissolve papain and filter)

- Papain (Worthington): 200ul (to final concentration of 20U/ml)
- nHBSS: 10ml
- L-cysteine (Sigma): 2mg (to 0.2mg/ml)

Protocol

1. Coat MatTek dishes with 1ml per dish of 0.1mg/ml Poly-D-Lysine (Sigma) for at least 1.5hrs to overnight at 37°C.
2. Wash dishes twice with culture grade water and let dry.
3. Prepare papain solution and leave at 37°C for 30 minutes.
4. Dissect hippocampi from at least 3 P0 mouse brains using a stereo microscope. Collect tissue in ice cold Neuronal HBSS (nHBSS).
5. Transfer isolated hippocampi into a fresh cold nHBSS containing dish.
6. Cut tissue into $\approx 1\text{mm}^3$ pieces and transfer into a 15ml Falcon tube with 10ml ice cold nHBSS and let sediment on ice.

Note: From this step on everything is done under a sterile hood.

7. Aspirate medium and wash 2-3 times with 10ml fresh ice cold nHBSS.
8. Add DNase to papain solution and filter sterilize. Incubate tissue prep with papain solution for 30 minutes at 37°C on a rocking platform.
9. Aspirate the enzyme solution and wash twice with plating medium and then twice with nHBSS.
10. Allow debris to settle for several minutes and collect supernatant.

Warning: Tissue will be softer after papain solution, handle with care to avoid dissociating the cells.

11. Resuspend samples in 2ml cold nHBSS. Gently dissociate neurons with a P1000 filter tip by pipetting up and down for 10-12 times.

Warning: Avoid generating any bubbles.

12. Count the cells.
13. For imaging, seed 75,000 neuronal cells as a drop (usually around 100ul) in the PDL-coated coverslip of MatTek dishes in Plating Medium.
14. After 3hr to overnight incubation at 37°C and 5% CO₂, change the plating medium to neuronal medium.
15. Remove 500ul of media and add 1ml of fresh neuronal media every 3-4 days.

Suggestion: Incubate the media at 37°C and 5% CO₂ before adding it to neurons.