

## Isolation of Nucleated Cells from Whole Blood

Steven B. Wells<sup>1</sup>, Peter A. Szabo<sup>2</sup>, Nora Lam<sup>2,3</sup>

<sup>1</sup>Department of Systems Biology, Columbia University Irving Medical Center, New York, NY 10032, USA

<sup>2</sup>Department of Microbiology and Immunology, Columbia University Irving Medical Center, New York, NY 10032, USA

<sup>3</sup>Department of Pathology and Cell Biology, Columbia University Irving Medical Center, New York, NY 10032

### Keywords

blood, CD45, lymphocytes, myeloid, isolation, density gradient, ficoll, immune, 10x, scRNAseq, flow cytometry, WBC, leukocyte, single cell suspension, T cell

### 1. Abstract

This protocol describes a method for the isolation of pan-lymphocytes and pan-myeloid cells from human whole blood. By providing defined media formulations, volumes at each step, and a defined dilution factor for density centrifugation, it yields consistent single-cell suspensions across samples.

### 2. Materials

Dulbecco's Phosphate Buffered Saline (DPBS) (Fisher Scientific, Cat. No.: 14-190-144)

Penicillin-Streptomycin-Glutamine (100X) (Fisher Scientific, Cat. No.: 10-378-016)

50mL Centrifuge Tube (Fisher Scientific, Cat. No.: 12-565-271)

Iscove's Modified Dulbecco's Medium (IMDM) (Fisher Scientific, Cat. No.: 12-440-053)

Fetal Bovine Serum (FBS) (Fisher Scientific, Cat. No.: 10-099-14)

EDTA 0.5M pH 8.0 (Fisher Scientific, Cat. No.: 15-575-020)

Ficoll-Paque™ PLUS Media (Fisher Scientific, Cat. No.: 45-001-749)

Mr. Frosty (Fisher Scientific, Cat. No.: 51000001)

Cryostor CS10 (Fisher Scientific, Cat. No.: NC9930384)

Cryogenic Vials (Fisher Scientific, Cat. No.: 09-761-71)

5mL Falcon™ Round-Bottom Polypropylene Tubes (Fisher Scientific, Cat. No.: 14-959-11A)

Solution 13 AO/DAPI (Chemometec, Cat. No.: 910-3013)

NC-Slide A8 (Chemometec, Cat. No.: 942-0003)

Falcon™ Plastic Disposable Transfer Pipets (Fisher Scientific, Cat. No.: 13-680-50)

### 3. Equipment

Centrifuge

Cell Counter - NC-3000

### 4. Protocol

#### 4.1. Preparing Mediums and Buffers

4.1.1. Create the following **IMDM-FBS-PSQ Media** in a 500mL bottle of IMDM by using the table below:

<i>Component</i>	<i>Volume (mL)</i>	<i>Starting Conc.</i>	<i>Final Conc.*</i>
IMDM	500		-
Penicillin-Streptomycin-Glutamine	5	100X	1X
FBS	50	100%	10%

*Table 1.*

*\*Final Concentration is approximate.*

4.1.2. Create the following **DPBS-FBS-EDTA Solution** in a bottle of DPBS by using the table below:

<i>Component</i>	<i>Volume (mL)</i>	<i>Starting Conc.</i>	<i>Final Conc.*</i>
DPBS	500	-	-
FBS	25	100%	5%
EDTA	1	0.5M	1mM

*Table 2.*

*\*Final Concentration is approximate.*

#### 4.2. Preparation of Blood

4.2.1. Record the total volume of whole blood to be processed.

\_\_\_\_\_ mL

4.2.2. Spin the whole blood 400 x g for 10 minutes in the anti-coagulant tubes, remove the plasma layer, and distribute to cryovials – up to 2mL/vial.

Record the total volume of plasma: \_\_\_\_\_ mL and the number of vials: \_\_\_\_\_.

4.2.3. Replace the plasma volume removed from the whole blood with DPBS-FBS-EDTA Solution.

4.2.4. Divide the whole blood into 10mL aliquots and distribute to separate 50mL tubes.

4.2.5. Dilute the whole blood using 4 volumes or 40mL DPBS-FBS-EDTA Solution; invert to mix.

**NOTE:** This is the optimum dilution to maximize cell recovery.

#### 4.3. Ficoll-Paque

4.3.1. Layer the blood/DPBS-FBS-EDTA Solution mixture from the 50mL tubes 25mL at a time in separate 50mL tubes on top of 15mL of Ficoll-Paque Media PLUS.

**NOTE:** For any remaining volume, add DPBS-FBS-EDTA Solution to bring the volume to 25mL, and layer as described in this step.

- 4.3.2. Spin for 20 minutes, 1200 x g at 20°C with 4 acceleration and 0 brake, evenly distribute the tubes across the entire rotor to prevent wobbling (use all four buckets if possible as opposed to just two).
- 4.3.3. Remove the mononuclear cell layer from each tube with a transfer pipet to 50mL tubes - mononuclear layers may be combined at this step to reduce the number of tubes to spin. Add cold DPBS-FBS-EDTA Solution to a final volume of 50mL and centrifuge the cell suspensions for 10 minutes at 400 x g, 4°C.
- 4.3.4. Remove the supernatant and re-suspend the cell pellet in 50mL cold DPBS-FBS-EDTA Solution and centrifuge the cell suspension for 10 minutes at 120 x g, 4°C.
- 4.3.5. Remove the supernatant and re-suspend the cell pellet in cold 10mL IMDM-FBS-PSQ Media.

#### 4.4. Cell Count

- 4.4.1. Count cells, and viability by using the NC-3000 cell counter. Calculate total viable cells and record below:

cell number: \_\_\_\_\_ cells/mL, \_\_\_\_\_ % viable

final volume: \_\_\_\_\_ mL

$$cell\ number\left(\frac{cells}{mL}\right) * viability(\%) * final\ volume(mL) = total\ viable\ cells$$

Total Viable Cells: \_\_\_\_\_

#### 4.5. Freeze-down and QC

- 4.5.1. **(Optional QC)** Aliquot  $2 \times 10^6$  cells to a 5mL Falcon tube and place on ice for subsequent flow cytometric analysis.
- 4.5.2. Aliquot cells for analysis or experimentation, and then freeze down remaining cells in up to  $2 \times 10^7$  aliquots using Cryostor CS10 Medium, a Mr. Frosty, and a -80°C freezer (1-1.5mL aliquots, round down to the nearest 20 million cells and discard/freeze/use any left over cells). Record the number of vials frozen: \_\_\_\_\_.