

CVNSP3A-c011 small scale purification

PAGE24-00789

Author: **Wang, Korvus**

Date Started: **2024-May-01**

Experiment Started:

Projects: **Purification;ASAP**

Related Pages: **PAGE23-00304;PAGE24-00534**

Referenced by:

Tags:

Expression and information

CVNSP3A-c011

Expression done by Ellie

CVNSP3A-e004

CVNSP3A-p004

Final pellet ~60g

MHHHHHHSSGVDLGTEENLYFQ/////SMVNSFSGYLKLTDNVYIKNADIVEEAKVKPTVVVNAANVYLKHGGGVAGALNKATNNA
MQVESDDYIATNGPLKVGGS CVLSGHNLA KHCLHVVGP NVNKGEDIQLLSAYENFNQHEVLLAPLLSAGIFGADPIHSLRVCVDT
VRTNVYLAVFDKNLYDKLVSSFL

Post-cleavage:

Molecular weight: 18157.79

Ext. coefficient 10430

Abs 0.1% (=1 g/l) 0.574, assuming all Cys residues are reduced

Purification

Buffers

Lysis buffer - 50 mM HEPES pH 7.5, 500 mM NaCl, 5% glycerol, 0.5mM TCEP

Wash buffer 1 - 50 mM HEPES pH 7.5, 500 mM NaCl, 5% glycerol, 20mM imidazole, 0.5mM TCEP

Elution Buffer (IMAC) – 50 mM HEPES pH 7.5, 500 mM NaCl, 5% glycerol, 500 mM imidazole, 50mM ArgGlu, 0.5mM TCEP

Gel Filtration Buffer (SEC) – 50 mM HEPES pH 7.5, 500 mM NaCl, 5% glycerol, 0.5mM TCEP

Cell Lysis

1. Add lysis buffer to pellet until total volume is 400mL. Supplement with 1:4000 dilution of benzonase, 0.5mg/mL lysozyme, 1mM final MgSO₄, and 1 PIC tablet.
2. Mix lysis mixture until homogenous.
3. Sonicated on ice at 40% amplitude for 10 minutes (2 seconds on 4 seconds off).
4. Clarified lysate by centrifugation at 18,000rpm, 4°C for 1 hour.

IMAC

1. Wash and equilibrate 2mL bed volume of Ni sepharose resin on gravity flow column, first with distilled water, then with wash buffer. (use cut pipette tip)
2. Add equilibrated resin to clarified lysate. Incubate in cold room on rotation for 1hr.
3. Put stopper on column. Re-suspend resin with small volumes of the clarified supernatant, and pour into 500mL v-bottom centrifuge bottles.

4. Incubate resin and supernatant on rotator for ~30min. Pour mixture back on gravity flow column and rinse remaining resin from bottle with ~30mL wash buffer. Load all onto column and collect FT.
5. Wash resin 2x with 50mL wash buffer each time. Collect both washes as wash 1 and 2.
6. Elute with 7.5mL elution buffer, 10mL incubation. Collect both elutions as elution 1 and 2.

Desalting and Cleavage

Imidazole was removed from eluate using HiPrep Desalting column into gel filtration buffer. His-2Strep-TEV was added in 1:20 molar ratio and incubated in cold room overnight.


rIMAC

Imidazole was added to final concentration of 10mM to the cleavage mix. Cleavage mix passed over washed IMAC resin 3 times. FT collected.

Washed with 10mL wash buffer.

Elution with 10mL elution buffer to see what stuck to the resin.

IMAC and cleavage results




Gel filtration

FT and wash was concentrated in Vivaspin 15 10kDa MWCO concentrators to ~1mL
ACCIDENTALLY SPILLED HALF


Solid aggregates pelleted by centrifugation at 27,000xg, 10 mins, 4oC

Sample injected onto superdex 75 16/600 on AKTA pure, 1mL/min

s75 chromatogram



SEC results




Final sample

Fraction F9-G8 were pooled and concentrated in Vivaspin 15 10kDa MWCO concentrators.

Final concentratrion: 58.17 mg/mL, 8x50uL

Total yield: 46.5mg

Title missing - double click to edit



MS

