





	Catalog No.	Size	
 Package contents	11766050	50 reaction (with ezDNase enzyme)	 Kit contents
	11766500	500 reaction (with ezDNase enzyme)	
	11756050	50 reaction	
	11756500	500 reaction	


 **Storage conditions** Store all contents at -20°C (non-frost-free)

 **Required materials**


- Template: RNA

 **Timing**

- Preparation time: 5 minutes
- Total incubation time: 25–27 minutes

 **Product description**

- The Invitrogen™ SuperScript™ IV VILO™ Master Mix with ezDNase enzyme is an optimized solution for first strand cDNA synthesis in two-step RT-PCR applications.
- The 5X master mix includes SuperScript™ IV Reverse Transcriptase, a proprietary recombinant RNase inhibitor, helper proteins, stabilizer proteins, oligo (dT)₁₈, random hexamer primers, MgCl₂, and dNTPs.
- The SuperScript™ IV VILO™ No RT Control contains all the components of the SuperScript™ IV VILO™ Master Mix except the reverse transcriptase enzyme. The No RT Control is used in a RT minus reaction to verify the absence of genomic DNA (gDNA) contamination in the RNA sample.
- ezDNase enzyme (Cat. No. 11766051) is a novel double-strand specific thermolabile DNase that is used to remove gDNA contamination from template RNA prior to the RT reaction. The enzyme is available as part of the kit or as a standalone product.

 **Online resources** Visit our **product page** for additional information and protocols. For support, visit thermofisher.com/support.

Guidelines for RNA preparation

- Use high-quality, intact RNA for accurate quantification. RNA must be devoid of RNase contamination and handled using aseptic conditions.
- Isolate total RNA with **TRIzol™ Reagent**, the **PureLink™ RNA Mini Kit**, or the **MagMAX™-96 Total RNA Isolation Kit**.
- Determine RNA quality using a bioanalyzer or by agarose gel electrophoresis.

Guidelines for reverse transcription

- Use up to 2.5 µg of total RNA as starting material in a 20-µL reaction.
- For GC-rich or structurally complex RNA templates, increasing the RT incubation temperature up to 65°C may improve cDNA synthesis results.
- To verify the absence of gDNA contamination in the template RNA, perform a no RT control reaction by substituting SuperScript™ IV VILO™ Master Mix with the SuperScript™ IV VILO™ No RT Control.

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.thermofisher.com/us/en/home/global/terms-and-conditions.html.

If you have any questions, please contact Life Technologies at www.thermofisher.com/support.

Important licensing information









These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

Disclaimer: TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.







Corporate entity: Life Technologies | Carlsbad, CA 92008 USA | Toll Free in USA 1.800.955.6288

©2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

Reverse transcription protocol for SuperScript™ IV VIL0™ Master Mix with ezDNase enzyme

Step	Action	Procedure details												
1 	Prepare gDNA digestion reaction mix (on ice)	<p>For each RT reaction or No RT Control reaction, prepare a 10 µL gDNA digestion reaction mix in a RNase-free tube on ice with the following components:</p> <table border="1"> <thead> <tr> <th>Component</th> <th>Volume</th> </tr> </thead> <tbody> <tr> <td>10X ezDNase Buffer</td> <td>1 µL</td> </tr> <tr> <td>ezDNase enzyme</td> <td>1 µL</td> </tr> <tr> <td>Template RNA (1 pg to 2.5 µg total RNA)</td> <td>varies</td> </tr> <tr> <td>Nuclease-free Water</td> <td>to 10 µL</td> </tr> </tbody> </table>	Component	Volume	10X ezDNase Buffer	1 µL	ezDNase enzyme	1 µL	Template RNA (1 pg to 2.5 µg total RNA)	varies	Nuclease-free Water	to 10 µL		
Component	Volume													
10X ezDNase Buffer	1 µL													
ezDNase enzyme	1 µL													
Template RNA (1 pg to 2.5 µg total RNA)	varies													
Nuclease-free Water	to 10 µL													
2 	Digest gDNA	Gently mix and incubate at 37°C for 2 minutes. Briefly centrifuge the reaction and place on ice.												
3 	Prepare RT and No RT Control reaction mixes (on ice)	<p>Add the following components to the tube containing the 10-µL reaction mix from step 2 on ice.</p> <table border="1"> <thead> <tr> <th>Component</th> <th>RT reaction</th> <th>No RT Control reaction</th> </tr> </thead> <tbody> <tr> <td>SuperScript™ IV VIL0™ Master Mix</td> <td>4 µL</td> <td>—</td> </tr> <tr> <td>SuperScript™ IV VIL0™ No RT Control</td> <td>—</td> <td>4 µL</td> </tr> <tr> <td>Nuclease-free Water</td> <td>6 µL</td> <td>6 µL</td> </tr> </tbody> </table>	Component	RT reaction	No RT Control reaction	SuperScript™ IV VIL0™ Master Mix	4 µL	—	SuperScript™ IV VIL0™ No RT Control	—	4 µL	Nuclease-free Water	6 µL	6 µL
Component	RT reaction	No RT Control reaction												
SuperScript™ IV VIL0™ Master Mix	4 µL	—												
SuperScript™ IV VIL0™ No RT Control	—	4 µL												
Nuclease-free Water	6 µL	6 µL												
4 	Anneal primers	Gently mix and incubate at 25°C for 10 minutes.												
5 	Reverse transcribe RNA	Incubate at 50°C for 10 minutes.												
6 	Inactivate enzyme	Incubate at 85°C for 5 minutes.												
7 	qPCR amplification	<p>Use the diluted or undiluted cDNA for qPCR or store at -20°C for up to one week, or -70°C for long term storage.</p> <p> Click to see Guidelines for optimizing qPCR amplification.</p>												

Reverse transcription protocol for SuperScript™ IV VILO™ Master Mix (without ezDNase enzyme treatment)

Step	Action	Procedure details															
1 	Prepare RT and No RT Control reaction mixes (on ice)	Add the following components to an empty RNase-free tube on ice. <table border="1" data-bbox="758 207 2003 427"> <thead> <tr> <th>Component</th> <th>RT reaction</th> <th>No RT Control reaction</th> </tr> </thead> <tbody> <tr> <td>SuperScript™ IV VILO™ Master Mix</td> <td>4 µL</td> <td>—</td> </tr> <tr> <td>SuperScript™ IV VILO™ No RT Control</td> <td>—</td> <td>4 µL</td> </tr> <tr> <td>Template RNA (1 pg to 2.5 µg total RNA)</td> <td>varies</td> <td>varies</td> </tr> <tr> <td>Nuclease-free Water</td> <td>to 20 µL</td> <td>to 20 µL</td> </tr> </tbody> </table>	Component	RT reaction	No RT Control reaction	SuperScript™ IV VILO™ Master Mix	4 µL	—	SuperScript™ IV VILO™ No RT Control	—	4 µL	Template RNA (1 pg to 2.5 µg total RNA)	varies	varies	Nuclease-free Water	to 20 µL	to 20 µL
Component	RT reaction	No RT Control reaction															
SuperScript™ IV VILO™ Master Mix	4 µL	—															
SuperScript™ IV VILO™ No RT Control	—	4 µL															
Template RNA (1 pg to 2.5 µg total RNA)	varies	varies															
Nuclease-free Water	to 20 µL	to 20 µL															
2 	Anneal primers	Gently mix and incubate at 25°C for 10 minutes.															
3 	Reverse transcribe RNA	Incubate at 50°C for 10 minutes.															
4 	Inactivate enzyme	Incubate at 85°C for 5 minutes.															
5 	qPCR amplification	Use the diluted or undiluted cDNA for qPCR or store at –20°C for up to one week, or –70°C for long term storage.  Click to see Guidelines for optimizing qPCR amplification.															