

Passing Placental Organoids generated from full-term tissue

Reagents, Solutions and Materials prepared in advance:

- a) **Pre-cool** blunt 200 µl pipette tips.
- b) Pre-warm multi-well TC plate (this protocol uses 24-well TC plate, cat# 3526, Costar) and Stem Pro Accutase (Gibco, Cat # A11105-01) supplemented with 10 μM Y-27632 (Sigma, Y0503-1MG; 100 × dilution from stock solution, this is Rock inhibitor used to prevent the stem cells from anoikis during the passaging process).
- c) **Pre-thaw** Matrigel on the ice for at least 2 hrs, we usually thaw it o/n on the ice.
- d) **Prepare** 20% (vol/vol) FBS medium and basal media as needed.
- 1. Remove old complete growth medium (TOM).
- 2. Add 500 ul of fresh basal medium to each well.
- 3. Pre-coat a wide orifice 1 ml pipette tip (Finntip 1000, Thermo Fisher, 9405160) using FBScontaining media and use this tip to gently scrape off the Matrigel domes (including organoids) without touching the bottom of the well.
- 4. Using the same wide orifice tip, carefully transfer the released mixture of Matrigel and organoids into a 15 ml conical centrifuge tube and briefly pipette several times.
- 5. Centrifuge at 600 RPM, RT for FIRST 6 minutes.
- Carefully remove supernatant as much as possible using a 1 ml pipette and then remove remaining media using a 200uL tip if necessary; do not use glass Pasteur pipette or vacuum to aspirate.
- 7. Add 1 ml of pre-warmed dissociation reagent
 - a. StemPro Accutase (Life Technologies, A11105-01) or TrypLE express Life Technologies, 12605-028). Pre-warm prior to use).
 - b. Add ROCK inhibitor to dissociation agent; for 1 ml of dissociation agent, add 10 ul of inhibitor.
- 8. Incubate in 37 °C water bath for 6-10 min and swirl from time to time during this incubation.
- 9. Centrifuge at 600 RPM, RT for SECOND 6 min.
 - a. Do additional centrifugation at 600 RPM (do not increase RPM) if Matrigel-organoids mixtures are still floating in supernatant.
- 10. Remove supernatant as much as possible by using 1 ml pipette and then remove remaining media using a 200uL tip if necessary; do not use glass Pasteur pipette to aspirate.
- 11. Add 200 ul of basal medium.
- 12. Use autopipette to disturb/resuspend pellet (pipetting times depend on organoid size)
 - a. Set autopipette to full force (level 8)
 - b. Pipette 200x
 - c. Check suspension
 - d. Pipette 50x
 - e. Check suspension
 - f. Pipette 30x if not evenly disrupted
- 13. Remove but do not discard autopipette tip.
- 14. Flush inside of autopipette tip using 1 ml of basal medium into above solution.
- 15. Centrifuge at 600 RPM, RT for 6 minutes.
- 16. Remove supernatant using 1 ml pipette; do not use glass Pasteur pipette to aspirate, then put the 15 ml conical tube with pellet into ice.
- 17. Resuspend the pellet with pre-thawed Matrigel using pre-cooled blunt 200 μl pipette tips (Fisher 02-707-134), the amount of Matrigel: 40 X number of wells desired.



- a. For example, for 6 wells, add 240-250 ul of Matrigel.
- b. Matrigel should be kept on ice.
- c. Do not discard pipette tip; transfer to 20-200ul pipette.
- 18. Carefully dispense 40 ul aliquot of Matrigel-organoid suspension into pre-warm 24-well plate using cold pipette tip
 - a. Do not touch the bottom of the plate.
 - b. Slowly and carefully lift up pipette tip as dispensing.
 - c. Do not push pipette tip fully down as this will introduce bubbles.
- 19. Place 24-well plate in 37 °C incubator for 2 minutes to allow Matrigel to pre-polymerize
- 20. Flip the plate and incubate additional 8 minutes to fully polymerize and evenly distribute the organoids fragments through Matrigel.
- 21. During the polymerization process, prepare TOM with Y-27632 with 200 × dilution
 - a. Need 500 ul of medium/well.
 - b. For example, for 6 wells, we need 3 ml of medium with 15 ul of ROCK inhibitor.
- 22. Submerge the polymerized Matrigel domes with 500 µl TOM per well, and culture them in 37 °C humidified CO₂ incubator.
- 23. Observe daily and renew the TOM every 48-72 hrs.

Media recipe:

Trophoblast organoid medium (TOM)

Ingredient	Volume(µl)	Final Concentration
100 × N2 (Life Technologies, 17502-048(500	1×
50 × B27 (Life Technologies, 17504-044)	1000	1×
500 × Primocin (InvivoGen, ant- pm-1)	100	100 µg/ml
80 × NAC (Sigma, A9165)	625	1.25 mM
100 × L-glutamine (Life Technologies, 35050-061)	500	2 mM
10000 × A83-01 (Tocris, 2939)	5	500 nM
10000 × CHIR99021 (Tocris, 4423)	5	1.5 µM
2000 × recombinant hEGF (Gibco, PHG0314)	25	50 ng/ml
2000 × recombinant R-spondin1 (R & D systems, 4645-RS-100)	25	80 ng/ml
2000 × recombinant hFGF2 (Peprotech, 100-18C)	25	100 ng/ml
2000 × recombinant hHGF (Peprotech, 100-39)	25	50 ng/ml
100 × Nicotinamide (NTM) (Sigma, N0636-100G	500	10 mM
<mark>500 × Y-27632</mark> (Sigma, Y0503- 1MG)	250	5 µM ↑

The following is the recipe of preparing 50 mL TOMs (# 6, 4th version)



2000 × PGE2 (R & D systems, 22- 961-0)	25	2.5 µM
FBS (heat inactivated) (Cytiva HyClone, SH30070.03(5 mL	10% (vol/vol)
Advanced DMEM/F12 (Life Technologies, 12634-010)	Adjust to 50 mL	N/A

Annotation: First add about 35 mL Advanced DMEM/F12 to the 50 mL centrifuge tube, then add the above supplements, adjust the final volume to 50 mL with Advanced DMEM/F12. Use the full medium within 1 month. The red highlighted are supplemented components post optimization for full-term placental tissue.