

Cloning outline/ protocol

*This protocol applies for mutagenesis (small insertions, base changes, deletions) using the InFusion Cloning kit only

1. Run gradient PCR to linearize and amplify vector
 - Objective: to figure out which melting temperature (T_m) works best for each primer set
 - Gradient PCR: allows to use a range of T_m s at the same time. One of our PCR machines can do this. I'll show you how to do it.
 - Choose a T_m range: $\sim -2^\circ\text{C}$ of lowest T_m & $+ 2^\circ\text{C}$ of highest T_m
 - Using wild type vector as DNA template
 - Using specific forward and reverse primers that insert mutations of interest -> these primers are unique for each mutant
 - Using Q5 High-Fidelity DNA polymerase
 - Don't forget to calculate the appropriate elongation time. For the Q5 polymerase method this is: 30secs/kb of expected PCR product
 - After setting up the gradient PCR protocol in the PCR machine, don't forget to take a picture of and/or write down T_m s included in the gradient for future reference
2. Run gradient PCR product samples in 1% agarose gel
 - Sample prep: 8ul PCR product + 2ul 6X loading dye (for a 1X final concentration)
 - Don't forget to include a ladder! Recommended ladder: 1kb Plus DNA ladder
 - Run: 1XTBE, 120V, ~ 30 mins
 - Image gel (will show you how) & save picture in usb. Also include in notebook.
3. Large scale PCR
 - Choose the T_m that worked best for each set of primers. Don't forget to write down the chosen T_m for future reference.
 - Run a large scale PCR reaction at the chosen T_m .
 - Using wild type vector as DNA template
 - Using specific forward and reverse primers that insert mutations of interest -> these primers are unique for each mutant
 - Using Q5 High-Fidelity DNA polymerase
 - Don't forget to calculate the appropriate elongation time. For the Q5 polymerase method this is: 30secs/kb

4. PCR cleanup
 - Pull together all PCR aliquots for each mutant and do a PCR cleanup using the Clonetech kit.
 - Follow kit's protocol
 - **Elute in 20ul WATER**

5. Run PCR purified large scale PCR product samples in 1% agarose gel
 - Sample prep: all of the PCR product after PCR cleanup (~20ul) + 20ul 6X loading dye
 - Gel prep: use large well comb (to fit all the sample)
 - Don't forget to include a ladder! Recommended ladder: 1kb Plus DNA ladder
 - Run: 1XTBE, 120V, ~30mins
 - Image gel (will show you how) & save picture in usb. Also include in notebook.
 - **SAVE GEL FOR NEXT STEP**

6. Gel extraction
 - Extract bands of interest using the Zymo Gel Purification kit.
 - Follow kit's protocol
 - Elute in 15ul WATER

7. NanoDrop gel extraction product
 - Objective: to measure the concentration of the linearized DNA product and to confirm the purity of the sample
 - Concentration should be >20ng/ul
 - 260/280 ratio should be ~1.8
 - 260/230 ratio should be ~2.0-2.2

8. InFusion reaction
 - Objective: to re-circularize plasmid
 - Follow protocol on spreadsheet and in InFusion website

9. Transformation of Stellar Competent Cells (E. coli strain)
 - Follow manufacturer's protocol (pdf sent to you)

10. Inoculation of LB cultures with colonies from transformation

- I will teach you how to do this

11. Plasmid prep using Zymo kit

- Objective: to extract re-circularized plasmid from bacteria
- Follow kit's protocol
- Elute in 30ul WATER

12. Sequencing of new mutant plasmids

- Objective: to confirm that the mutations were inserted
- See sequencing spreadsheet for calculations
- Will show you how to prepare the sequencing order and drop off the samples