

# One-pot native barcoding protocol

## Scope:

This is our 'one-pot ligation' protocol for Oxford Nanopore native barcoded ligation libraries using shearing. We have seen no reduction in performance compared to standard libraries yet it can be made faster by using the Ultra II ligation module which is compatible with the Ultra II end repair/dA-tailing module removing a clean-up step. It can be used with or without the optional FFPE DNA repair step. If you have the time I would recommend using the double incubation times in **blue**, if you are in a hurry the times in **red** are a good compromise between speed and efficiency.

## Required:

g-TUBEs (optional)  
SQK-LSK108 1D ligation kit  
Native barcoding kit  
Ultra II End Repair/dA-Tailing Module  
Ultra II Ligation Module  
FFPE DNA Repair Mix (optional)  
Ampure XP beads  
80% ethanol  
EB (10 mM Tris-HCl pH 8)

## Protocol:

Dilute 1 ug gDNA to 50 µl with EB and pipette into Covaris g-TUBE  
Centrifuge g-TUBE at 5000 rpm for 2 x 30 sec spins, invert tube and repeat  
Transfer newly sheared DNA into a clean microcentrifuge tube  
Remove 1 µl and assess concentration by Qubit  
Set up either the following end-prep reaction for each sample:

Sheared DNA (500 ng)	25 µl
Ultra II End Prep Reaction Buffer	3.5 µl
Ultra II End Prep Enzyme Mix	1.5 µl
<b>Total</b>	<b>30 µl</b>

Or to include FFPE DNA repair set up the following combined reaction:

DNA (500 ng)	24 µl
Ultra II End Prep Reaction Buffer	1.75 µl

FFPE DNA Repair Buffer	1.75 µl
Ultra II End Prep Enzyme Mix	1.5 µl
FFPE DNA Repair Mix	1 µl
<b>Total</b>	<b>30 µl</b>

Incubate at RT for **5 mins** or **10 mins** then 65°C for **5 mins** or **10 mins**

Place on ice for 30 secs

Add the following directly to the previous reactions:

NBXX barcode	2.5 µl
Ultra II Ligation Master Mix	20 µl
Ligation Enhancer	1 µl
<b>Total</b>	<b>53.5 µl</b>

Incubate at RT for **10 mins** or **20 mins**, 70°C for 5 mins then place on ice

Pool all barcoded fragments together into a clean 1.5 ml Eppendorf tube

Add 26.75 µl Ampure XP beads per sample

Incubate for **5 mins** or **10 mins**

Place on a magnet rack until clear

Remove supernatant

Add 200 µl 80% ethanol to the tube still on the magnetic rack

Incubate 30 secs

Remove supernatant

Repeat last three steps

Spin down and remove residual 70% ethanol

Air dry 1 min

Resuspend in 31 µl EB

Incubate off the magnetic rack for **5 mins** or **10 mins**

Replace on magnetic rack

Remove 1 µl and assess concentration by Qubit

Set up the following adapter ligation reaction:

Cleaned-up barcoded fragments (~3 µg)	30.0 µl
BAM 1D	20.0 µl
Ultra II Ligation module	40.0 µl
Ultra II Ligation enhancer	1.0 µl
<b>Total volume</b>	<b>91.0 µl</b>

Incubate at RT for 10 mins or 20 mins

Add 45.5 µl Ampure XP beads

Incubate for 5 mins or 10 mins

Place on a magnetic rack until clear

Remove supernatant

Add 150 µl ABB and resuspend by flicking (CAUTION: do not use 80% ethanol)

Place on magnetic rack until clear

Remove supernatant

Repeat ABB wash

Spin down and remove residual ABB

Add 12 µl ELB and resuspend by flicking

Incubate at RT for 5 mins or 10 mins

Place on magnetic rack

In a new tube prepare library dilution for sequencing:

	Without LLB	With LLB
RBF	35 µl	35 µl
Nuclease-free water	28 µl	2.5 µl
LLB	-	25.5 µl
Library	12 µl	12 µl
<b>Total</b>	<b>75 µl</b>	<b>75 µl</b>

Mix well

Remove 1 µl and assess concentration by Qubit (wait until beads have settled before measuring)