

Bacmid minipreps

- (1) Grow a single colony in 3 mL of LB containing 50 $\mu\text{g}/\text{mL}$ kanamycin, 7 $\mu\text{g}/\text{mL}$ gentamicin, 10 $\mu\text{g}/\text{mL}$ tetracycline overnight. It's best to start these cultures in the morning, since the cells grow very slowly.
- (2) Make a glycerol stock by mixing 1 mL of the cells with 120 μL of 80% glycerol. Store in the $-80\text{ }^{\circ}\text{C}$ freezer.
- (3) Transfer the remaining culture to a 2-mL microcentrifuge tube and spin at maximum speed for 1 min.
- (4) Remove all of the supernatant using the vacuum aspirator.
- (5) Resuspend the cells in 300 μL of Qiagen Buffer P1. Do not use P1 containing the blue lysis indicator dye.
- (6) Add 300 μL of Qiagen Buffer P2, then incubate at room temperature for 5 min.
- (7) Slowly add 300 μL of Qiagen Buffer N3, mixing gently during addition. A thick white precipitate of protein and genomic DNA will form. Place the sample on ice for 5-10 min.
- (8) Centrifuge for 10 min at maximum speed. During this step, label another tube and add 800 μL isopropanol to it.
- (9) Gently transfer the supernatant into the tube containing the isopropanol. Avoid transferring any of the white precipitate. Mix gently by inverting the tube a few times and place on ice for 5-10 min. Note: The sample can be stored at $-20\text{ }^{\circ}\text{C}$ overnight (or over a weekend) to get higher yield.
- (10) Centrifuge for 15 min at maximum speed.
- (11) Remove the supernatant and add 500 μL of 70% ethanol. Invert the tube several times to wash the pellet. Centrifuge for 10 min at maximum speed. Repeat the wash step.
- (12) Carefully remove all of the supernatant using the vacuum aspirator. To generate gentle suction, use either a gel-loading pipet tip or a 10- μL tip as nozzle.
- (13) Air-dry the pellet (about 5-10 min). Do not overdry.
- (14) Resuspend the DNA in 40 μL of TE buffer. Do not attempt to add buffer before the ethanol has dried out, or the DNA will not go into solution.
- (15) We generally store the sample at $4\text{ }^{\circ}\text{C}$, until it is used for transfection. After that, it can be transferred to $-20\text{ }^{\circ}\text{C}$ for long-term storage. Avoid repeated freeze-thaw cycles.