

***Bacillus Licheniformis* Enzyme Mix (1 mL per 23 mg tissue):**

100 μ L b. lich 100 mg/mL (10 mg/mL final) (Sigma, P5380)

1 μ L 0.5 M EDTA (Sigma, A8806)

899 μ L DPBS (no Ca, Mg) ThermoFisher (cat. #14190)

Preparing enzymes:

The enzyme is made up in DPBS (#14190). It is aliquoted and stored at -80 °C at 100 mg/mL in 100 μ L aliquots..

Reagents

Enzymes, trypsin inhibitor, BSA and DNase are made up in DPBS (no Ca, no Mg) from Thermo Fisher (14190).

Bovine Serum Albumin - Sigma (A8806).

Hypothermosol FRS

Required supplies:

2 mL dounce homogenizer – Bellco (1984-10002)

Centrifuge for 1.5 mL, 15 mL conicals

Pipettes and pipet tips

15 ml Conicals (MLS)

1.5 mL tubes (MLS)

30 μ M filters - Miltenyi (130-098-458)

Petri dishes (MLS)

Razor blades (MLS)

Ice bucket w/ice

Hemocytometers - InCyto Neubauer Improved (DHC-NO1-5)

Cell yield:

Skin: 3000 cells/mg (69,750 cells total from 23 mg tissue)

Protocol:

1. After euthanizing mouse, remove hair using Nair: dab with Nair, wait 30 secs, wipe with wet paper towel.
2. Isolate tissue and place in ice-cold hypothermosol.
3. Scrape off underlying layer of fatty / connective tissue using scalpel before proceeding.
4. Mince skin tissue thoroughly on petri dish on ice for 3-4 min on ice into 1-mm³ pieces using razor blade while manipulating tissue with forceps – you will need to use grinding motion and vigorously break up tissue.
5. Place 23 mg minced tissue into 1 mL B. Lich enzyme cocktail. Incubate on ice.
6. Shake every min; triturate every 2 min with p1000 w/tip cut (start triturating at 2 min).

7. After 20 mins of triturating on ice, use pipet to transfer digest mix to 2 mL dounce homogenizer. Use 10 strokes of Pestle A every 2 min (4 series total, 8 min). Digest mix should become turbid.
8. Transfer back to 1.5 mL tube using 1 mL serological pipet. Mix thoroughly and allow to settle on ice 2 min.
9. Save 70% (700 μ L) of supernatant, leaving chunks at the bottom of the tube; apply to 30 μ M filter on 15 mL conical. Rinse filter w/5 mL ice-cold PBS/BSA 0.04%. Save flow through on ice and keep filter on tube for 2nd layer.
10. Add additional 1 mL b. Lich enzyme mix to residual tissue chunks.
11. Triturate every 2 min, shake every min while incubating on ice for 20 additional mins. (50 min. total digest time).
12. Transfer entire volume to same 30 μ M filter on 15 mL conical. Rinse with additional 5 mL ice-cold PBS/BA 0.04%.
13. Centrifuge at 300 g for 5 min at 4 °C. Remove supernatant & re-suspend in 100 μ L PBS/BSA 0.04%. Examine using hemocytometer with trypan blue.
14. Adjust concentration to 1,000 cells/ μ L for Chromium or 100 cells/ μ L for DropSeq.