

D68EV3C large scale purification

PAGE24-00037

Author: **Wang, Korvus**

Date Started: **2024-Jan-10**

Experiment Started:

Projects: **Purification;ASAP**

Related Pages: **PAGE23-01458**

Referenced by:

Tags:

Notes and Thoughts from the last large scale...

Potential ways to remove DNA:

- PEI

need to check if there are clusters of negatively charged region in the protein - matlab bioinf toolbox

Doesn't look like it

Is the protein DNA binding?

No, 3C is a protease

Only concern is target is slightly negatively charged overall (PI=7.21) at HEPES buffer pH (7.4)

- DEAD52

Ask Mike if he knows where it is

- Benz soak

After loading onto column, wash with binding buffer. Soak resin in buffer with 1:4000 benz, 2mM MgCl₂, 150-200mM NaCl depending on what the protein can take. Leave capped in cold room for about 1hr, then wash with high salt buffer.

Other alterations:

Maybe raise the pH of the buffer to around 8 for IMAC, and back down to 7.4 during gel filtration

This might help with the precipitation issue as well

Protein has two free cysteines, which may cause homodimer formation during purification (need more TCEP?)

D683C sequence analysis

```
1  MGPGFDFAQA  IMKKNTVIAR  TEKGEFTMLG  VYDRVAVIPT  HASVGEIIIYI  NDVETRVLDA
61  CALRDLTDTN  LEITIVKLDR  NQKFRDIRHF  LPRCEDDYND  AVLSVHTSKF  PNMYIPVGQV
121 TNYGFLNLGG  TPTHRILMYN  FPTRAGQCGG  VVTTTGKVIG  IHVGGNGAQG  FAAMLLHSYF
181  TDTQKHHHHHH  H
```

SwissModel showing cysteine locations



IMAC

Expression ID: D68EV3CPROA-e006 (Nathan)

Purification ID: D68EV3CPROA-p007

D68EV3CPROA-c001

MGPGFDFAQAIMKKNTVIARTEKGEFTMLGVYDRVAVIPTHASVGEIYINDVETRVLDACALRDLTDTNLEITIVKLDNRNQKFR
DIRHFLPRCEDDYNDVAVLSVHTSKFPNMYIPVGQVTNYGFLNLGGTPTHRILMYNFPTRAGQCGGVTTTGKVIKIHVGGNGA
QGFAAMLLHSYFTDTQKHHHHHH

21283.3 Da

10430 mM-1cm-1

PI=7.21

Buffers:

Lysis buffer - 50 mM HEPES pH 8.0, 500 mM NaCl, 5% glycerol, 0.5mM TCEP, 20mM imidazole

Wash buffer 1 - 50 mM HEPES pH 8.0, 100 mM NaCl, 5% glycerol, 0.5mM TCEP, 2mM MgCl, 1:2000 benzonase

Wash buffer 2 - 50 mM HEPES pH 8.0, 500 mM NaCl, 5% glycerol, 0.5mM TCEP, 30mM imidazole

Elution Buffer - 50 mM HEPES pH 8.0, 500 mM NaCl, 5% glycerol, 0.5mM TCEP, 500mM imidazole

Gel Filtration Buffer (SEC) - 50 mM HEPES pH 7.5, 500 mM NaCl, 5% glycerol, 0.5mM TCEP

Expression was carried out by Nathan Wright in bubble columns. 10L culture total. 244g pellet harvested.

1. 244g of pellet was thawed and resuspended with enough lysis buffer so that the final volume is 2L. Supplement with 1:4000 dilution of benzonase, 1mg/mL lysozyme, and 2mM MgCl. Incubate for 30min at RT.
1. Sonicated on ice at 50% amplitude for a total of 7-minute sonication time (4 seconds on 12 seconds off) with thick probe. Sonication done in 2 batches.
2. Clarified lysate by centrifugation at 38,000rpm, 4°C for 1 hour. Used JLA16.250 rotor, done in 2 batches. Supernatant poured into 2L beaker.

IMAC

1. Wash and equilibrate 30mL bed volume of Ni Sepharose resin on large gravity flow column, first with distilled water, then with wash buffer.
2. Put stopper on column. Re-suspend resin with small volumes of the lysis buffer, and pour into beaker with lysate. Stir in cold room gently for 30min.
3. Pour mixture back onto gravity flow column and let flow through.
4. Add 200mL of wash buffer 1. Cap the column and leave rotating for 30min to digest remaining DNA. Drain wash buffer 1 and add 200mL wash buffer 2. Allow to flow through.
5. Elute with 30mL elution buffer, 3 elutions carried out
Elution 2 started to precipitate mid-elution!! perhaps because concentration was too high? Immediately diluted with 20mL of base buffer and the precipitants appear to have re-solubilised.
Made all the elutions up to 50mL to prevent them from crashing overnight.

68

E1: 1.16 mg/mL


E2: 6.22 mg/mL

E3: 3.24 mg/mL

A260 looks good across all elutions. DNA issue resolved.

Surprisingly lower than expected, since elution 2 was crashing out?

IAMC results



Sample concentration for gel filtration

Elutions pooled and diluted to 500mL with base buffer to reduce the final [imidazole] after concentration
 Attempted to concentrate with Vivaspin 50 (10kDa MWCO) tangential flow concentrators.


Sample spurted everywhere at one point

Also the peristaltic pump generates a lot of heat and kind of cooked some of the protein.

Ended up using 3 cassettes and had to bury them in ice. Still a thick layer of jelly-like precipitated protein on top.

After reaching ~90mL, stopped concentration. Spun samples down at 20,000xg in JLA25. Ran sample and actually looks like most of what crashed out was the contaminants?

Concentrated sample



Gel filtration run

Decided to use XK 50/100 Superdex 200 pg column. 2L bed volume and loads up to 30mL sample. 500mg sample load can achieve analytical resolution.

MAKE SURE TO FLANGE ALL TUBINGS DUE TO HIGH PRESSURE

Run 30mL first to check if sample crashes at pH 7.5.

Loaded sample using 50mL superloop, prepared using the following video as guide

<https://www.youtube.com/watch?v=EavIht55bq0>

Collect 35mL fractions with falcon tubes in the fraction collector (6 cassettes). Start collection after 0.4CV (anything before 0.3 is void)

All looks good. No pressure change which indicated no significant precipitation. Sample came out near the end. Sample falcon tubes ran out near the end but is fine cuz all the peaks already came out.

Further concentrated the remaining samples down to ~45mL. Clarified by centrifugation at 20,000xg and injected via 0.2µM syringe filter.

Also ran fine. Peak profile looks similar to previous run.

Run 1 Fractions:

3B3: 1.74 mg/mL


4A1: 2.25 mg/mL

4A2: 0.58 mg/mL


Run 2 Fractions:
3B2: 2.68 mg/mL
3B3: 5.26 mg/mL
4A1: 2.73 mg/mL
4A2: 0.56 mg/mL

Fractions ran on gel and concentrated in Amicon 10kDa MWCO concentrators.


Run 1



Run 2



SEC results



Final sample

Final sample:

22.54 mg/mL

241x100uL

543.2 mg total

NEED FINAL SAMPLE GEL

Previous purification (p004) got 128.3mg from 2L.

Expected from 10L: $128.3 \times 5 = 641.5$ mg

Achieved 84.5% of expected yield

probably lost some during initial concentration

Conclusion: scale-up is fairly comparable to original shake-flask method..

D68EV3CPROA-p007

Good mass!!

