Lysosomal flux assay (as reported by Giacomo Monzio Compagnoni)

### PURPOSE

To assess the autophagic flux in cells, by evaluating LC3II and p62 amount before and after bafilomycin treatment

## MATERIALS AND REAGENTS

- Western blot running tank
- iBlot transfer system
- LI-COR Clx scanner
- bench centrifuge
- 4-12 % Bis-Tris gel
- iBlot nitrocellulose stacks
- LDS 4X buffer
- Protease/phosphatase inhibitors
- Bafilomycin
- PBS
- Odyssey blocking buffer
- tween20
- LC3 antibody (Cell Signaling 2775)
- p62 antibody (Millipore MABN130)
- Actin antibody (Sigma A2066)
- Secondary antibodies (800 anti Rb, 800 anti Ms, 680 anti Rb)

## PROTOCOL

## **Bafilomycin treatment:**

- 1. Change medium to complete medium + bafilomycin (eg 50nM, 100 nM, 200nM)
- 2. Treat control wells with complete medium + DMSO (eg 500X, to match bafilomycin dilution) (bafilomycin aliquots are resuspended in DMSO)
- 3. Keep at 37° until designated collection timepoint(s) (e.g. 12 h)
- 4. Collect cell pellet

## Protein extraction:

- 5. Keep samples on ice throughout extraction
- 6. Dilute 4X blue LDS buffer (Cat. no. B0007, Life Tech) in water plus protease inhibitors to get 1X LDS buffer
- 7. Resuspend each pellet in 1X LDS buffer( $100\mu$ L, but reduce or increase the volume according to the pellet size)
- 8. Sonicate twice for 15 seconds at 50% power, keeping sample on ice
- 9. Boil for 5 minutes at 100°C; return directly to ice
- 10. Centrifuge for 10 minutes at 850xg

- 11. Retain the supernatant
- 12. Perform BCA assay on protein samples; dilute the standards and the blank in 1X LDS buffer
- 13. Store lysates at -80°C

# Western blot:

- 14. Calculate the volume of each sample containing 30  $\mu$ g of proteins; add 1X LDS to bring volume to 22.5  $\mu$ L; add 2.5  $\mu$ L Thermo Fisher reducing reagent per sample (if using a 10-wells gel)
- 15. Boil samples at 100°C for 5min
- 16. Load 25  $\mu$ L per well into 10-well mini Bis-Tris 4-12% gel(s)
- 17. Dilute protein ladder in 1X LDS (e.g. 4  $\mu$ l protein ladder + 16  $\mu$ l 1X LDS buffer)
- 18. Run in MES buffer (1X), 30 minutes, at 200V
- 19. Cut off wells and bottom of gel; move gel directly to transfer stack
- 20. Transfer using iBlot P0 program to nitrocellulose membrane (Cat. no. IB23002, Thermo)
- 21. Do not touch membrane with gloves-use forceps and razor
- 22. Cut off edges (perimeter of gel)
- 23. Cut across sample section between 2<sup>nd</sup> and 3<sup>rd</sup> (100kDa) ladder bands from top (or somewhere else above 75 KDa)
- 24. Cut across samples at 25kDa (lower red) band (middle of the band)
- 25. Rehydrate membrane in PBS for 5min on orbital shaker
- 26. Block 1 hour at room temp in Licor Odyssey Buffer PBS (Cat. no. 927-40000, Licor)
- 27. Prepare primary antibody solutions in Odyssey plus 0.1% Tween:
  - HMW: rabbit-anti-actin 1:1200 (-20°C) + mouse-anti-p62 1:1000 (labeled "SQSTM1", at 4°C)
  - LMW: rabbit-anti-LC3 1:1000
- 28. Recover blocking solution to be used for secondary antibodies—keep on ice or at 4°C
- 29. Incubate in primary antibody solutions for 2 hours at room temp on orbital shaker or at 4°C overnight
- 30. Wash 4x 5 minutes with PBS-T (0.05% Tween)
- 31. Incubate secondary antibodies at 1:10,000 dilution in Licor Odyssey plus 0.1% Tween, 1 hour at room temp in black box or aluminium foil:
  - o HMW: 680-anti-rabbit, 800-anti-mouse
  - o LMW: 800-anti-rabbit
- 32. Wash 4x 5 minutes with PBS-T (0.05%), in black box or aluminium foil
- 33. Wash 1x 5 minutes with PBS
- 34. Change to fresh PBS
- 35. Immediately before acquisition, dry membranes on kimwipe
- 36. Reassemble membrane and scan with Licor Clx scanner
  - Flip membranes so that lower left corner is in upper left

- 37. Analysis: LC3 II normalized over Actin (or LC3 I), then divided by baseline (DMSO condition)
  - LC3 II is larger than LC3 I, but charge makes it run faster: ratio is lower band divided by upper band (or Actin)
  - Normalize p62 over Actin to corroborate LC3

•