

CVMPROA-c001 p003 SARSMpro in pNIC transformation, expression and purification

PAGE24-00585

Author: **Wang, Korvus**

Date Started: **2024-Mar-22**

Experiment Started:

Projects: **Expression;Purification;ASAP**

Related Pages: **PAGE23-00603;PAGE23-00782**

Referenced by:

Tags:

Transformation

Plasmid was purchased from Twist and dissolved in 10uL of nuclease-free water. Transformation carried out per the standard protocol.

CVMPROA-c001

Expression

22/03/2024

6x1L Invitrogen TB + Kan + 4mL glycerol after autoclave + Kan50
1mL of ON from transformation was used to inoculate per L of TB. Grown at 37oC 200rpm until OD~1.8, induced with 0.5mM IPTG
Temperature reduced to 18oC and speed to 180rpm. Grown overnight.

Harvest next day. Total pellet weight = 88g

IMAC, cleavage and rIMAC

MGSSHHHHHHGSDSEVNQEAKPEVKPEVKPETHINLKVSDGSSEIFFKIKKTTPLRRLMEAFKRQ GKEMDSLRFlyDGIRIQA
DQTPEDLDMEDNDIIEAHREQIGGSGFRKMAFPSGKVEGCMVQVTCGTTTLNGLWLDDVVYCPRHVICTSEDMLNPNYEDLL
IRKSNHNFLVQAGNVQLRVIGHSMQNCVLKLVDTANPKTPKYKFVRIQPGQTFSVLACYNGSPSGVYQCAMRPNFTIKGSFL
NGSCGSVGFNIDYDCVSFCYMHMELPTGVHAGTDLEGNFYGPFVDRQTAQAAGTDTTITVNVLAWLYAAVINGDRWFLNR
FTTTLNDFNLVAMKYNYEPLTQDHVDILGPLSAQTGIAVLDMCASLKELLQNGMNGRTILGSALLEDEFTPFDVVRQCSGVTFQ

Purification ID: CVMPROA-p001 (labelled p003 on bag)

Buffers:

Lysis buffer - 50 mM HEPES pH 7.4, 150 mM NaCl, 5% glycerol, 0.5mM TCEP

Wash buffer - 50 mM HEPES pH 7.4, 150 mM NaCl, 5% glycerol, 0.5mM TCEP

Elution Buffer - 50 mM HEPES pH 7.4, 150 mM NaCl, 5% glycerol, 0.5mM TCEP, 500mM imidazole

Gel Filtration Buffer (SEC) - 50 mM HEPES pH 7.4, 150 mM NaCl, 5% glycerol, 0.5mM TCEP

no imidazole in binding and wash step due to previous difficulties with target not binding to IMAC resin

1. Lysis buffer supplement with 1:4000 dilution of benzonase, 0.5 mg/mL lysozyme. Incubate for 30min at RT.
2. Sonicated on ice at 50% amplitude for a total of 7-minute sonication time (4 seconds on 12 seconds off) with thick probe.
3. Clarified lysate by centrifugation at 18,000rpm, 4°C for 1 hour. Used JLA16.250 rotor. Supernatant poured into clean beaker.

IMAC

1. Wash and equilibrate 5mL bed volume of Ni Sepharose resin on gravity flow column, first with distilled water, then with wash buffer.
2. Resuspend resin with lysis buffer and add to beaker containing clarified supernatant.
3. Stir gently in cold room for 30mins. Pour only gravity flow column and allow to flow through.
4. Wash resin with 25mL wash buffer twice.
5. Elute with 7.5mL elution buffer, 2 elutions carried out



TEV cleavage

1. Desalt elution with HiLoad 10/26 desalting column on ATKA. Desalt into lysis buffer.
2. Pool desalted protein from the fractions
3. Add SENP1 in 1:300 OD ratio. Leave to incubate in cold room.


rIMAC and Gel filtration

1. in morning, IMAC resin washed with some lysis buffer to remove imidazole.
2. cleavage mix passed through the resin 1 time (flow rate was a bit slow this time)
3. Wash resin with 10mL wash buffer
4. Elute to see what stuck to the resin
5. rIMAC FT concentrated in 10kDa MWCO yellow vivaspin concentrators to final volume of ~5mL
6. Injected onto Superdex s200 16/60 column and run in SEC buffer at 1mL/min
7. After SEC, peaks were analysed by SDS-PAGE

IMAC results



SEC result




Final sample

Fraction 4C5 - 4D10 were pooled and concentrated in the yellow vivaspin 10kDa MWCO concentrators

Final concentration: 26.96 mg/mL, 47*100 + 50 uL

Total yield: 128.1 mg

Title missing - double click to edit



final sample

