

Isolating Placental Organoids from Full-term human placental tissue

DOI: dx.doi.org/10.17504/protocols.io.5jyl8943dv2w/v1

Reagents needed:

Wash buffer: HAMS/F12 medium (Thermo Fisher, 11765054) supplemented with 1x pen/strep

Digestion stop solution: 20% FBS/wash medium

0.2% Trypsin-250/0.02% EDTA/PBS solution:

0.3g glucose, 12g NaCl, 0.3g KCl, 1.725g disodium hydrogen orthophosphate, 0.3g potassium dihydrogen orthophosphate, 2g trypsin (Alfa Aesar, J63993-09), 0.2g EDTA (Sigma-Aldrich, E9884-100G) made up in 1L of water. Filter, aliquot, and store for up to 6mo at -20°C.

Collagenase V digest solution:

1.0 mg/mL Collagenase V (Sigma-Aldrich, C9263-100MG) in Ham's F12/K with 10% FBS. Make fresh for each use.

<u>Matrigel stocks</u>: Stocks should be thawed overnight at 4°C.

Tissue dissection and digestion

Pre-warm 0.2% Trypsin/0.02% EDTA and Collagenase V solutions at 37°C

- 1. Dissect placental chorionic villi. Carefully remove decidua.
- 2. Generate very small fragments of tissue.
- 3. Place tissue fragments in glass bottle containing wash buffer and a stir bar and wash extensively on a magnetic stirrer (VWR, 76447-030) at room temperature as much as possible (at least 5 times). Continue to wash until wash buffer is clear.
- 4. Allow tissue to settle to bottom of glass bottom by gravity and then remove wash buffer using a vacuum aspirator.
- 5. Pour 25-75mL (depending on quantity of fragments) of pre-warmed 0.2% Trypsin/0.02% EDTA solution into the same glass bottle containing a small stir bar. Place bottle on a heated (37°C) stir plate with gentle stirring (set at ~ 80 rpm). Incubate for 8min.
- 6. Filter the suspension using muslin gauze placed over a funnel. Add digestion stop solution to neutralize trypsin. Save the flow through as you will use this in a subsequent step. Be sure to save the undigested tissue, which will be used in the next step.
- Retrieve the undigested tissue from the muslin gauze and place in 12mL of collagenase V solution in a glass bottle containing a stir bar. Place the bottle on a heated stir plate and stir with gentle agitation (~ 80 rpm) for 8min.
- Remove the glass bottle from the stir plate and manually disrupt tissue using a narrow orifice 10mL serological pipette (VWR, 89130-898). Forcefully pipette up and down ~10 times to break up tissue prior to filtration.



- 9. Filter the suspension as described in step 6. Remaining undigested tissue on the gauze can be discarded.
- 10. Pool the cells collected in steps 6 and 9.
- 11. Centrifuge at 400g for 5min.
- 12. Resuspend the pellet in 5mL of Advanced DMEM/F12 medium (Life Technologies 12634-010) and transfer to a 15mL conical.
- 13. Centrifuge at 600g for 5-6min to pellet.
- 14. Carefully remove as much of the supernatant as possible. *Do not disrupt the pellet*. It is not recommended to use vacuum aspiration for this step being careful not to disrupt the pellet.
- 15. Resuspend the pellet with pre-thawed Matrigel (Corning 356231) using blunt 200 μl pipette tips (Fisher 02-707-134). The volume will depend on how many wells will be plated. 40uL of Matrigel is used per well. Use pre-chilled large orifice pipette tips to prevent Matrigel polymerization due to temperature.
- 16. Carefully dispense 40uL of cells/Matrigel to the middle of each well of a pre-warmed 24-well plate to create a "dome".
- 17. Carefully place the plate in the 37°C incubator for 3min to allow for pre-polymerization.
- 18. Flip the plate upside down to disperse cells throughout the Matrigel dome. Incubate for 8min.
- 19. Add 500uL of growth medium containing ROCK inhibitor Y-27632 (to prevent stem cell death) and culture for 72-96h in this medium. After this time, add fresh medium without Y-27632.

Media recipe:

Trophoblast organoid medium (TOM)

Ingredient	Volume(µl)	Final Concentration
100 × N2 (Life Technologies, 17502-048(500	1×
50 × B27 (Life Technologies, 17504-044)	1000	1×
500 × Primocin (InvivoGen, ant- pm-1)	100	100 µg/ml
80 × NAC (Sigma, A9165)	625	1.25 mM
100 × L-glutamine (Life Technologies, 35050-061)	500	2 mM
10000 × A83-01 (Tocris, 2939)	5	500 nM
10000 × CHIR99021 (Tocris, 4423)	5	1.5 µM
2000 × recombinant hEGF (Gibco, PHG0314)	25	50 ng/ml
2000 × recombinant R-spondin1 (R & D systems, 4645-RS-100)	25	80 ng/ml
2000 × recombinant hFGF2 (Peprotech, 100-18C)	25	100 ng/ml
2000 × recombinant hHGF (Peprotech, 100-39)	25	50 ng/ml
100 × Nicotinamide (NTM) (Sigma, N0636-100G	500	10 mM

The following is the recipe of preparing 50 mL TOMs (# 6, 4th version)



500 × Y-27632 (Sigma, Y0503- 1MG)	250	5 µM ↑
2000 × PGE2 (R & D systems, 22- 961-0)	25	2.5 µM
FBS (heat inactivated) (Cytiva HyClone, SH30070.03(5 mL	10% (vol/vol)
Advanced DMEM/F12 (Life Technologies, 12634-010)	Adjust to 50 mL	N/A

Annotation: First add about 35 mL Advanced DMEM/F12 to the 50 mL centrifuge tube, then add the above supplements, adjust the final volume to 50 mL with Advanced DMEM/F12. Use the full medium within 1 month. The red highlighted are supplemented components post optimization for full-term placental tissue.