

Direct cDNA Sequencing (SQK-DCS109)

Version: DCS_9090_v109_revL_14Aug2019
 Last update: 22/12/2020



Flow Cell Number:

DNA Samples:

Before start checklist		
Materials	Consumables	Equipment
<input type="checkbox"/> 100 ng PolyA+ RNA, or 70-200 ng already-prepared cDNA <hr/> <input type="checkbox"/> Direct cDNA Sequencing Kit (SQK-DCS109) <hr/> <input type="checkbox"/> Flow Cell Priming Kit (EXP-FLP002)	<input type="checkbox"/> Agencourt AMPure XP beads <hr/> <input type="checkbox"/> NEBNext End repair / dA-tailing Module (E7546) <hr/> <input type="checkbox"/> NEB Blunt/TA Ligase Master Mix (M0367) <hr/> <input type="checkbox"/> 1.5 ml Eppendorf DNA LoBind tubes <hr/> <input type="checkbox"/> 0.2 ml thin-walled PCR tubes <hr/> <input type="checkbox"/> Nuclease-free water (e.g. ThermoFisher, cat # AM9937) <hr/> <input type="checkbox"/> Freshly prepared 70% ethanol in nuclease-free water <hr/> <input type="checkbox"/> 10 mM dNTP solution (e.g. NEB N0447) <hr/> <input type="checkbox"/> LongAmp Taq 2X Master Mix (e.g. NEB M0287) <hr/> <input type="checkbox"/> Maxima H Minus Reverse Transcriptase (200 U/μl) with 5x RT Buffer (ThermoFisher, cat # EP0751) <hr/> <input type="checkbox"/> RNaseOUT™, 40 U/μl (Life Technologies, 10777019) <hr/> <input type="checkbox"/> RiboShredder (Epicentre, RS12500), or RNase Cocktail Enzyme Mix (ThermoFisher, AM2286)	<input type="checkbox"/> Hula mixer (gentle rotator mixer) <hr/> <input type="checkbox"/> Magnetic separator, suitable for 1.5 ml Eppendorf tubes <hr/> <input type="checkbox"/> Microfuge <hr/> <input type="checkbox"/> Vortex mixer <hr/> <input type="checkbox"/> Thermal cycler <hr/> <input type="checkbox"/> Ice bucket with ice <hr/> <input type="checkbox"/> Timer <hr/> <input type="checkbox"/> Pre-chilled freezer block at -20° C for 200 μl tubes (e.g. Eppendorf 022510509) <hr/> <input type="checkbox"/> Multichannel pipette capable of 20-200 μl <hr/> <input type="checkbox"/> Pipettes and pipette tips P2, P10, P20, P100, P200, P1000
INSTRUCTIONS		NOTES/OBSERVATIONS
Reverse transcription and strand-switching		
IMPORTANT		
<input type="checkbox"/> If you have already prepared your cDNA, use 70-200 ng cDNA and start from the End-prep step.		
Prepare the RNA in Nuclease-free water <input type="checkbox"/> Transfer 100 ng PolyA+ RNA into a 1.5 ml Eppendorf DNA LoBind tube <input type="checkbox"/> Adjust the volume to up to 7.5 μl with Nuclease-free water <input type="checkbox"/> Mix by flicking the tube to avoid unwanted shearing <input type="checkbox"/> Spin down briefly in a microfuge		

Direct cDNA Sequencing (SQK-DCS109)

Version: DCS_9090_v109_revL_14Aug2019
 Last update: 22/12/2020



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Prepare the following reaction in a 0.2 ml PCR tube:</p> <ul style="list-style-type: none"> <input type="checkbox"/> x µl poly A+ RNA, 100 ng <input type="checkbox"/> 2.5 µl VNP <input type="checkbox"/> 1 µl 10 mM dNTPs <input type="checkbox"/> 7.5-x µl RNase-free water <p><input type="checkbox"/> Mix gently by flicking the tube, and spin down.</p> <p><input type="checkbox"/> Incubate at 65° C for 5 minutes and then snap cool on a pre-chilled freezer block.</p> <p>In a separate tube, mix together the following:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 4 µl 5x RT Buffer <input type="checkbox"/> 1 µl RNaseOUT <input type="checkbox"/> 1 µl Nuclease-free water <input type="checkbox"/> 2 µl Strand-Switching Primer (SSP) <p><input type="checkbox"/> Mix gently by flicking the tube, and spin down.</p> <p><input type="checkbox"/> Add the strand-switching buffer to the snap-cooled, annealed mRNA, mix by flicking the tube and spin down.</p> <p><input type="checkbox"/> Incubate at 42° C for 2 minutes.</p> <p><input type="checkbox"/> Add 1 µl of Maxima H Minus Reverse Transcriptase. The total volume is now 20 µl.</p> <p><input type="checkbox"/> Mix gently by flicking the tube, and spin down.</p> <p>Incubate using the following protocol:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Reverse transcription and strand-switching 90 mins @ 42° C (1 cycle) <input type="checkbox"/> Heat inactivation 5 mins @ 85° C (1 cycle) <input type="checkbox"/> Hold @ 4° C 	
<p>RNA degradation and second strand synthesis</p>	
<ul style="list-style-type: none"> <input type="checkbox"/> Add 1 µl RiboShredder or RNase Cocktail Enzyme Mix (ThermoFisher, AM2286) to the reverse transcription reaction. <input type="checkbox"/> Incubate the reaction for 10 minutes at 37° C. <input type="checkbox"/> Resuspend the AMPure XP beads by vortexing. <input type="checkbox"/> Transfer the sample to a clean 1.5 ml Eppendorf DNA LoBind tube. <input type="checkbox"/> Add 17 µl of resuspended AMPure XP beads to the reaction and mix by flicking the tube. <input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT. <input type="checkbox"/> Prepare 500 µl of fresh 70% ethanol in Nuclease-free water. <input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant. 	

Direct cDNA Sequencing (SQK-DCS109)

Version: DCS_9090_v109_revL_14Aug2019
 Last update: 22/12/2020



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><input type="checkbox"/> Keep the tube on the magnet and wash the beads with 200 µl of freshly prepared 70% ethanol without disturbing the pellet. Remove the ethanol using a pipette and discard.</p> <p><input type="checkbox"/> Repeat the previous step.</p> <p><input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual ethanol. Allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking.</p> <p><input type="checkbox"/> Remove the tube from the magnetic rack and resuspend pellet in 20 µl Nuclease-free water.</p> <p><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 10 minutes at RT.</p> <p><input type="checkbox"/> Pellet beads on magnet until the eluate is clear and colourless.</p> <p><input type="checkbox"/> Remove and retain 20 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube.</p> <p>Prepare the following reaction in a 0.2 ml thin-walled PCR tube:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 25 µl 2x LongAmp Taq Master Mix <input type="checkbox"/> 2 µl PR2 Primer (PR2) <input type="checkbox"/> 20 µl Reverse-transcribed sample from above <input type="checkbox"/> 3 µl Nuclease-free water <p>Incubate using the following protocol:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 94 °C 1 mins 1 <input type="checkbox"/> 50 °C 1 mins 1 <input type="checkbox"/> 65 °C 15 mins 1 <input type="checkbox"/> 4 °C ∞ <p><input type="checkbox"/> Resuspend the AMPure XP beads by vortexing.</p> <p><input type="checkbox"/> Transfer the sample to a clean 1.5 ml Eppendorf DNA LoBind tube.</p> <p><input type="checkbox"/> Add 40 µl of resuspended AMPure XP beads to the reaction and mix by flicking the tube.</p> <p><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT.</p> <p><input type="checkbox"/> Prepare 500 µl of fresh 70% ethanol in Nuclease-free water.</p> <p><input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant.</p> <p><input type="checkbox"/> Keep the tube on the magnet and wash the beads with 200 µl of freshly prepared 70% ethanol without disturbing the pellet. Remove the ethanol using a pipette and discard.</p> <p><input type="checkbox"/> Repeat the previous step.</p> <p><input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual ethanol. Allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking.</p> <p><input type="checkbox"/> Remove the tube from the magnetic rack and resuspend pellet in 21 µl Nuclease-free water.</p> <p><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 10 minutes at RT.</p>	

Direct cDNA Sequencing (SQK-DCS109)

Version: DCS_9090_v109_revL_14Aug2019
 Last update: 22/12/2020



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<ul style="list-style-type: none"> <input type="checkbox"/> Pellet beads on magnet until the eluate is clear and colourless. <input type="checkbox"/> Remove and retain 21 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube. <input type="checkbox"/> Analyse 1 µl of the strand-switched DNA for size, quantity and quality. 	
<p>End-prep</p>	
<p>IMPORTANT</p> <p><input type="checkbox"/> If you have prepared your own cDNA instead of performing reverse transcription using the Direct cDNA Sequencing Kit, please start this step with 70-200 ng cDNA in 20 µl Nuclease-free water.</p> <p>Perform end repair and dA-tailing of fragmented DNA as follows:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 20 µl cDNA sample <input type="checkbox"/> 30 µl Nuclease-free water <input type="checkbox"/> 7 µl Ultra II End-prep reaction buffer <input type="checkbox"/> 3 µl Ultra II End-prep enzyme mix <ul style="list-style-type: none"> <input type="checkbox"/> Mix gently by pipetting and spin down. <input type="checkbox"/> Using a thermal cycler, incubate at 20°C for 5 mins and 65°C for 5 mins. <input type="checkbox"/> Resuspend the AMPure XP beads by vortexing. <input type="checkbox"/> Transfer the sample to a 1.5 ml DNA LoBind Eppendorf tube. <input type="checkbox"/> Add 60 µl of resuspended AMPure XP beads to the end-prep reaction and mix by pipetting. <input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT. <input type="checkbox"/> Prepare 500 µl of fresh 70% ethanol in Nuclease-free water. <input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant. <input type="checkbox"/> Keep the tube on the magnet and wash the beads with 200 µl of freshly prepared 70% ethanol without disturbing the pellet. Remove the ethanol using a pipette and discard. <input type="checkbox"/> Repeat the previous step. <input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual ethanol. Allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking. <input type="checkbox"/> Remove the tube from the magnetic rack and resuspend pellet in 30 µl Nuclease-free water. Incubate for 2 minutes at RT. <input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless. <input type="checkbox"/> Remove and retain 30 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube. 	
<p>Take forward 30 µl of end-prepped cDNA into adapter ligation.</p>	

Direct cDNA Sequencing (SQK-DCS109)

Version: DCS_9090_v109_revL_14Aug2019
 Last update: 22/12/2020



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Adapter ligation</p> <p>Check the contents of each tube are clear of any precipitate and are thoroughly mixed before setting up the reaction.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Mix the contents of each tube by flicking <input type="checkbox"/> Check that there is no precipitate present (DTT in the Blunt/TA Master Mix can sometimes form a precipitate) <input type="checkbox"/> Spin down briefly before accurately pipetting the contents in the reaction <p>Taking the end-prepped DNA, perform adapter ligation as follows, mixing by flicking the tube between each sequential addition.</p> <ul style="list-style-type: none"> <input type="checkbox"/> 30 µl End-prepped DNA <input type="checkbox"/> 5 µl Adapter Mix <input type="checkbox"/> 50 µl Blunt/TA Ligation Master Mix <input type="checkbox"/> 15 µl Nuclease-free water <ul style="list-style-type: none"> <input type="checkbox"/> Mix gently by flicking the tube, and spin down. <input type="checkbox"/> Incubate the reaction for 10 minutes at RT. 	
<p>Adapted and tethered DNA library.</p>	
<p>AMPure XP bead binding</p> <ul style="list-style-type: none"> <input type="checkbox"/> Resuspend the AMPure XP beads by vortexing. <input type="checkbox"/> Add 40 µl of resuspended AMPure XP beads to the adapter ligation reaction from the previous step and mix by pipetting. <input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT. <input type="checkbox"/> Place on magnetic rack, allow beads to pellet and pipette off supernatant. <input type="checkbox"/> Add 200 µl of Wash Buffer (WSB) to the beads. Resuspend the beads by pipetting up and down. Return the tube to the magnetic rack, allow beads to pellet and pipette off the supernatant. <input type="checkbox"/> Repeat the previous step. <input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual supernatant. <input type="checkbox"/> Remove the tube from the magnetic rack and resuspend pellet in 25 µl of Elution Buffer (EB). <input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 10 minutes at RT. <input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless. <input type="checkbox"/> Remove and retain 25 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube. <input type="checkbox"/> Quantify 1 µl of eluted cDNA using a Qubit fluorometer - recovery aim ~60 fmol. 	

Direct cDNA Sequencing (SQK-DCS109)

Version: DCS_9090_v109_revL_14Aug2019
 Last update: 22/12/2020



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>IMPORTANT</p> <p>Please be aware that the flow cell's pore occupancy could be compromised when loading lower amounts of cDNA. Please use the table below as a guide:</p> <p><input type="checkbox"/> Please check the Mass to Molarity table in the protocol</p>	
<p>The prepared library is used for loading onto the flow cell. Store the library on ice until ready to load.</p>	
<p>Priming and loading the flow cell</p>	
<p><input type="checkbox"/> Thaw the Sequencing Buffer (SQB), Loading Beads (LB), Flush Tether (FLT) and one tube of Flush Buffer (FB) at RT.</p> <p><input type="checkbox"/> Mix the Sequencing Buffer (SQB), Flush Buffer (FB) and Flush Tether (FLT) tubes by vortexing and spin down at RT.</p> <p><input type="checkbox"/> To prepare the flow cell priming mix, add 30 µl of thawed and mixed Flush Tether (FLT) directly to 1 tube of thawed and mixed Flush Buffer (FB), and mix by vortexing.</p> <p><input type="checkbox"/> Load the flow cell(s) into the docking ports within the PromethION.</p> <p>Prime the flow cell using the following steps, taking care to avoid the introduction of air bubbles.</p> <p><input type="checkbox"/> Turn the valve to expose the inlet port (Port 1)</p> <p><input type="checkbox"/> A small tract of air will be visible beyond the inlet port. Draw back a small volume to remove any air bubbles (a few µl):</p> <ol style="list-style-type: none"> Set a P1000 pipette to 200 µl Insert the tip into the inlet port Turn the wheel until the dial shows 220-230 µl, or until you can see a small volume of buffer entering the pipette tip <p><input type="checkbox"/> Using a P1000 pipette, flush 500 µl of the Priming Mix into the inlet port of the flow cell, avoiding the introduction of air bubbles</p> <p><input type="checkbox"/> Wait five minutes. During this time you can prepare your library for loading, as described in the next steps</p> <p><input type="checkbox"/> Repeat the priming step with another 500 µl flush</p> <p><input type="checkbox"/> Thoroughly mix the contents of the Loading Beads (LB) tubes by vortexing.</p> <p>In a new tube, prepare the library for loading as follows:</p> <p><input type="checkbox"/> 75 µl SQB</p> <p><input type="checkbox"/> 51 µl LB</p> <p><input type="checkbox"/> 24 µl DNA library</p> <p>Load your sample</p> <p><input type="checkbox"/> Load 150 µl of your sample through the inlet port</p> <p><input type="checkbox"/> Close the valve to seal the inlet port and close the PromethION lid when ready</p> <p><input type="checkbox"/> Wait a minimum of 10 minutes after loading the flow cells onto the PromethION before initiating any experiments. This will help to increase the sequencing output.</p>	
<p>Ending the experiment</p>	
<p><input type="checkbox"/> After your sequencing experiment is complete, if you would like to reuse the flow cell, please follow the Wash Kit instructions and store the washed flow cell at 2-8°C, OR</p> <p><input type="checkbox"/> Follow the returns procedure by washing out the flow cell ready to send back to Oxford Nanopore.</p>	

Direct cDNA Sequencing (SQK-DCS109)

Version: DCS_9090_v109_revL_14Aug2019
Last update: 22/12/2020



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>IMPORTANT</p> <p><input type="checkbox"/> If you encounter issues or have questions about your sequencing experiment, please refer to the Troubleshooting Guide that can be found in the online version of this protocol.</p>	