DOH Workshop Protocol Part 1: Manual Purification of High-Molecular Weight Genomic DNA from Gram-Negative Bacteria (MagAttract HMW DNA)

This protocol enables extraction of high-molecular-weight DNA from Gram-negative bacterial cultures and is an edited version of the Qiagen MagAttract HMW DNA Handbook protocol. Edited by Vesa Qarkaxhija and Bryan Wee

Starting Materials

□ 1.5 ml of overnight bacterial cell culture.

Equipment

- \square P1000 pipette (Micropipette with 100 1000 µl range)
- **□** P200 pipette (Micropipette with 20– 200 µl range)
- **D** P20 pipette (Micropipette with $2 20 \mu l$ range)
- **D** P10 pipette (Micropipette with $0.5 10 \mu l$ range)
- □ Thermomixer (for heating, cooling and mixing 1.5ml tubes)
- $\hfill\square$ Benchtop centrifuge for 1.5 μl (up to 5000 Gs) OR Bento Lab
- Minicentrifuge (for brief spin downs of 1.5ml eppendorfs and 0.2ml PCR tubes) OR Bento lab with improvised adaptor made from 0.5ml eppendorf within a 1.5ml eppendorf with caps cut off.
- □ PCR thermoblock (21°C 80°C required) OR Bento Lab
- DNA fluorometer (Promega Quantus or Themofisher Qubit)
- Vortex mixer

Consumables

- $\hfill \hfill \hfill$
- **D** P200 filter pipette tips (with 20– 200 μl range)
- **D** P20 filter pipette tips (with $2 20 \mu l$ range)
- **D** P10 filter pipette tip (with $0.5 10 \mu l$ range)
- □ 1.5 ml Eppendorf DNA LoBind tubes
- □ 0.2 ml thin-walled PCR tubes
- □ Absolute ethanol (>96%)
- Distilled water (1400 ml per sample)
- HMW Magattract Kit
 - Nuclease Free Water
 - Buffer MW1 (Ethanol must be added)
 - Buffer PE (Ethanol must be added)
 - MagAttract Suspension G
 - Proteinase K
 - RNAase A (100mg/ml)
 - o Buffer ATL

Things to do before starting:

- If precipitate has formed in Buffer ATL, dissolve by incubating at 37°C with occasional shaking. The presence of precipitate can result in incomplete lysis.
- Prepare (but do not start) the thermomixer to 56°C at 900rpm for 30min for the lysis step.

Important points throughout protocol:

- Be gentle when pipetting, too much fast pipetting/ extended vortexing can shear the DNA which will result in poorer sequencing outcomes.
- Do not let the beads dry after removal of supernatant. This can be avoided by keeping Eppendorf tubes closed if beads are not submerged. Note: When bead pellet is moist, it appears shiny. As they start to dry, the shine reduces, and cracks start to form.
- **SHAKE WELL** When working with beads ensure they are thoroughly mixed before using. This can be achieved by vortexing for at least 10 seconds before use and vortexing between use to prevent beads settling.

Start of protocol:

1.	Penet 1.5mm of the Grann negative bacterial cell culture by spinning at	
	5000 Gs / 5 kG for 3 min.	
	Note 1: When removing tube from centrifuge, be careful not to shake	
	the tube too much which can dislodge the pellet.	
2.	Remove and discard supernatant without disturbing the pellet. Use a	
	P200/P20 pipette if there is a small amount of supernatant remaining.	
	Note: The cell pellet can be stored at −30 to −15°C or −90 to −60°C for	
	future use, or can be used immediately	
3.	Resuspend the bacterial pellet in 180 μ l Buffer ATL, add 20 μ l	
	Proteinase K and mix by flicking/stirring the tube.	FRAGILE
4.	Incubate for 30 min on thermomixer at 56°C shaking at 900 rpm.	
Things	to do while waiting:	
Things •	to do while waiting: Ensure that Buffers MW1 and PE were prepared according to	
Things •	to do while waiting: Ensure that Buffers MW1 and PE were prepared according to instructions on them (i.e. adding appropriate amount of Ethanol).	
Things • 5.	to do while waiting: Ensure that Buffers MW1 and PE were prepared according to instructions on them (i.e. adding appropriate amount of Ethanol). Add 4 μl RNase A to the sample, mix by pulse vortexing or by tapping	
Things • 5.	to do while waiting: Ensure that Buffers MW1 and PE were prepared according to instructions on them (i.e. adding appropriate amount of Ethanol). Add 4 μl RNase A to the sample, mix by pulse vortexing or by tapping the tube several times, and incubate for 2 min at room temperature.	FRASILE
Things • 5. 6.	to do while waiting: Ensure that Buffers MW1 and PE were prepared according to instructions on them (i.e. adding appropriate amount of Ethanol). Add 4 μl RNase A to the sample, mix by pulse vortexing or by tapping the tube several times, and incubate for 2 min at room temperature. Set the thermomixer to 21°C to give it time to cool for next steps	FRAGILE
Things	to do while waiting: Ensure that Buffers MW1 and PE were prepared according to instructions on them (i.e. adding appropriate amount of Ethanol). Add 4 μl RNase A to the sample, mix by pulse vortexing or by tapping the tube several times, and incubate for 2 min at room temperature. Set the thermomixer to 21°C to give it time to cool for next steps Fully mix MagAttract Suspension G. Vortex the MagAttract Suspension	FRANLE
Things • 5. 6. 7.	to do while waiting: Ensure that Buffers MW1 and PE were prepared according to instructions on them (i.e. adding appropriate amount of Ethanol). Add 4 μl RNase A to the sample, mix by pulse vortexing or by tapping the tube several times, and incubate for 2 min at room temperature. Set the thermomixer to 21°C to give it time to cool for next steps Fully mix MagAttract Suspension G. Vortex the MagAttract Suspension G vigorously until thoroughly mixed ~10 seconds	
Things • 5. 6. 7.	to do while waiting: Ensure that Buffers MW1 and PE were prepared according to instructions on them (i.e. adding appropriate amount of Ethanol). Add 4 μl RNase A to the sample, mix by pulse vortexing or by tapping the tube several times, and incubate for 2 min at room temperature. Set the thermomixer to 21°C to give it time to cool for next steps Fully mix MagAttract Suspension G. Vortex the MagAttract Suspension G vigorously until thoroughly mixed ~10 seconds Note: It is important that beads are fully mixed and re-mixed between	FRAGILE SHAKE WELL
Things • 5. 6. 7.	to do while waiting: Ensure that Buffers MW1 and PE were prepared according to instructions on them (i.e. adding appropriate amount of Ethanol). Add 4 μl RNase A to the sample, mix by pulse vortexing or by tapping the tube several times, and incubate for 2 min at room temperature. Set the thermomixer to 21°C to give it time to cool for next steps Fully mix MagAttract Suspension G. Vortex the MagAttract Suspension G vigorously until thoroughly mixed ~10 seconds Note: It is important that beads are fully mixed and re-mixed between every few samples.	HAKE WELL

8.	Add 15 μ l MagAttract Suspension G and 280 μ l Buffer MB to the sample. Mix by pulse vortexing (1-2 second quick presses x 3).	FRAGILE
9.	Place the tube onto the thermomixer and incubate at 21°C for 3 min at 1400 rpm.	
10.	Place tubes onto magnetic rack and wait until beads pellet (~1min).	
11.	While leaving tubes on the magnet, carefully remove supernatant without disturbing the beads using a P1000.	FRAGILE
12.	Remove tube from magnetic rack and add 700 μ l Buffer MW1 to the tube, flick the tube until bead pellet is unstuck from tube wall.	MW1 Wash 1
13.	Place the tube onto the thermomixer and incubate at 21°C for 3 min at 1400 rpm.	
14.	Place tubes onto magnetic rack and wait until beads pellet (~1min).	
15.	While leaving tubes on the magnet, carefully remove supernatant without disturbing the beads using a P1000.	FRAGILE
16.	For the second wash, remove tube from magnetic rack and add 700 μ l Buffer MW1 to the tube, flick the tube until bead pellet is unstuck from tube wall.	MW1 Wash 2
17.	Place the tube onto the thermomixer and incubate at 21°C for 3 min at 1400 rpm.	
18.	Place tubes onto magnetic rack and wait until beads pellet (~1min).	
19.	While leaving tubes on the magnet, carefully remove supernatant without disturbing the beads using a P1000.	FRAGILE
20.	Remove tube from magnetic rack and add 700 µl Buffer PE to the tube, flick the tube until bead pellet is unstuck from tube wall.	PE Wash 1
21.	Place the tube onto the thermomixer and incubate at 21°C for 3 min at 1400 rpm.	
22.	Place tubes onto magnetic rack and wait until beads pellet (~1min).	
23.	While leaving tubes on the magnet, carefully remove supernatant without disturbing the beads using a P1000.	FRAGILE
24.	For the second wash, remove tube from magnetic rack and add 700 μ l Buffer PE to the tube, flick the tube until bead pellet is unstuck from tube wall.	PE Wash 2
25.	Place the tube onto the thermomixer and incubate at 21°C for 3 min at 1400 rpm.	
 26.	Place tubes onto magnetic rack and wait until beads pellet (~1min).	
27.	While leaving tubes on the magnet, carefully remove supernatant without disturbing the beads using a P1000. Then use a P200 to remove any traces of Buffer PE.	FRAGILE
28.	While keeping the tube on the magnetic rack, slowly pipette 700µl distilled water by pipetting down the side of the tube opposite the pellet, without disturbing any of the beads/pellet. Leave for 30-60 seconds and then remove the supernatant without disturbing any of the beads/pellet.	FRAGILE
29.	Repeat step 28. Note: <i>If necessary, use a P20 to remove any remaining distilled water.</i>	FRAGILE
30.	Remove the tube from the magnetic rack and add 50μ l of nuclease free water. Resuspend the beads by slowly pipetting or stirring with the pipette tip.	FRAGILE

Note: <i>Be as gentle as possible while ensuring that pellet is resuspended.</i>	
 Place the tube onto the thermomixer and incubate at 21°C for 3 min at 1400 rpm. 	
 32. Pellet beads on magnet and slowly pipette DNA eluate into a new lobind tube. Note: Save the tube with magnetic beads in case of incomplete DNA elution, so you can repeat. 	FRAGILE

Part 1.1: Promega Pronex protocol

This protocol is to purify extracted dsDNA, removing contaminants (e.g., buffers, proteins, salts, etc.) and low molecular weight DNA (e.g., dsDNA adapters, ssDNA oligonucleotides and nucleotides) and is an edited version of the Promega 6.A.Protocol. Edited by Vesa Qarkaxhija and Bryan Wee

Materials

□ Extracted genomic DNA

Equipment

- $\hfill\square$ P1000 pipette (Micropipette with 100 1000 μl range)
- **□** P200 pipette (Micropipette with 20– 200 µl range)
- **D** P20 pipette (Micropipette with $2 20 \mu l$ range)
- **D** P10 pipette (Micropipette with $0.5 10 \mu l$ range)
- DNA fluorometer (Promega Quantus or Themofisher Qubit)
- Vortex mixer

Consumables

- $\hfill\square$ P1000 filter pipette tips (with 100 1000 μl range)
- $\hfill\square$ P200 filter pipette tips (with 20– 200 μl range)
- **D** P20 filter pipette tips (with $2 20 \mu l$ range)
- \square P10 filter pipette tip (with 0.5 10 µl range)
- □ 1.5 ml Eppendorf DNA LoBind tubes
- **Qubit dsDNA HS Assay Kit OR Promega QuantiFluor® ONE Dye**
- ❑ Qubit[™] Assay Tubes
- □ Absolute ethanol (>96%)
- □ ProNex[®] Size-Selective Purification System
 - Pronex beads
 - Wash buffer (Ethanol must be added)

Things to do before starting:

• Ensure that the Pronex Wash Buffer is prepared according to instructions on them (i.e. adding appropriate amount of Ethanol).

Important points throughout protocol:

- Be gentle when pipetting, too much fast pipetting/ extended vortexing can shear the DNA which will result in poorer sequencing outcomes.
- Do not let the beads dry after removal of supernatant. This can be avoided by keeping Eppendorf tubes closed if beads are not submerged. Note: When bead pellet is moist, it appears shiny. As they start to dry, the shine reduces, and cracks start to form.



• **SHAKE WELL** When working with beads ensure they are thoroughly mixed before using. This can be achieved by vortexing for at least 10 seconds before use and vortexing between use to prevent beads settling.

Start of protocol:

1.	Resuspend the Pronex beads by vortexing for 10 seconds or longer	SHAKE WELL
2.	Into your extracted DNA tube pipette 80µl of Pronex beads and mix into	
	the sample by slowly pipetting 10 times.	
	Note: If sticky clumps of bead-bound DNA form, be careful not to take	
	any beads either in the pipette tip or on the outside of the pipette tip.	
3.	Leave at room temperature for 10 min.	
4.	Place the sample on a magnetic rack for 2 min (or longer if necessary)	
	until the solution becomes clear and the beads form a pellet.	
5.	While leaving tubes on the magnet, carefully remove and discard	
	supernatant without disturbing the beads using a P1000.	
6.	While the tube is on the magnetic rack, add 200µl of Pronex Wash	Wash Step 1
	Buffer without flushing directly onto the pellet. If $200\mu l$ is not enough	
	to submerge the pellet, use more Wash Buffer.	
7.	Allow to incubate for 30-60 seconds	
8.	While leaving tubes on the magnet, carefully remove and discard Wash	
	Buffer without disturbing the beads using a P1000.	
9.	For the second wash: Keeping the tube on the magnetic stand, carefully	Wash Step 2
	add 200 μ l of Pronex Wash Buffer without flushing directly onto the	
	pellet. If 200μ l is not enough to submerge the pellet, use more Wash	
	Buffer.	
10.	Allow to incubate for 30-60 seconds	
11.	While leaving tubes on the magnet, carefully remove and discard Wash	
	Buffer without disturbing the beads using a P1000.	
12.	Allow the sample to air dry with lids open for 2-5 min watching it until	کن
	the pellet is no longer shiny.	×
13.	Remove the sample from the magnetic stand.	

 14. Remove the tube from the magnetic rack and add 32µl of nuclease free water. Resuspend the beads by slowly pipetting or stirring with the pipette tip. Note: Be as gentle as possible while ensuring that pellet is resuspended. 	FRAGILE
15. Leave for 5 min at room temperature to elute the DNA	
 Pellet beads on magnet for 1 minute until solution becomes clear and slowly pipette DNA eluate into a new LoBind tube. Note: Save the tube with Pronex beads in case of incomplete DNA elution, so you can repeat. 	FRAGLE
17. Quantify 1µl of DNA elute using the Qubit™ dsDNA HS Assay Kit or QuantiFluor® ONE dsDNA System.	

Troubleshooting:

1. My DNA yield is low, do I have to redo the extraction?

There are a few things you can do to increase yield before having to redo the extraction. Try these solutions:

- Try re-eluting the DNA from the remaining Pronex beads from step 16 of the **Pronex protocol**. Repeat steps 14–16.
- Repeat Pronex protocol
- Try re-eluting the DNA from the remaining MagAttract Suspension G beads from step 32 of the **DNA extraction protocol**. Repeat steps 30, and re-elute by placing the tube onto the thermomixer and incubate at 21°C for **6 min** at 1400 rpm.

DOH Workshop Protocol Part 2: Library preparation for Rapid Sequencing DNA V14 Barcoding kit (SQK-RBK114.24) with Pronex modification

This protocol performs creation of Nanopore sequencing libraries for the MinION flow cell and is an edited version of the SQK-RBK114.24 protocol. Edited by Vesa Qarkaxhija and Bryan Wee

Material

D gDNA (~200ng- 300ng of DNA, maximum 10µl). If DNA is too concentrated, dilute with NFW to 10µl.

Equipment

- **D** P1000 pipette (Micropipette with $100 1000 \mu$ l range)
- **D** P200 pipette (Micropipette with 20– 200 μ l range)
- **D** P20 pipette (Micropipette with $2 20 \mu l$ range)
- **D** P10 pipette (Micropipette with $0.5 10 \mu l$ range)
- PCR thermoblock (21°C 80°C required) OR Bento Lab
- DNA fluorometer (Promega Quantus or Themofisher Qubit)

Consumables

- **D** P1000 filter pipette tips (with $100 1000 \mu$ l range)
- **D** P200 filter pipette tips (with 20– 200 μl range)
- **D** P20 filter pipette tips (with $2 20 \mu l$ range)
- **D** P10 filter pipette tip (with $0.5 10 \mu l$ range)
- □ 0.2 ml thin-walled PCR tubes x Number of Samples
- □ Nuclease free water (NFW)
- □ 1.5 ml Eppendorf DNA LoBind tubes
- □ ProNex[®] Size-Selective Purification System
 - Pronex beads
 - Wash buffer
- □ Nanopore Rapid sequencing V14 Amplicon sequencing (SQK-RBK114.24)
 - Rapid Adapter (RA)
 - Adapter Buffer (ADB)
 - AMPure XP Beads (AXP) (not used in this protocol)
 - Elution Buffer (EB)
 - Sequencing Buffer (SB)
 - Library Beads (LIB)
 - Flow Cell Flush (FCF)
 - Flow Cell Tether (FCT)
 - Rapid Barcodes x 24 (RB01-24)

□ Ice bucket with Ice

Important throughout the protocol:



Be gentle when pipetting, too much fast pipetting/ extended vortexing can shear the DNA which will result in poorer sequencing outcomes.

• Do not let the beads dry after removal of supernatant. This can be avoided by keeping Eppendorf tubes closed if beads are not submerged. Note: When bead pellet is moist, it appears shiny. As they start to dry, the shine reduces, and cracks start to form._



Be careful when pipetting into a Nanopore flow cell to avoid introducing bubbles into the channels. Bubbles can damage Nanopore array. Practice pipetting on a used flow cell.



• SHAFE WELL When working with beads ensure they are thoroughly mixed before using. This can be achieved by vortexing for at least 10 seconds before use and vortexing between use to prevent beads settling.

Things to do before starting:

1.	Calculate the dilution of your starting genomic DNA from the	
	extraction step to make up 10µl with a maximum concentration of	
	20 - 30ng/μl (200-300 ng total DNA). (For example, if your DNA from	
	the previous protocol was 100ng/µl, you would only need 2µl of	
	your sample, diluted in 8μl of NFW).	
2.	Assign each sample to a Barcode (1-24) and note this down	
3.	Program the thermal cycler to incubate at 30°C for 2 min then 80°C	
	for another 2 min. Do not start it yet.	

Start of protocol:

4.	Prepare one 0.2 ml thin-walled PCR tube for each sample from the previous step. Label the top of the tube with the barcode number	
5.	 In each 0.2 ml thin-walled PCR tubes: Pipette an appropriate amount of your sample (1-10μl) and add NFW (1-10μl), if necessary, to get 200-300ng of DNA in a 10μl total volume. 1.5 μl of your chosen Rapid Barcode (RB01-24). (1 barcode per sample) Note: Spin barcodes down for 2-3 seconds before use 	
6.	Mix by gently by stirring or pipetting until thoroughly mixed. If there are any bubbles present, spin down briefly (2-3 seconds) to ensure all liquid is at the bottom. Note: <i>Barcodes will be a thicker liquid, visually check they have been mixed in.</i>	FRAGILE

7.	Incubate the tubes in the thermal cycler (PCR machine) at 30°C for 2	
0	Priofly place the tubes on ise to see	
0.	Briefly place the tubes of ite to cool.	
9.	Pool all your samples into a 1.5ml clean Eppendorf DNA LoBind tube.	
	Note down the total volume after pooling. It should be 11.5µi	FRAGILE
Only 1	normal is required to carry out the following stops	
	person is required to carry out the following steps.	
10.	Resuspend the Pronex beads by vortexing for 10 seconds or longer.	
11.	Use a 1:1 ratio of sample to Pronex beads and mix into the sample	
	by slowly pipetting 10 times.	FRAGILE
	Note: If sticky clumps of bead-bound DNA form, be careful not to	
	take any beads either in the pipette tip or on the outside of the	
	pipette tip.	
12.	Incubate the sample at room temperature for 10 min.	
 To do v	vhile waiting	
•	Take out Elution Buffer (EB) to thaw on ice	
13.	Place the sample on a magnetic stand for 2 min until the solution	
	becomes clear and the beads form a pellet on one side of the tube.	
14.	While leaving the tube on the magnet, carefully remove and discard	
	supernatant without disturbing the beads.	
15.	Wash 1: While still on the magnetic stand, carefully add 200 μ l of	
	Pronex Wash Buffer without flushing directly onto the pellet. If	
	200µl is not enough to submerge the pellet, use more Wash Buffer.	
16.	Allow to incubate for 30-60 seconds	
17	Milette land the Aules and the second strength and fully second attended	
17.	while leaving the tube on the magnet, carefully remove and discard	
10	Wash 2: While still on the magnetic stand, carefully add 200ul of	
10.	Proper Wash 2. While still on the magnetic stand, carefully add 200µ1 of	
	200ul is not enough to submerge the pellet use more Wash Buffer	
19	Allow to incubate for 30-60 seconds	
20	While leaving the tube on the magnet, carefully remove and discard	
20.	supernatant without disturbing the heads	
21	Allow the sample to air dry with lids onen for 2-5 min watching it	<u></u>
21.	until the pellet is no longer shiny.	-20:-
22.	Remove the tube from the magnetic rack and add 15ul of Elution	
	Buffer (included in the Nanopore Rapid Barcoding kit). Resuspend	FRAGILE
	the beads by slowly pipetting or stirring with the pipette tip.	
	Note: Be as gentle as possible while ensuring that pellet is resuspended	
23.	Incubate the samples at room temperature for 10 min to elute the	
	DNA.	
24.	Return the tube to the magnetic stand for 1 minute until the	
	solution becomes clear and the beads form a pellet.	
25.	Store the Pronex beads in the fridge.	

26. Transfer 11μ l of the eluate into a clean 1.5ml Eppendorf DNA LoBind	
tube.	FRAGILE
27. Take enother 1.1 of the alute from the tube on the meanatic stand	AND CONTRACTOR
27. Take another 1 μ i of the eluce from the tube on the magnetic stand	
for quantification on a fluorometer (Qubit or Quantus). The	
remaining beads can be kept in a closed tube on ice, for re-elution, if	
necessary.	
28. In a new 1.5 ml Eppendorf DNA LoBind tube, mix the following:	
 1.5µl Rapid Adapter (RA) 	
 3.5μl Adapter Buffer (ADB) 	
29. Add 1l of this $PA + ADR$ mixture to the DNA	
20. Mix conthe building and anin down briefly (2.2 cocondo)	
30. IVIX gently by flicking and spin down briefly (2-3 seconds)	
31. Incubate this for 30 min at room temperature.	
Preparing the flowcell	
32 Remove the following Nanonore Ranid Kit (RRK-114 24) items from	
the 20 °C freezer spin down and store on ice	
the -20°C freezer, spin down and store on ice.	
• SB (Sequencing Burler)	
 LIB (LIDrary Beads) 	
• FCT (Flow Cell Tether)	
• FCF (Flow Cell Flush)	
 Bovine Serum Albumin (BSA) at 50mg/ml 	
33 Prenare the flow cell Priming Mix in a fresh 1 5ml Ennendorf DNA	
LoBind Mix by inverting the tube	
0 1,170µ FCF	
24. Remove the flow cell you want to use and slide it under the metal	
steriove the now cell you want to use and side it under the metal	
correct contact with the thermal and electrical connections	
Correct contact with the thermal and electrical connections.	
 MK1B: Plug in the MiniON MK1B to a laptop with Minknow 	
software	
• Mk1C: Turn on the MinION Mk1C.	
35. Complete a flow cell check to assess the number of pores available	
On the now cell.	
36. Rotate the flow cell priming port cover clockwise to open the	
priming port.	~
37. After opening the priming port there will be a small air bubble under	
the cover that needs to be removed.	
\circ Set a P1000 pipette to 200 µl	
 Insert the tip into the priming port 	
\circ Turn the adjustment wheel slowly, pausing every few µls,	
until the pipette shows 220 μ l - 230 μ l to draw a total of 20-	
30 μ l out of the priming port, or until you can see a small	
volume of liquid entering the pipette tip.	
Note 1: There may be a small delay before the liquid comes out of	
the port into the pipette tip. Do not draw out more than 30 μ l.	
Note 2: Check that there is a continuous flow of buffer from the	



