

Positive Control Serum (Titer):

Negative Control Serum:

Conjugate assessed:

Conjugate Lot No:

Conjugate dilution	Control sera	16	32	64	128	256	512	1024	2048	4096
1/50	Pos.									
1/50	Neg.									
1/100	Pos.									
1/100	Neg.									
1/200	Pos.									
1/200	Neg.									
1/400	Pos.									
1/400	Neg.									
1/800	Pos.									
1/800	Neg.									

RESULTS:

Plateau Endpoint Titer of Positive Control Serum: _____

Conjugate Endpoint Titer: _____

Optimal Conjugate Dilution*: _____

*Taken to be within one twofold dilution from the conjugate endpoint titer, as long as negative control serum is not reactive at the lowest sample dilutions using this conjugate dilution.