

Supplementary information

Viral crosslinking and solid-phase purification enables discovery of ribonucleoprotein complexes on incoming RNA virus genomes

In the format provided by the authors and unedited

Supplementary Figure 1



Supplementary Figure 1. CLASP with RNA marker input.

1% Agarose gel of RNAs pulled down with CLASP performed with the RNA marker (lane 1) as input. Lanes 2-5 demonstrate that increasing the bead buffer ratio (bead buffer : sample buffer) results in purification of lower molecular weight RNAs.