

COBRA Protocol Tray Map: Hua protocol.txt

1	2	3	4	5	6	7	8	9	10	11	12
								Rinse in buffer	Rinse in buffer	Rinse in buffer	
Plate 1											

1	2	3	4	5	6	7	8	9	10	11	12
OsO4 in buffer	Potassium ferricya		Rinse in water	Rinse in water	Rinse in water		OsO4 in water		Rinse in water	Rinse in water	Rinse in water
Plate 2											

1	2	3	4	5	6	7	8	9	10	11	12
							Rinse in water	Rinse in water	Rinse in water		TCH
Plate 3											

1	2	3	4	5	6	7	8	9	10	11	12
Watson's lead		Rinse in water	Rinse in water	Rinse in water							
Plate 4											

1	2	3	4	5	6	7	8	9	10	11	12
100% Resin	100% Resin	100% Resin	75% Resin	50% Resin	100% Acetone	100% Acetone	100% Acetone	95% Acetone	85% Acetone	75% Acetone	50% Acetone
Plate 5											

1	2	3	4	5	6	7	8	9	10	11	12
Uranyl acetate		Rinse in water	Rinse in water	Rinse in water							
Plate 6											